

Brewer

MODEL 1500 BLOOD DRAWING CHAIR

Assembly and Operation Manual



<i>Important Information</i>	2	<i>Operation</i>	5
<i>Safety Information</i>	2	<i>Maintenance</i>	5
<i>Specifications</i>	3	<i>Options, Accessories</i>	6
<i>Assembly</i>	3-4	<i>Limited Warranty</i>	6

Brewer
more possibilities

The Brewer Company
13901 Main Street
Menomonee Falls, WI 53051
P 1.888.Brewer.1
F 262.251.1786
www.brewercompany.com

IMPORTANT INFORMATION

General

The Brewer Blood Drawing Chair is designed to provide positioning and support of patients during the blood drawing process conducted by qualified medical professionals. Read and understand all safety information, operating instructions and maintenance requirements contained in this manual prior to operating the chair. Become familiar with all of the chair's functions before using it. The chair is constructed of tubular steel, welded for strength and longevity.

Service

If you require assistance with the installation or operation of your Brewer Blood Drawing Chair, call the Brewer Customer Service Department at **1-888-Brewer1**. Our trained staff will attempt to assist you in correcting the problem directly over the phone. If service is required, a factory authorized technician will be sent to your location. Please fill in the following information for use when calling the Brewer Company or your distributor with questions regarding your unit. See Figure 1 for model and serial number location.


Date of Purchase _____
Model Number _____
Authorized Dealer Name _____
Dealer Phone Number _____

SAFETY INFORMATION

Safety Information


The primary concern of The Brewer Company is that the equipment is operated and maintained with the safety of the patient and healthcare staff in mind. To ensure safe and reliable operation:

- Read and understand all instructions in this manual before attempting to install or operate the chair.
- Ensure that appropriate personnel are informed on the manual contents. This is the responsibility of the purchaser.
- Ensure that this manual is located near the chair, or if possible permanently affixed to the chair.

Failure to follow instructions with this symbol  could result in serious personal injury and/or damage to the unit.

 **DANGER** indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

 **WARNING** indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

 **CAUTION** indicates an imminently or potential hazardous situation which, if not avoided, may result in serious injury and/or equipment damage.

SPECIFICATIONS

Weight of Chair.....	44 lbs.	Hip Clearance.....	20 in.
(With Carton).....	52 lbs.	Seat Cushion Dimensions.....	18 in. x 18 in.
Chair Weight Capacity.....	350 lbs.	Back Cushion Dimensions.....	18 in. x 12.5 in.
Chair Height (Overall).....	37 in.		
Chair Depth (Overall).....	23.5 in.		
Width:			
With Both Worksurfaces attached.....	43.25 in.		
Armrests Only.....	29.25 in.		
Inside Arm Width.....	19.25 in.		
Seat Height.....	22.0 in.		

ASSEMBLY

Overview

Perform the following sequence when setting up the chair.

1. Unpack the Chair.
2. Assemble the Chair.
3. Position Chair in desired location.

Note: Inspect all boxes and contents for damage. Report any damage to the carrier immediately.

Note: Save all packaging in case re-shipment is required.

⚠ CAUTION

To avoid damaging the Chair's upholstery or painted surfaces, DO NOT use a knife or any other sharp object to open the packaging.

⚠ CAUTION


The chair weighs approximately 44 lbs. Two or more people should assist in removing the chair from the shipping carton. Also, use proper lifting techniques when lifting. Failure to do so could result in serious injury.




Chair shown in packed configuration. Armrests are reversed for shipping

ASSEMBLY (continued)

1. Remove Armrests from Chair and set aside. (The Armrests are installed at the factory in a shipping configuration that is opposite of their intended position). See figure on previous page.
2. Insert (4) 1/4-20 x 1.75" screws and washers through seat supports from underneath chair into Seat Bottom and tighten. Seat Bottom should be positioned with Notch and Label facing the front of the chair. (See inset picture to locate front of Seat Bottom)
3. Insert (4) 1/4-20 x 0.75" screws and washers through Tabs on upper frame, from back of chair into Seat Back and tighten.
4. Position Armrests back into tubing, install Hand Knob and tighten.




No Work Surface Option: To remove the Work Surface portion of the Armrest, remove (4) screws, slide Work Surface and its support plate away from Armrest, then reinstall the (4) screws into the Armrest.



OPERATION

Adjusting the Armrest

1. While supporting the Armrest to be adjusted, loosen the Hand Knob.
2. Move the Armrest vertically to the desired position and tighten Hand Knob.
3. To remove the Armrest completely from the chair loosen the Hand Knob until the armrest lifts freely out of the square tube.



MAINTENANCE

Preventative Maintenance

- Inspect mechanical functions to ensure satisfactory operation.
- Check fasteners to make sure they are present and secured tightly.

Cleaning

For Light Soiling:

1. A solution of 10% household liquid dish soap with warm water applied with a soft damp cloth will remove most soiling.
2. If necessary, use a solution of liquid cleanser and water applied with a soft bristle brush. Wipe away the residue with a warm damp cloth.

For more difficult stains not removed by the above method:

1. Dampen a soft white cloth with a solution of household bleach (sodium hypochlorite); 10% bleach, 90% water.
2. Rub gently.
3. Rinse with a water dampened cloth to remove residual bleach.

All surfaces should be cleaned within your facility's guidelines concerning antimicrobials and blood borne pathogens.

Painted Metal Surfaces

Wipe all painted metal surfaces with a damp cloth or sponge at least once a week.

OPTIONS, ACCESSORIES, AND COMMON SERVICE PARTS

Description	Order No.	User/Restrictions
Replacement Upholstery Seat Back Seat Bottom	1515-XX 1550-XX	Intended Use-Replacement upholstery with hardware for 1500 Blood Draw Chair in 18 Standard colors.
Special Upholstery Seat Back Seat Bottom	1515-SP 1550-SP	Intended Use-Allows the purchase of specialty vinyls in different colors. (Additional charges may apply)
Work Surface	101244	Work surface adjacent to Armrest.
Armrest	101245	Right or Left hand patient Armrest.
Hand Knob	08082	Intended use-Allows Armrest to be adjusted.
Drawer Option	101059	Optional storage drawer for blood drawing supplies.
Tube Plugs	100506	Tube caps for bottoms of chair legs (6)

LIMITED WARRANTY

**The Brewer Company
GENERAL TERMS AND CONDITIONS**

Warranty: The Brewer Company warrants its Blood Draw Chair to be free from defects in parts and workmanship under normal use and service for a period of five (5) years from date of shipment. The Brewer Company will not be responsible for any Product failure due to abuse, misuse, modification or improper use or for any use which exceeds the published capacity of the Product. Products returned by prepaid freight for inspection and found defective will, at the option of The Brewer Company, be repaired or replaced at no charge, but no claim for outside labor or other charges will be allowed. Products must not be returned without proper written authorization from The Brewer Company. Requests for authorization must be in writing and accompanied by the original purchase order, The Brewer Company invoice number and a copy of the invoice for the Product. THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES AND REMEDIES WHATSOEVER, INCLUDING BUT NOT LIMITED TO IMPLIED WARRANTIES OF MERCHANTABILITY AND/OR FITNESS FOR A PARTICULAR PURPOSE. NO AGENT, EMPLOYEE, OR REPRESENTATIVE OF THE BREWER COMPANY HAS ANY AUTHORITY TO MAKE ANY AFFIRMATION, REPRESENTATION, OR WARRANTY NOT SET FORTH IN THESE TERMS AND CONDITIONS CONCERNING ANY PRODUCTS OF THE BREWER COMPANY. THE BREWER COMPANY SHALL HAVE NO LIABILITY WHATSOEVER FOR DAMAGES CAUSED BY TRANSPORTATION, ACCIDENTS, FIRE, UNAUTHORIZED ALTERATION, OR NORMAL WEAR OR ABUSE, NOR SHALL THE BREWER COMPANY HAVE ANY LIABILITY WHATSOEVER FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES, including without limitation, lost profits or any such damages arising from (a) the design, manufacture, sale, delivery, installation, repair, operation or use of products of The Brewer Company or any part thereof, (b) any actual or alleged failure or defect in products of The Brewer Company or any part thereof, or (c) any actual or alleged breach or non-performance by The Brewer Company of this limited warranty.

Freight Claims: Upon receipt, merchandise should be carefully examined to ascertain that proper amount has been received in good condition. Any claim for shortage or damage must be made with delivering carrier within five (5) days from receipt of goods. We do not assume any responsibility for loss or damage in transit, and compensation for such loss must be obtained from the carrier.

NOTES