

LiftMate™

Brewer Low/High Patient Lift

Lift with your finger, not your back.



The Brewer LiftMate™ is an affordable, low/high mobile patient lift designed to decrease the injury risk associated with patients getting onto and off of standard exam tables. The first product of its kind, Brewer LiftMate reduces patient falls and clinical staff injury caused by manually lifting patients.

Safe patient handling has never been so easy or affordable.

- Helps prevent patient falls
- Improves clinical staff safety
- Removes patient mobility barriers
- Reduces capital investment
- Mobile, maneuverable and easy to use



Brewer

Furnishing solutions for safe practices.

www.brewerliftmate.com

LiftMate™

Brewer Low/High Patient Lift

Durability

- Three-year warranty

Product Safety

- 500 lb. weight capacity
- Large non-slip platform surface
- Two emergency stop buttons
- Grab bars for patient stability and positioning
- Four auto-lock rubber feet for stability
- Four retractable swivel casters prevent rolling during use

Function

- Convenient hand pendant for easy operation
- Powered by standard outlet plug or rechargeable battery
- Easy-to-clean design
- Fits through standard doorways

Patient Safety

- Smooth movement for enhanced patient security
- Quiet actuators for patient comfort
- Large platform with ample room for patient positioning
- Easy access for walkers
- Designed to prevent patient falls

Caregiver Benefits

- Easy-to-use mobile device
- Reduces need to assist patients onto exam tables
- Protects employees from workplace hazards

Dimensions

- Height: 38" - 47"
- Depth (outer dimensions): 21.5"
- Width (outer dimensions): 34.75"
- Depth (platform): 18"
- Width (platform): 21"
- Width (between grab bars): 27.75"

Electrical

- 115 VAC, 60Hz
- Current in: 2.0 amp max
- Class II double insulated
- Safe, extra-low voltage actuators

Battery

- Lead Acid Sealed Battery
- 24 volt
- 1.2 amp-hour capacity
- Rechargeable
- Light/Audible signal to indicate low battery charge
- Replacement Part No.: 2102150

Certifications

All Brewer products conform to U.S. and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

