



Brewer

BLOOD DRAWING CHAIR

BLOOD DRAWING CHAIR SERIES

Comfort. Adjustability. Durability.

Durability

- 5 year warranty

Safety

- 350 lb. weight capacity*
- Work surfaces covered in a bacteriological, chemical, and fungal resistant material
- Meets BIFMA (Business and Institutional Furniture Manufacturer's Association) standards

Function

- Innovative angled armrests allow optimal arm position for blood drawing
- Durable, powder-coated welded steel tube frame
- Contoured, seamless upholstery for easy cleaning
- Two standard add-on work surfaces to accommodate prep materials, paperwork, and vials. (10" D x 7" W)
- Work surfaces are easily removable to accommodate space constraints
- Taper-lock armrest knob adjustment system ensures smooth movement

Patient Comfort

- 2" thick foam for patient comfort
- Seat alignment positions patients securely on chair

Aesthetics

- 15 standard upholstery color options
- Custom colors available

Chair Dimensions with upholstery

- Height: 37"
- Length: 23.5"
- Width: 20"; 29.25" (with armrests); 43.25" (with both work surfaces)
- Floor to top of Seat: 22"
- Chair Depth: 23.5"

Armrest Function

Angle Degree: 15°

Height Adjustment: 29.5" - 38"

Upholstery Dimensions

Seat cushion: 18" L x 18" W

Backrest cushion: 18" L x 12.5" W

*under normal use

Two standard add-on work surfaces accommodate prep materials, paperwork, and vials



OPTIONAL BLOOD DRAWING CHAIR DRAWER

Function

- Easily attaches to either side of chair
- Work surfaces covered in a bacteriological, chemical, and fungal resistant material
- Interior plastic lining for easy cleaning

Optional Drawer Dimensions

- Height: 4"
- Width: 8"
- Length: 18"



#101059

All Brewer blood drawing chairs come with a 5-year warranty and are available in 15 standard colors.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
N88 W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Front cover: 1500 in Emerald

Brewer

07.2012
00020