



Brewer

STEP STOOLS

STEP STOOLS

SERIES Stability with exceptional durability.



Step Stool with Hand Rail

- Chrome plated 18 gauge steel for durability
- Skid resistant rubber mat for safety
- Dimensions: 8.75" H x 11.25" W x 14.25" L
- Handrail is .88" chrome plated steel
- Height: 35.75" from top of rail to floor
- Weight capacity: 350 lbs.
- Some assembly required

#11220

accessories

#99121 replacement crutch tips (set of 4)



Step Stool

- Chrome plated 18 gauge steel for durability
- Skid resistant rubber mat for safety
- Dimensions: 8.75" H x 11.25" W x 14.25" L
- Weight capacity: 350 lbs.

#11200-1
1 stool per box

#11200-2
2 stools per box

accessories

#99121 replacement crutch tips (set of 4)

Two-Step, Step Stool with Hand Rail

- Chrome plated 16 gauge steel for durability
- Skid resistant rubber mat for safety
- Dimensions (first step): 9" H x 12" W x 16" L
- Dimensions (second step): 7" H x 12" W x 16" L
- Handrail is 1" chrome plated steel (can put on either side)
- Height: 43" from top at highest point of rail to floor
- Weight capacity: 600 lbs.
- Some assembly required

#31200
step stool

#31220
step stool w/hand rail

accessories

#50016-1 handrail



Bariatric Step Stool with Hand Rail

- Chrome plated 16 gauge steel for durability
- Skid resistant rubber mat for safety
- Dimensions: 8.75" H x 11.25" W x 14.25" L
- Handrail is 1" chrome plated steel
- Height: 34.5" from top of rail to floor
- Weight capacity: 600 lbs.
- Some assembly required

#21200
step stool

#21220
step stool w/hand rail

accessories

#99121 replacement crutch tips (set of 4)



All Brewer step stools come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
N88 W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Brewer