

BrewerAccess High-Low

EXAM TABLE
700



Advanced patient care for the health of your practice.

- 700 lb weight capacity!
- Unique safety features
- Maximum storage and workspace
- Easy, safe patient movement



Explore the Access High-Low family here. →



Chair-To-Table Design

- Smoothly moves patient from consultation to examination
- 18" low height for easy and safe patient transfers
- Return to chair feature easily takes table back to chair position in one simple step



Standard Patient Transfer Grab Bars

- Provides patients a stable and secure grab point to use when getting on and off table
- Conveniently swings out for lateral wheelchair transfers
- Easily removable



Pass-Through Work Surface

- Pulls out on either side
- Saves space by eliminating the need for additional stands or trays
- Perfect for charts, instruments, procedure preparation, or portable computer



Abundant Storage Capacity

- Patented pass-through side drawers allow left or right access
- Additional spacious front drawer
- Convenient covered paper roll storage holds up to six rolls



Comfortable, Ergonomic Upholstery

- Seamless, medicinal stain resistant
- Available in 15 standard colors and UltraLeather™ plush upholstery
- Cal-133 fire code compliant upholstery available



Reliable, Durable Design

- Industry-leading standard 3-year warranty
- Easy to operate and service
- UL Approved

Brewer

Furnishing solutions for safe practices.

BrewerAccess High-Low EXAM TABLE 700

Advanced Safety

- Safety grab bars for patient stability and positioning
- Return to chair feature
- Anti-tamper on/off switch prevents accidental movement
- Sensors prevent injuries during table movement
- 18" (45.7 cm) low height, 28" (71.1 cm) when flat
- Easy on / off for patients
- Powered backrest adjusts for easy patient positioning

Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density support upholstery
- Standard adjustable ergonomic stirrups
- Toe board / leg extension / stainless steel debris tray

Abundant Storage

- Pass-through work surface on side of table
- Patented pass-through side drawers accessible from right or left side of the table
- Front drawer: 19.8" x 19.2" x 3.8" (50.3 cm x 48.8 cm x 9.7 cm)
- Side drawers: 22.6" x 20.2" x 7" (57.4 cm x 51.3 cm x 17.8 cm)
- Covered storage for up to 6 additional paper rolls
- Paper roll storage capacity: 23" x 23" x 3.25" (58.4 cm x 58.4 cm x 8.3 cm)
- Slide through holder accepts 21" x 3" (53.3 cm x 7.6 cm) and 18" x 3" (45.7 cm x 7.6 cm) paper rolls

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

- UL 60601-1



Table Specifications

- 700 lb (317 kg) weight capacity
- Standard patient assist grab bars
- Standard adjustable stirrups
- Standard foot control
- Standard 3-year warranty
- Height: 18" - 37" (45.7 cm - 94 cm)
- Length: 57" extends to 78" (144.8 cm - 198.1 cm)
- Upholstery width: 28" (71.1 cm)
- Leg extension length: 21" (53.3 cm)
- Goes flat at 28" (71.1 cm)
- 115VAC, 60 Hz, 5.8 amps max
- Safe, low voltage actuators
- Safe, single power cord design, detachable
- Easy to clean, seamless upholstery
- 15 standard upholstery color options
- Cal-133 fire code compliant upholstery available in 4 colors
- Premium, plush upholstery available in 5 colors
- Custom colors also available

6801 model only:

- Hospital grade outlet standard
- 115VAC, 60 Hz, 9.3 amps max
- 5° pelvic tilt / 60-Watt drawer heater in front top drawer

models

- | | |
|-------|---|
| #6800 | Access High-Low 700 Exam Table with power back |
| #6801 | Access High-Low 700 Exam Table with power back, pelvic tilt, drawer heater and outlet |

accessories

- | | |
|---------|--|
| #98464 | stainless steel treatment pan |
| #99502 | one universal socket clamp |
| #99504* | deluxe articulating knee crutches |
| #100253 | urology drain pan system |
| #100407 | front drawer dividers (can order up to 7) |
| #100540 | two support adapters for knee crutches (plant installed) |
| #100550 | two support adapters for knee crutches (field installed) |
| #102031 | Brewer Welch Allyn® IV bracket (plant installed) |
| #102032 | Brewer Welch Allyn® IV bracket (field installed) |

* Needs two #99502 and one of either #100540 or #100550

The Brewer Company is an ISO 13485:2003 certified company.

All Brewer Access™ High-Low Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company LLC
N88W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 | www.brewercompany.com

Brewer