

# Brewer

# Access High-Low

## EXAM TABLE

### ADVANCED PATIENT CARE FOR THE HEALTH OF YOUR PRACTICE.

- Increase your patient satisfaction rankings by preserving patient dignity and comfort during exams and procedures.
- Keep pace with the ever changing patient population requirements by accommodating a wider range of patients with a single table.
- More standard features. **Pivoting safety grab bars, laptop work surface, and pass-through drawers** equate to more storage and increased flexibility.
- Reduce your total cost of ownership with our industry-leading **three year warranty** for peace of mind.



# Access High-Low

## EXAM TABLE

Advance staff safety and care for more of your patients.



### ABUNDANT STORAGE

- Patented pass-through side drawers eliminates left or right setup issues
- Pass-through laptop surface
- Convenient table paper storage under backrest

### ADVANCED SAFETY

- 450 & 700 lb weight capacities
- Anit-tamper on/off switch
- Adjustable, ergonomic stirrups integrated in table
- UL 60601 certified

### SUPERIOR DESIGN

- Takes patient safely and automatically from consultation to examination
- 18" low height for easy and safe patient transfers
- Made in USA

### GRAB BAR TRANSFERS

- Stable and secure grab point for wheelchair and other transfers
- Conveniently swings out and are easily removable

### ERGONOMIC UPHOLSTERY

- Seamless & easy to clean
- Multi-density, bolstered support upholstery
- 15 standard colors available

#### models

#6000	Access High-Low with pneumatic back
#6001	Access High-Low with pneumatic back, pelvic tilt, and front drawer warmer
#6500	Access High-Low with power back
#6501	Access High-Low with power back, pelvic tilt, front drawer warmer and outlet
#6800	Access High-Low 700 with power back
#6801	Access High-Low 700 with power back, pelvic tilt, front drawer warmer and outlet

#### accessories

#99464	stainless steel treatment pan
#99502	one universal socket clamp
#99504	deluxe articulating knee crutches (pair) *
#100253	urology drain pan system
#100407	front drawer dividers (can order up to 7)
#100540	two knee crutch support adapters (plant install)
#100550	two knee crutch support adapters (field install)
#101800	electrical outlet (6001 only, plant install)
#102031	Brewer Welch Allyn® IV bracket (plant install)
#102032	Brewer Welch Allyn® IV bracket (field install)

\* Needs two #99502 and one of either #100540 or #100550

