


## Backrest Upright Bar Knob Differences



Glossy Finish

Zinc or Nickel Plated Screw


Coarse Threads

Narrower Screw (slightly less than 3/8")

### Type A

\*For 21340B Millenium Surgeon Stool

\*Effective 21-Feb-2013. Subject to change. Contact Brewer Customer Service to make sure this knob has not changed.



Matte Finish

Black Oxide Coated Screw

Fine Threads

Wider Screw (3/8")

### Type B

\*For 21520B Millenium Lab Stool, 21435B Millenium Task Stool, 11001B Century Stool, and 22400B & 22500B Value Stool

\*This is for all iterations of stools with backrests.

For example:

11001B

11001BV11001BLCD, etc.