

Brewer

ACCESS[®] HIGH-LOW PLUS

POWER EXAM TABLE



*Access High-Low PLUS
shown in Blue Fog*

Discover the Access High-Low PLUS Power Exam Table, fully equipped to meet every clinical need. Features such as ComfortGrip™ Patient Assist Handles, expansive storage, and adjustable ergonomic stirrups enhance clinical efficiency and patient care. Featuring a lighter, modern cabinet color that enhances the contemporary aesthetic of any clinic. The optional EasyGlide™ Mobility System streamlines exam table movement, facilitating effortless room reconfiguration and enabling easy cleaning. Embrace unmatched versatility and efficiency with this all-in-one solution.



ComfortGrip™ Patient Assist Handles



EasyGlide™ Mobility System



Industry Best Storage

- ADA Compliant seat height 18"-37"
- Length: 60"; extendable to 78"
- Upholstery width: 28"
- Front storage drawer: 19^{3/4}" x 19^{1/4}" x 4"
- Side storage drawer: 22^{1/2}" x 20^{1/4}" x 8"
- Paper storage capacity: 23" x 23^{1/4}" x 3^{1/4}"
- Industry leading 3 safety zones
- ComfortGrip™ Patient Assist Handles included
- Mayo Pass-Through Work Station included
- 6 cubic ft. of storage, 3 times more than others

FULLY FEATURED, PERFECTLY ENHANCED



Models	Weight Capacity	ADA Compliant 18" Seat Height	Power Base	Power Back	Manual Back	Patient Assist Handles	Stirrups	Pelvic Tilt & Drawer Warmer	Outlet	Work Surface	Mobility
6000PLUS	450	•	•		•	•	•			•	optional
6001PLUS	450	•	•		•	•	•	•		•	optional
6500PLUS	450	•	•	•		•	•			•	optional
6501PLUS	450	•	•	•		•	•	•	•	•	optional
6800PLUS	700	•	•	•		•	•			•	optional
6801PLUS	700	•	•	•		•	•	•	•	•	optional

PolyCloud Contour Upholstery

Maximizes patient comfort with a unique blend of adaptive foam and ergonomic design.

Up to 700 lb. Weight Capacity

Ensuring safety for patients of all sizes.

Brewer's Storage Hub

Industry's best storage - 5 cubic ft. of storage streamlines supplies management and boosts clinical efficiency.

ComfortGrip™ Patient Assist Handles

Facilitate easy transfers and adjustments, enhancing patient independence and safety during examinations.

Enhanced handle gripping surface is designed to ADA 2010 Standards - 505.7.2 Non-Circular Cross section.

Patented Pass-Through Drawers

Streamline care with easy-access, dual-side pass-through drawers—maximizing convenience and functionality.

Ergonomic Stirrups

Support women's health exams and treatments comfortably.

Mayo Pass-Through Work Surface

Conveniently holds laptops and supplies, enhancing clinical productivity.

ADA Compliant Seat Height

ADA compliant seat height lowers to 18 in. for simpler wheelchair transfers and patient accessibility.

EasyGlide™ Mobility System

Provides a single-step, whisper-quiet mobility solution for effortless room reconfiguration and improved cleaning effectiveness.



Available Colors

Need a closer look? Contact us to receive physical samples of our vinyl or base materials directly to you!

