

Brewer

FLEX™ HIGH-LOW PLUS

POWER EXAM TABLE



*FLEX High-Low PLUS
shown in Blue Fog*

Upgrade to the FLEX High-Low PLUS Power Exam Table, where advanced features meet affordability, enhancing your practice without a premium price. Designed for cost-effectiveness, this versatile table allows for customization with optional features like the innovative EasyGlide™ Mobility System and ComfortGrip™ Patient Assist Handles, ensuring superior mobility, functionality, and accessibility tailored to your clinical needs and budget. The lighter, modern cabinet color enhances the contemporary aesthetic of any clinic.



ComfortGrip™ Patient Assist Handles



EasyGlide™ Mobility System



Ergonomic Stirrups

- ADA Compliant seat height 18"-37"
- Length: 60"; extendable to 78"
- Upholstery width: 28"
- Industry leading 3 safety zones
- Industry best 3 year warranty
- Front storage drawer: 19³/₄" x 19¹/₄" x 4"
- 0.9 cubic ft. of storage
- Optional ComfortGrip™ Patient Assist Handles
- Optional Ergonomic Stirrups
- Optional Mayo Pass-Through Work Surface

INTUITIVELY DESIGNED, SEAMLESSLY EVOLVED



Models	Weight Capacity	ADA Compliant 18" Seat Height	Power Base	Power Back	Manual Back	Assist Handles	Stirrups	Pelvic Tilt & Drawer Warmer	Outlet	Work Station	Mobility
5700PLUS	700	•	•		•	optional	optional			optional	optional
5701PLUS	700	•	•		•	optional	•	•	•	optional	optional
5800PLUS	700	•	•	•		optional	optional			optional	optional
5801PLUS	700	•	•	•		optional	•	•	•	optional	optional

PolyCloud Contour Upholstery

Maximizes patient comfort with a unique blend of adaptive foam and ergonomic design.

ComfortGrip™ Patient Assist Handles

Facilitate easy transfers and adjustments, enhancing patient independence and safety during examinations.

Enhanced handle gripping surface is designed to ADA 2010 Standards - 506.7.2 Non-Circular Cross section.

ADA Compliant Seat Height

ADA compliant seat height lowers to 18 in. for simpler wheelchair transfers and patient accessibility.

Ergonomic Stirrups

Support women's health exams and treatments comfortably.

Mayo Pass-Through Work Surface

Conveniently holds laptops and supplies, enhancing clinical productivity.

EasyGlide™ Mobility System

Provides a single-step, whisper-quiet mobility solution for effortless room reconfigurations and improved cleaning effectiveness.



Available Colors

Need a closer look? Contact us to receive physical samples of our vinyl or base materials directly to you!

