

Brewer 135DSS

Dynamic Saddle Stool

Optimize posture and production with Brewer's patented breakthrough – Dynamic Split-Saddle-Seat Technology.

Patented Split Saddle Seat improves comfort

- Eases the load on back muscles by distributing weight over “sits- bones”
- Wider hip abduction lessens the tendency to slump
- Helps cool the body, relieve pressure, and improve circulation

Unique circular 360° movement strengthens core muscles

- Dynamic movement offers full 360° tilt to maintain proper neutral posture while you work
- Lessens spinal muscle degeneration and back pain
- Core muscles become active and support the spine when sitting, which helps them become stronger over time

Ergonomic saddle-seat design enhances posture

- Naturally aligns vertebrae and discs
- Relieves shoulder and neck tension

Engineered to support deeper breathing

- More oxygen promotes circulation and helps reduce fatigue
- Natural posture elongates the torso, more room for lung expansion
- Improved mobility



135DSS

All 135DSS Saddle Stools feature:

- Unique split-saddle seat design
- 23” dia. cast aluminum base w/ large casters
- 250lb weight capacity
- Standard 5-year warranty

	adjustable height range	shipping weight	seamless upholstery	recommended user height
135DSS	22.5” - 29.75”	28.4	Yes	5’3” - 5’11”
135DSS-L	18.75” - 23.75”	27.8	Yes	Under 5’3”
135DSS-H	24” - 33.75”	29	Yes	Over 5’11”



A unique angle on comfort.

The Brewer 135DSS is designed to provide the most beneficial seating position possible.

Simply adjust the seat height so your upper and lower legs create a 135° angle while both feet are flat on the ground. Tailbone should be almost flush with the back of the saddle. This position naturally aligns the spine to minimize strain on back muscles.

In studies, 87% of women and 80% of men have preferred the patented split-saddle-seat technology found in the innovative Brewer 135DSS to other saddle seating options.

Brewer

Manufactured with SALLI Technology Inside.

The Power To Advance. Safely. Financially. Competitively.

www.brewercompany.com | P: 1.888.273.9371

Brewer PS135

Polyurethane Saddle Stool

A unique angle on comfort

The Brewer PS135 is designed to provide the most comfortable seating position possible.

Simply adjust the seat height so your upper and lower legs create a 135° angle while both feet are flat on the ground. This position naturally straightens the spine to minimize strain on back muscles.



- Promotes better posture, relieving tension in the shoulders and neck
- Balances the pelvis in an upright position distributing the weight over the “sits-bones” and not the soft tissues
- Comfortable and durable soft black polyurethane seat
- Pneumatic height adjustment, three models available to accommodate differing user heights
- 21.5” diameter stable five-star black ABS plastic base
- Large dual-wheel casters, standard, for easy rolling
- Standard 5-year product warranty
- 250lb weight capacity
- Built in the USA, with US and imported parts, in accordance with US BIFMA quality standards



PS135

	adjustable height range	shipping weight (lbs)	recommended user height
--	-------------------------	-----------------------	-------------------------

PS135-1	18.5" - 23.5"	19	Under 5'3"
PS135-2	21.75" - 29"	19.8	5'3" - 5'11"
PS135-3	25" - 34"	20.6	Over 5'11"



Brewer

Manufactured with  SALLI Technology Inside.

The Power To Advance. Safely. Financially. Competitively.

www.brewercompany.com | P: 1.888.273.9371