

# Brewer

# Assist & AssistPRO

POWER PROCEDURE TABLE

ADVANCED PATIENT CARE FOR  
THE HEALTH OF YOUR PRACTICE.

- Preserve patient dignity and comfort during procedures, increasing your patient satisfaction rankings.
- Maximize throughput by having accessories and simple one-touch programming at clinician fingertips.
- ADA compliant to accommodate multiple patient populations.
- Reduce your total cost of ownership with our industry-leading **three year warranty** for peace of mind.



# Assist & AssistPRO

POWER PROCEDURE TABLE

Make procedures easier, safer and more comfortable for staff and patients.



## PROGRAMMABILITY

- AssistPRO - simple two-step programming process
- Memory for storing up to four preferred positions
- Visual and audible signals with auto-run for safety

## MORE FLEXIBILITY

- Supine, left lateral plane, lithotomy and proctology positions standard
- Patented 3-point pivot headrest system
- Extendable leg rest

## ADVANCED SAFETY

- 19" low height for easy wheelchair transfers
- 40" height and positioning for efficient procedures
- Pinch-free design and safety switches

## PATIENT COMFORT

- Seamless, multi-density, bolstered support
- ADA compliant to accommodate many patient populations

## SUPERIOR QUALITY

- Industry-leading **three year warranty**, standard
- UL 60601 certified
- **Made in USA**

### models

|       |  |
|-------|--|
| #7000 | Assist Power Procedure Table                 |
| #7500 | AssistPRO Programmable Power Procedure Table |

### standard features

|                                  |
|----------------------------------|
| Three-year warranty              |
| 450 pound weight capacity        |
| Height: 19" - 40"                |
| Length: 69", extends to 80"      |
| Width: 28"                       |
| Tilt range/Trendelenburg: to 30° |
| UL 60601-1 Certified             |

### accessories

|         |  |
|---------|--|
| #99497  | grab bars**                                    |
| #99499  | articulating armboard**                        |
| #99501  | side seat rails (pair)                         |
| #99502  | one universal socket clamp                     |
| #99503  | standard arm board**                           |
| #99504  | delux articulating knee crutches (pair) *      |
| #99505  | stainless steel treatment pan                  |
| #100253 | urology drain pan system                       |
| #102025 | Brewer Welch Allyn® IV bracket (plant install) |
| #102025 | Brewer Welch Allyn® IV bracket (field install) |

\* Needs one order of #99501

\*\* Needs one order of #99501 and two of #99502

