

Brewer

Access High-Low

EXAM TABLE

ADVANCED PATIENT CARE FOR
THE HEALTH OF YOUR PRACTICE.

- Increase your patient satisfaction rankings by preserving patient dignity and comfort during exams and procedures.
- Keep pace with the ever changing patient population requirements by accommodating a wider range of patients with a single table.
- More standard features. Pivoting safety grab bars, laptop work surface, and pass-through drawers equate to more storage and increased flexibility.
- Reduce your total cost of ownership with our industry-leading three year warranty for peace of mind.



Access High-Low

EXAM TABLE

Advance staff safety and care for more of your patients.

ABUNDANT STORAGE

- Patented pass-through side drawers eliminates left or right setup issues
- Pass-through laptop surface
- Convenient table paper storage under backrest



GRAB BAR TRANSFERS

- Stable and secure grab point for wheelchair and other transfers
- Conveniently swings out and are easily removable

ADVANCED SAFETY*

- 450 & 700 lb weight capacities
- Exclusive SafeGlide™ leg extension increases patient safety and satisfaction
- Adjustable, ergonomic stirrups integrated in table

SUPERIOR DESIGN

- Takes patient safely and automatically from consultation to examination
- 18" low height for easy and safe patient transfers
- Made in USA

ERGONOMIC UPHOLSTERY

- Seamless & easy to clean
- Multi-density, bolstered support upholstery
- 12 standard colors available

models

#6000	Access High-Low with pneumatic back
#6001	Access High-Low with pneumatic back, pelvic tilt, and front drawer warmer
#6500	Access High-Low with power back
#6501	Access High-Low with power back, pelvic tilt, front drawer warmer and outlet
#6800	Access High-Low 700 with power back
#6801	Access High-Low 700 with power back, pelvic tilt, front drawer warmer and outlet

accessories

#99464	stainless steel treatment pan
#99502	one universal socket clamp
#99504	deluxe articulating knee crutches (pair) **
#100253	urology drain pan system
#100407	front drawer dividers (can order up to 7)
#100540	two knee crutch support adapters (plant install)
#100550	two knee crutch support adapters (field install)
#101800	electrical outlet (6001 only, plant install)
#102031	Brewer Welch Allyn® IV bracket (plant install)
#102032	Brewer Welch Allyn® IV bracket (field install)

** Needs two #99502 and one of either #100540 or

*UL60601 certified



The Brewer Company LLC
N88W13901 Main Street, Ste 100
Menomonee Falls, WI 53051
1.888.Brewer.1
brewercompany.com

Brewer is an ISO 13485:2003 certified company.