

# Quality Products for Quality Care

# A single, convenient source for so many solutions.

Brewer understands that you need more than quality power tables to equip your exam and procedure rooms for success. That's why we offer a broad selection of ruggedly constructed clinical support solutions to complete your clinical care space. Each product is manufactured to Brewer's rigorous quality standards for design, materials and craftsmanship to provide the reliable, trouble-free service you need. Plus, all our clinical support products are backed by a standard 5-Year Warranty for added peace of mind.



# Hampers & Carts

# Durable designs you can depend on.

#### 18" Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Dimensions: 35.5" H x 18.25" outside diameter (90.2 cm x 46.3 cm).

#### Model 33330

#### Accessories

**31331** - 18" Rope style 10 oz. cotton bag **34000** - Bag holder

#### **Multi Purpose Cart**

Strong steel frame with durable white powder coat finish. Total weight capacity: 50 lbs (22.7 kg). Dimensions with handle:

37.5" H x 19.5" W x 29" D (95.3 cm x 49.5 cm x 73.6 cm). Baskets sold separately.

#### Models

63400 - Cart frame only

63021 - 12" White powder coat wire basket 63020-1 - 6" White powder coat wire basket (1 per carton)

**63020-2** - 6" white powder coat wire basket (2 per carton)

#### **Plastic Utility Cart**

Incudes three ABS plastic shelves,  $16'' \text{ W} \times 21.75'' \text{ L}$  (40.6 cm x 55.3 cm). Total weight capacity: 200 lbs (90.7 kg). Dimensions with handle:  $28'' \text{ L} \times 17.5'' \text{ W} \times 28'' \text{ H}$  (71.1 cm x 44.5 cm x 71.1 cm).

Model 63530

## 18" Tilt-Top Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Step-on pedal for hands-free use. Dimensions: 35.75" H x 18" outside diameter

 $(90.6 \text{ cm} \times 45.7 \text{ cm})$ 

#### Model 33334

#### Accessories

31331 - 18" Rope style 10 oz. cotton bag 34000 - Bag holder

## **Adjustable Folding Linen Hamper**

Holds 18'', 22'' or 25'' diameter bags (45.7 cm, 55.9 cm or 63.5 cm). Features durable white powder-coat finish. Folds for compact storage. Dimensions:  $37'' \text{ H} \times 20'' \times 20'' \text{ D}$   $(93.9 \text{ cm} \times 50.8 \text{ cm}) \times 50.8 \text{ cm}$ .

#### Model 33395

#### Accessories

31331 - 18" Rope style 10 oz. cotton bag 31341 - 25" Rope style 10 oz. cotton bag

# Tilt-Top Square Hamper

Open design for easy unloading. Step-on pedal for hands-free use. White powder-coat finish. Built-in bag holder supports heavy- duty bags. Dimensions:

 $38" \text{ H} \times 18.5" \text{ W} \times 19.5" \text{ D } (96.5 \text{ cm} \times 47 \text{ cm} \times 49.5 \text{ cm}).$ 

#### Model 11410

#### Accessories

31331 - 18" Rope style 10 oz. cotton bag

# Mobile Equipment Cart

Powder coated steel cabinet with integral handles, locking casters, & rubber bumpers to prevent wall damage. Has three small drawers &

one large drawer.

Dimensions: 34.4" H
x 24.2" W x
21.6" D (87.4 cm x
61.5 cm x 54.9 cm).

#### Model 63100

#### Accessories

100407 - Drawer dividers

# Stainless Steel All Purpose Cart

Features non-corrosive stainless steel construction & three 15.5" W x 24" L (39.4 cm x 61 cm) shelves. Total weight capacity: 300 lbs (136 kg). Dimensions with handle: 28" L x 16.5" W x 32" H (71.1 cm x 41.9

Model 63500

cm x 81.3 cm)



Visit brewercompany.com

2

# Step Stools

# Designed & built to support confident care.

#### Step Stool

Features rugged, steel construction with durable, black powdercoat finish and a skid-resistant rubber mat for safety. Dimensions: 9" H x 11.25" W x 14.25 L  $(22.8 \text{ cm} \times 28.6 \text{ cm} \times$ 36.2 cm). With handrail, total height is 35.25" (89.5 cm). Weight capacity: 350 lbs (158 kg).

#### Models

11200-1 - 1 Stool per box 11200-2 - 2 Stools per box 11220 - w/Hand rail

#### Accessories

99121 - Replacement crutch tips (set of 4)

#### Two-Step, Step Stool

Features heavy-duty, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety. Dimensions (1st step): 9" H x 12" W x 16 L (22.9 cm x

**Bariatric Step Stool** 

powder-coat finish

and a skid-resistant

rubber mat for safety.

Dimensions: 9.25" H x

11.25" W x 14.25 L

 $(23.5 \text{ cm} \times 28.6 \text{ cm} \times$ 

total height is

34" (86.4 cm).

Models

Accessories

Weight capacity:

600 lbs (272 kg).

21200 - Step stool

(set of 4)

21220 - Step stool w/hand rail

99121 - Replacement crutch tips

36.2 cm). With handrail,

Features heavy-duty, steel

construction with durable, black

 $30.5 \text{ cm} \times 40.6 \text{ cm}$ ). Dimensions (2nd step): 7" H x 12" W x

16 L (17.8 cm x  $30.5 \text{ cm} \times 40.6 \text{ cm}$ ). With handrail, total height is 44"

(111.8 cm).

Weight capacity: 600 lbs (272 kg).

Models **31200** - Step stool

31220 - Step stool w/hand rail

#### Accessories

50016-1 - Hand rail

# Instrument Stands, Tables & Kick Bucket

#### **Mobile Instrument Stand**

Adjustable height from 29"-47.5" (73.6 cm-120.7 cm) for comfortable use sitting or standing. Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm  $\times$  32.3 cm  $\times$  1.9 cm).

#### Model 43465

#### Accessories

2102568 - Replacement tray

#### **Mayo Instrument Stand**

U-shaped base minimizes interference during procedures. Adjustable height from 31"-48.5" (78.7 cm-123.2 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75"  $(48.8 \text{ cm} \times 32.3 \text{ cm})$  $\times 1.9 \text{ cm}$ ).

#### Model 43460

#### Accessories

2102568 - Replacement tray

#### **Mayo Instrument Stand**

"California" style base for easy positioning. Adjustable height from 34"-52" (86.4 cm-132.1 cm). Removable stainless steel tray:  $19.2'' \times 12.7'' \times 0.75''$  $(48.8 \text{ cm} \times 32.3 \text{ cm} \times 1.9 \text{ cm}).$ 

#### Model 11100

#### Accessories

2102568 - Replacement tray

## Stainless Steel Mayo Stand

Features non-corrosive stainless steel construction. Adjustable height from 31"-48.5" (78.7 cm-123.2 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75"  $(48.8 \text{ cm} \times 32.3 \text{ cm} \times 1.9 \text{ cm}).$ 

#### Model 48461

#### Accessories

2102568 - Replacement tray

#### "U" Base Table

"U" base fits easily in restricted spaces. Offers convenient, spring-loaded height adjustment from 28.5"-43.25" (72.4 cm-109.9 cm). Top measures  $15'' \times 30''$  (38.1 cm x 76.2 cm).

#### Models

11620 - w/Walnut grain laminate top & chrome base 11630 - w/Gray laminate top &

gray epoxy finish on base

#### Kick Bucket w/Frame

Features non-corrosive stainless steel construction. Heavy-duty design includes 360° rubber bumper on frame 13 quart size (3.25 gallon).

#### Model 48711

48712 - Kick bucket only 48714 - Kick bucket frame only

#### Low "H" Base Table

Designed for low care beds. Features spring-loaded height adjustment. With casters: 19.375"-26.5" (49.2 cm-67.3 cm). With glides: 17.625"-24.5" (44.8 cm-62.2 cm). Top measures 15" x

30" (38.1 cm x 76.2 cm)

#### Model

11640 - w/Walnut grain laminate top & chrome base

#### "H" Base Table

Offers convenient, spring-loaded height adjustment from 28.5"-43.25" (72.4 cm-109.9 cm). Top measures 15" × 30"

(38.1 cm x 76.2 cm)

#### Models

11600 - w/Gray laminate top & gray epoxy finish on base

11610 - w/Walnut grain laminate top & chrome base



# Pump Stands, IV Poles & Waste Cans

# Quality that stands the test of time.

#### **Large Infusion Pump Stand**

Our 25.5" diameter (64.8 cm) design delivers the most stable base available, handling up to four infusion pumps and/ or IV controllers (w/ accessories 43433 & 43430). Exclusive "No Loss" knob allows easy height adjustment from 55"-96.5" (139.7 cm-245.1 cm).

#### Models

43408 - 2 Ram's horn hook 43409 - 4 Ram's horn hook

## **Heavy Base Four-Leg IV Pole**

Features tough, four-leg epoxy-coated steel base with 22" (55.9 cm) diameter, wall-protecting rubber bumpers on the legs, and exclusive "No Loss" knob for easy height adjustment from 51.5"-93" (130.8 cm-236.2 cm).

#### Models

**43403** - 2 Ram's horn hook 43406 - 4 Ram's

horn hook

## **Short Wheel Infusion Pump Stand**

Stable, space-saving design features heavy center weight with only a 17.5" diameter (44.5 cm) and exclusive "No Loss" knob for easy height adjustment from 54"-96.5" (137.2 cm-245.1 cm).

#### Models

horn hook **43417** - 4 Ram's horn hook

43416 - 2 Ram's

#### **Aluminum Base Five-Leg IV Pole**

Includes five-leg cast aluminum base with 22.25" (56.5 cm) diameter, and exclusive "No Loss" knob for easy height adjustment from 52"-93" (132 cm-236.2 cm).

## Models

11350 - 2 Ram's horn hook 11360 - 4 Ram's horn hook

# Two-Hook Four-Leg IV Pole

Offers four-leg, epoxycoated steel base with 22" (55.9 cm) diameter and easy, twist-lock height adjustment from 51"-92" (129.5 cm-233.6 cm).

## Model

11300 - 2 Ram's horn hook

#### Accessories for 43416, 43417, 43408, 43409

The following accessories all attach with the #43433 clamp, which is sold separately.

## **Universal Clamp**

Offers fast, flexible attachment & positioning of accessories with set-screw locks.

Model 43433

# Oxygen Cylinder Holder

Holds "D" & "E" size tanks. Ideal for ambulatory patients.

Model 43429

## **Infusion Pump Frame**

Mount up to four pumps on a single stand. Easily mount, adjust and remove pump rack. Only for use with infusion pump stands #43408 and #43409.

#### Model 43430

\*Two 43433 clamps needed.

# Add-A-Tray

Removable stainless steel tray measures 9.75" x 13.5" x 0.75"  $(48.9 \text{ cm} \times 32.1 \text{ cm})$  $\times$  1.9 cm).

Model 43431



#### Large Round Steel Waste Can 32 Quart (8 Gallon) Models

Model 45269

45266 - White enamel finish

45267 - Red enamel finish **45268** - Beige

enamel finish **45269** - Stainless steel finish

# **Square Steel Waste Can**

32 Quart (8 Gallon) Models

35266 - White enamel finish 35267 - Red enamel finish

35268 - Beige enamel finish

Model 35271

#### 20 Quart (5 Gallon) Models

35269 - White enamel finish 35270 - Red enamel finish 35271 - Beige enamel finish

# **Rectangular Plastic Waste Can**

## 52 Quart (13 Gallon) Models

25266 - White plastic **25267** - Red plastic

25268 - Beige plastic



Model 25267

### 32 Quart (8 Gallon) Models

25269 - White plastic **25270** - Red plastic

25271 - Beige plastic

Brewer

Visit brewercompany.com

# Task Lighting

# Brilliant additions to any clinical care setting.

Brewer offers a broad range of solutions to meet your room lighting needs. All feature flexible, "no-drift," gooseneck arms for precise positioning.

#### **LED ECO-Series**

Combines premium performance with a low cost of ownership.

This eco-friendly LED light consumes significantly less energy than halogens and its 50,000-hour bulb life is twice that of other LED task lights available. Efficient power consumption also means it remains cooler to the touch.

Yet it provides an industry-best, bright white light (more than 20,000 lux / 880 Lumens with a 5''diameter beam at 18") to drive clinical efficiency. Color temperature is 4000°-6000° Kelvin.

#### Models 19100

200165 - Optional five-caster base 2102755 - Replacement bulb (50,000-hour bulb life)

#### **Incandescent Series**

Advanced reflector system provides shadow-free, crisp white light. Includes a 25-Watt incandescent bulb.

Features a non-metallic head to eliminate surface heat.

#### Model 17100

200165 - Optional five-caster base 001525 - 25-Watt replacement bulb (500-hour bulb life)

#### Halogen 35 Series

Delivers crisp, white light (1,500 Footcandles at 24" [61 cm]). High-intensity, low-voltage 35-Watt bulb included (Color temperature is 3000° Kelvin).

Features a non-metallic head to eliminate surface heat; dichroic film-coated reflectors filter out 80% of beam heat; and heat sink technology removes bulb socket heat to help keep the housing cooler to the touch.

#### Model 18100

200165 - Optional five-caster base 001635 - 35-Watt replacement bulb (4,000-hour bulb life)



#### **Chrome Series**

Accepts 10 to 40-Watt incandescent bulb (not included).

Features a durable chrome finish throughout.

Model 11500





The Brewer Company, LLC N88 W13901 Main Street, Suite 100 Menomonee Falls, WI 53051

1.888.273.9371 www.brewercompany.com

