




Brewer

INFUSION PUMP STANDS

INFUSION PUMP STANDS

SERIES Strength. Stability. Versatility.



Large Infusion Pump Stand

- The most stable base available (26" diameter)
- Handles up to four infusion pumps and/or IV controllers (w/43430 and 43433 accessories)
- Tough six-leg epoxy coated steel base
- Chrome plated top post assembly
- 3" rubber wheel, ball bearing casters
- Exclusive "No Loss" knob
- Height: 52.25" - 102.75"
- Load capacity: 45 lbs.

#43408
2 ram's horn hook

#43409
4 ram's horn hook



Short Wheel Infusion Pump Stand


- Heavy center weight affords stability while saving space (16" diameter)
- Tough five-leg epoxy coated steel base
- Chrome plated top post assembly
- 3" rubber wheel, ball bearing casters
- Exclusive "No Loss" knob
- Height: 54" - 90.5"
- Load capacity: 45 lbs.

#43416
2 ram's horn hook

#43417
4 ram's horn hook

ACCESSORIES


The following accessories all attach with the #43433 clamp, which is sold separately



Universal Clamp

- Fast, flexible attachment & positioning of accessories
- Set-screw locks accessory in place
- Secure grip; versatile adjustment


#43433



Infusion Pump Frame

- Mount up to four pumps on a single stand
- Easy mounting, removal & adjustment of pump rack
- Only suitable for use with #43408 and #43409 infusion pump stands

#43430



Steering Handle

- Allows easy & convenient mobility
- Non-slip black vinyl coating


#43428



Oxygen Cylinder Holder

- Holds "D" and "E" size tanks
- Excellent for ambulatory patients

#43429



Add-A-Tray

- Holds instruments or patient items
- Removable stainless steel tray
- Tray size: 9.75" x 13.5" x .75"

#43431

All Brewer infusion pump products come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
N88 W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Brewer