




Brewer

IV POLES

# IV POLES

SERIES Stable. Versatile. Reliable.




## Aluminum Base Five-Leg

- Five-leg, 25" cast aluminum base
- Chrome plated steel tubing
- 2" Dual wheel nylon casters
- Exclusive "No Loss" knob
- Top pole: .75" diameter
- Bottom pole: 1" diameter
- Height: 52.25" - 92.75"
- Load capacity: 45 lbs.
- Some assembly required

#11350  
2 ram's horn hook

#11360  
4 ram's horn hook




## Heavy Base Four-Leg

- Tough four-leg, epoxy coated 22" steel base
- Chrome plated steel tubing
- 3" Rubber wheel, ball bearing casters
- Exclusive "No Loss" knob
- Rubber bumpers on legs to protect walls
- Top pole: .75" diameter
- Bottom pole: 1" diameter
- Height: 51.75" - 93.5"
- Load capacity: 45 lbs.
- Some assembly required

#43403  
2 ram's horn hook

#43406  
4 ram's horn hook



## Two-Hook Four-Leg

- Four-leg, epoxy coated 22" steel base
- Chrome plated steel tubing
- Dual wheel nylon casters
- Easy twist lock height adjustment
- Height: 51.25" - 91.75"
- Top pole: .75" diameter
- Bottom pole: 1" diameter
- Load capacity: 45 lbs.
- Some assembly required

#11300  
2 ram's horn hook

### All Brewer IV poles come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

**Call your local Brewer sales representative to learn about our complete product line.**

The Brewer Company  
N88 W13901 Main Street, Menomonee Falls, WI 53051  
1.888.Brewer.1 [www.brewercompany.com](http://www.brewercompany.com)

# Brewer