



Brewer

UTILITY CARTS

UTILITY CARTS

SERIES

Sturdy. Mobile. Durable.



Multi Purpose Cart

- Strong chrome plated steel frame
- 3" rubber wheel, ball bearing casters
- Convenient height of 37.5"
- Coated steel basket supports
- Total capacity: 50 lbs.
- No assembly required
- Dimensions: 37.5" H x 19.5" W x 30" D
- Large Basket: 12.0" H x 14.5" W x 24.5" D
- Small Basket: 6.0" H x 14.5" W x 24.5" D
- Baskets sold separately

#63400	cart frame only
#63021	12" wire basket
#63020-1	6" wire basket (1 per carton)
#63020-2	6" wire basket (2 per carton)



Plastic Utility Cart

- 3 ABS plastic shelves (21.75" x 16")
- 12" between shelves
- Four, 3" non-locking, swivel polyurethane ball-bearing casters
- Total capacity: 200 lbs.
- Easy to assemble
- Color: light grey
- Dimensions with handle:
- 25.25" L x 17.5" W x 35" H

#63530 assembly required



Mobile Equipment Cart

- Great for the Primary Care Market, used in exam rooms
- Powder coated steel cabinet
- Four black, locking casters
- Rubber bumpers prevent wall damage
- Contoured top covered in a bacteriological, chemical, and fungal resistant material
- Two handles integrated directly into the cart
- Cart size: 34.4" H x 24.2" W x 21.6" D
- 3 small drawers: 4.0" H x 19.7" W x 19.2" D
- 1 large drawer: 10.0" H x 19.7" W x 19.2" D
- Inside dimension of top: 23.0" W x 20.375" D
- Top weight capacity: 100lbs
- Drawer weight capacity: 15lbs

#63100 assembled

accessories

#100407 drawer dividers



Stainless Steel All Purpose Cart

- Non-corrosive stainless steel construction
- 3 stainless steel shelves (15.5"W x 24"L)
- Roomy 12" between shelves
- Quietly rides on 3.5" rubber wheel casters
- Total weight capacity: 300 lbs.
- Dimensions with handle: 28" L x 15.5" W x 33" H

#63500 assembled

All Brewer carts come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
N88 W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Brewer

07.2012
00009.10