



Brewer

INSTRUMENT STANDS

# INSTRUMENT STANDS

## SERIES

Ideal. Portable. Accessible.




### Mayo Instrument Stand

- "California" style base for easy positioning
- 1" chrome plated steel tubing
- Two rubber wheel, ball bearing casters
- Removable stainless steel tray:
- 19.25" x 12.62" x .75"
- Base: 20" wheels, 19" diameter
- Height: 34" - 53"
- Weight capacity: 5lbs

#11100

#### accessories

#09412 replacement tray,  
19.25" x 12.62" x .75"



### Mobile Instrument Stand

- Adjustable height for comfortable use while sitting or standing
- 1" chrome plated steel tubing
- Stable 25" five-leg cast aluminum base
- Removable stainless steel tray:
- 19.25" x 12.62" x .75"
- 2" Dual wheel nylon casters
- Height: 29.25" - 48.75"
- Weight capacity: 5lbs

#43465

#### accessories

#09412 replacement tray,  
19.25" x 12.62" x .75"



### Electrosurgical Unit Stand

- Sturdy and stable for easy room-to-room portability and storage
- Adapts easily to most electrosurgical units
- Angled mounting panel for easy viewing
- Includes convenient cord and treadle hooks
- Stable five-leg cast aluminum base
- Dual wheel nylon casters
- Base: 22"
- Height: 37.75"

#43577

#### accessories

#43431 optional add-a-tray  
9.75" x 13.5" x .75"

#43433 optional clamp



### Mayo Instrument Stand

- Designed for convenience during surgical procedures
- 1" chrome plated steel tubing
- U-shaped base minimizes interference with other equipment or patient
- Two rubber wheel, ball bearing casters
- Removable stainless steel tray: 19.25" x 12.62" x .75"
- Base: 19" wheels, 19.5" diameter
- Height: 32" - 50"
- Weight capacity: 5lbs

#43460

#### accessories

#09412 replacement tray,  
19.25" x 12.62" x .75"

### All Brewer instrument stands come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

**Call your local Brewer sales representative to learn about our complete product line.**

The Brewer Company  
N88 W13901 Main Street, Menomonee falls, WI 53051  
1.888.Brewer.1 [www.brewercompany.com](http://www.brewercompany.com)

# Brewer