



Brewer

OVERBED TABLES / STAINLESS STEEL

# OVERBED TABLES

Versatile. Convenient. Classic.

# STAINLESS STEEL

Easy to clean with exceptional durability



## "H" Base Tables

- Easy spring loaded height adjustment
- Four point H-base with two locking dual wheel nylon casters
- Height: 26.75" - 43.25"
- Base: 15.25" x 26.38"; Top: 15" x 30"
- Weight capacity: 50 lbs.
- Some assembly required

#11600 "H" base  
gray laminate top, gray epoxy finish on base

#11610 "H" base  
walnut grain laminate top, chrome base



## "U" Base Tables

- "U" base to fit easily in restricted space areas
- Easy spring loaded height adjustment
- Four point U-base with two locking dual wheel nylon casters
- Height: 26.75" - 43.25"
- Base: 15.25" x 26.38"; Top: 15" x 30"
- Weight capacity: 50 lbs.
- Some assembly required

#11620 "U" base  
walnut grain laminate top, chrome base

#11630 "U" base  
gray laminate top, gray epoxy finish on base



## Low "H" Base Table

- Convenient table solution for low care beds
- Easy spring loaded height adjustment
- Four point H-base with two locking dual wheel nylon casters
- 1.75" casters and plastic glides included
- Height: Casters 19.375" - 26.5"; Glides 17.625" - 24.5"
- Base: 15.25" x 26.38"; Top: 15" x 30"
- Ground clearance: .5" with glides; 2.25" with casters
- Weight capacity: 50 lbs.
- Some assembly required

#11640 low "H" base  
walnut grain laminate top, chrome base



## Stainless Steel Mayo Stand

- Non-corrosive stainless steel construction
- Removable stainless steel tray (19.25" x 12.63" x .75")
- Height: 32" - 50"
- Rubber wheel, ball bearing casters

#48461

accessories

#09412 replacement tray



## Stainless Steel All Purpose Cart

- Non-corrosive stainless steel construction
- 3 stainless steel shelves (15.5"W x 24"L)
- Roomy 12" between shelves
- Quietly rides on 3.5" rubber wheel casters
- Total weight capacity: 300 lbs.
- Dimensions with handle: 28" L x 15.5" W x 33" H

#63500 assembled



## Kick Bucket w/Frame

- Non-corrosive stainless steel construction
- Heavy duty, exceptionally durable
- 360° rubber bumper on frame
- 13 quart size, 12" diameter, 10" deep
- 2" rubber casters

#48711

#48712  
kick bucket only

#48714  
kick bucket frame only

**All Brewer overbed tables and stainless steel products come with a 5-year warranty.**

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

**Call your local Brewer sales representative to learn about our complete product line.**

The Brewer Company  
N88 W13901 Main Street, Menomonee falls, WI 53051  
1.888.Brewer.1 www.brewercompany.com

**Brewer**