



**Which Control Plate and Backrest Upright Bar Do You Have?
(for all Medical Stools except the foot operated Millennium stools)**

Straight Tube



End View




Angle View



Top View

Old Style Control Plate 101598, 2008 thru 13-Feb-2013


Double Wall
Formed Tube
(Dome Shaped)



End View




Angle View



Top View


New Style Control Plate 2102167, from 13-Feb-2013 on



100956

100957

Old Style Backrest Upright Bar, 2008 thru 13-Feb-2013



2102168

2102169

New Style Upright Bar, from 13-Feb-2013 on