

BrewerAssistPRO

POWER PROCEDURE TABLE

BrewerAssist

POWER PROCEDURE TABLE



Advanced patient care for the health of your practice.

- 450 lb weight capacity
- 3-year standard warranty
- Safe, extra-low voltage actuators
- Infinitely adjustable to meet any procedure needs



Explore the Assist Power Table family here. →



Supine



Left Lateral Plane



Lithotomy




Proctology



Easy wheelchair patient transfer at 19"



Convenient table paper storage



Supports up to 450 pounds



Comfortable, ergonomic upholstery; easy-to-clean surface



Exclusive, easy-to-adjust three-point pivot headrest system



Footrest extends an additional 5"



Standard intuitive foot controls



Handheld controls standard on AssistPRO

Brewer

Furnishing solutions for safe practices.

AssistPRO Programmability

- Intuitive hand pendant and foot control included
- Four preferred positions can be programmed into the table's memory
- Simple two-step programming process
- Visual and audible signals create a user-friendly interface
- Auto-run toggle for safety

Advanced Safety

- Safety release actuators
- Safe extra-low voltage actuators
- 19" (48.3 cm) low height
- 40" (101.6 cm) maximum height
- Leg rest extends an additional 5" (12.7 cm)
- Adjustable three-point pivot headrest
- Easy-to-clean seamless upholstery
- 18" (45.7 cm) or 21" (53.3 cm) paper dispenser

Patient Comfort

- Easy lateral wheelchair transfers
- Infinitely adjustable headrest
- Extendable leg rest
- Standard adjustable ergonomic stirrups
- Multi-density, bolstered support upholstery

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

- UL 60601-1



Table Specifications

- 450 lb (204 kg) weight capacity
- Standard 3-year warranty
- Standard adjustable stirrups
- Standard foot control
- Hand pendant standard on AssistPRO (7500 models)
- Height: 19" - 40" (48.3 cm - 101.6 cm)
- Length: 69" (75.3 cm) when flat
- Length with headrest & leg rest extended: 80" (203.2 cm)
- Length in upright position: 40" (101.6 cm)
- Width: 28" (71.1 cm)
- Leg rest: 16.6" x 11.25" (42.2 cm x 27.9 cm)
- Headrest: 13" x 11" (33 cm x 27.9 cm)
- Tilt range / Trendelenburg: to 30°
- Base: 24" W x 36" L (61 cm x 91.4 cm)
- Power cord: 70" (177.8 cm)
- 115VAC, 60 Hz
- 15 standard upholstery color options
- Cal-133 fire code compliant upholstery available in 4 colors
- Premium, plush upholstery available in 5 colors
- Custom colors available

models

#7500	Brewer AssistPRO Power Procedure Table
#7000	Brewer Assist Power Procedure Table

accessories

#99497*	grab bars
#99499*	articulating armboard
#99500	hand pendant (standard on Brewer AssistPRO)
#99501	side seat rails (qty. 2)
#99502	one universal accessory clamp
#99503*	standard arm board
#99504**	deluxe articulating knee crutches
#99505	stainless steel treatment pan
#100253	urology drain pan system
#102024	Brewer Welch Allyn® IV bracket (plant installed)
#102025	Brewer Welch Allyn® IV bracket (field installed)

* Needs one order of #99501

** Needs one order of #99501 and two of #99502

The Brewer Company is an ISO 13485:2003 certified company.

All Brewer Assist™ and AssistPRO™ Power Procedure Tables come with a standard 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
 N88W13901 Main Street, Menomonee Falls, WI 53051
 1.888.Brewer.1 | www.brewercompany.com