

BrewerBasic

EXAM TABLE




Abundant Storage Capacity

- Largest storage capacity available
- Front Drawers: 19.7" x 19.2" x 4"
- Side Drawers: 22.7" x 20.2" x 7"
- Side drawers available on right or left side (patient's perspective)



Sturdy Step

- Sturdy, retractable step
- Slip resistant mat for patient safety
- 20.5" x 11" x 8.6"
- Tucks away easily



Comfortable, Ergonomic Upholstery

- Seamless, medicinal stain resistant
- Available in 15 standard colors
- Cal-133 fire code upholstery available
- Simple pneumatic lift backrest



Ergonomic Stirrups

- Extendable, adjustable
- Conveniently stored in the table
- Manually extended through infinite positions and four lateral positions



Large Drawer Volume

- Front drawer: 1,512.96 in³ each
- Side drawer: 3,209.78 in³ each
- Total table volume: 9,445.48 in³

Brewer

Furnishing solutions for safe practices.

BrewerBasic

EXAM TABLE

Advanced Safety

- Sturdy, retractable patient step with mat
- Simple pneumatic back for easy patient positioning

Patient Comfort

- Adjustable ergonomic stirrups
- Leg extension / plastic debris tray
- Multi-density, bolstered support upholstery

Abundant Storage

- Largest storage capacity available
- Side drawers available on right or left side (from patient's perspective)
- Front drawer: 19.8" x 19.2" x 3.8"
(50.3 cm x 48.8 cm x 9.7 cm)
- Side drawers: 22.6" x 20.2" x 7"
(57.4 cm x 51.3 cm x 17.8 cm)

Table Dimensions

- 325 lb (147 kg) weight capacity
- Standard 3-year warranty
- Height: 32" (81.3 cm)
- Length: 57.5" extends to 71" (146.1 cm - 180.3 cm)
- Width: 26.5" (67.3 cm)
- Leg extension length: 15" (38.1 cm)
- Step dimensions: 20.5" x 11" x 8.6"
(52.1 cm x 27.9 cm x 21.8 cm)
- Optional 5° pelvic tilt / 60-Watt front drawer heater
- Easy-to-clean seamless upholstery
- 15 standard upholstery color options
- Cal-133 fire code compliant upholstery available
- Custom colors also available

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

- UL 60601-1



4001-L in Azure Blue

models

- | | |
|-------|---|
| #4000 | Basic Exam Table with pneumatic back |
| #4001 | Basic Exam Table with pneumatic back, pelvic tilt and drawer heater |

accessories

- | | |
|---------|--|
| #98464 | stainless steel treatment pan |
| #100253 | urology drain pan system |
| #100407 | front drawer dividers (can order up to 7) |
| #102022 | Brewer Welch Allyn® IV bracket (plant installed) |
| #102023 | Brewer Welch Allyn® IV bracket (field installed) |

The Brewer Company is an ISO 13485:2003 certified company.

All Brewer Basic Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company LLC
N88W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 | www.brewercompany.com

Brewer