

# Brewer

# Access

## EXAM TABLE

### ADVANCED PATIENT CARE FOR THE HEALTH OF YOUR PRACTICE.

- Reduce the risk of profit-robbing patient and staff injury with large step and easy accessibility.
- Preserve exam room space and save clinicians' time with up to **30% more drawer area** to maximize tight exam room spaces.
- Improve your patient's satisfaction and help infection control efforts with Access' superior upholstery.
- Reduce your total cost of ownership with our industry-leading **five year warranty** for peace of mind.



# Access

## EXAM TABLE

Drive safety, efficiency and patient satisfaction in the examination room.



### MAXIMUM STORAGE

- The most standard storage space - **up to 30% more drawer volume**
- Convenient table paper storage under backrest

### ERGONOMICS & COMFORT

- Seamless, multi-layer upholstery standard
- Simple pneumatic backrest
- Adjustable, integrated ergonomic stirrups

### ADVANCED VERSATILITY

- Patented pass-through side drawers allows easier access to equipment from the left or right side of the exam table

### SUPERIOR QUALITY

- Industry-leading standard **five year warranty** included
- UL 60601 certified
- **Made in the USA**

### SAFEST ACCESSIBILITY

- **Largest step surface, 17.7" x 20.4"**, reduces the risk of patient falls
- **Secure, four leg design** provides stability & support

#### models

#5000	Access Exam Table with pneumatic back
#5001	Access Exam Table with pneumatic back, pelvic tilt and drawer heater

#### standard features

Five-year warranty
500 pound weight capacity
Large, stable, four-leg patient step
Patented pass-through side drawers
Hospital-grade duplex electrical outlets
Simple, pneumatic lift backrest
Covered extra paper roll storage
Dimensions: 32" high x 28" wide
Length: 57.3" extends to 71"
UL 60601-1 Certified

#### accessories

#98464	stainless steel treatment pan
#100253	urology drain pan system
#100407	front drawer dividers (can order up to 7)
#102022	Brewer Welch Allyn® IV bracket (plant install)
#102023	Brewer Welch Allyn® IV bracket (field install)

