

# Brewer

## Design Guide



For  
Architects,  
Designers,  
Builders,  
Planners

The Power to Advance

Founded by E. F. Brewer in 1947, the Brewer Company leads the way in the development of innovative solution through our comprehensive line of power procedure and exam tables as well as the industry's first mobile patient lift, the Brewer LiftMate™.

Our process-oriented solutions are designed to help you safely meet the patient handling needs of individuals with mobility issues as well as improve operational efficiencies while more effectively protecting patients and staff.

Whether it's an aging population, obese patients, pregnant women or those with an injury or physical disability, Brewer provides the practical, common-sense approaches to today's do-more-with-less demands that you won't find anywhere else.

What's more, our tables and patient lift products are proudly made in the USA. Our extensive line of solutions is manufactured to Brewer's rigorous quality standards and backed by best-in-class standard warranties. It's all part of our commitment to provide you with The Power to Advance. Safely. Financially. Competitively.



## Table of Contents

INTRODUCTION	2
PRODUCT OVERVIEW/FLOOR PLANS	
Assist Power Procedure Table	4
Access High-Low Power Exam Table	8
Access Examination Table	12
Basic Examination Table	16
LiftMate Low/High Mobile Patient Lift	20
OTHER BREWER PRODUCT LINES	
Task Lighting Series	24
Exam Stools and Medical Seating	26
Additional Support Products	28
BREWER COLOR OPTIONS	
	30



# BrewerAssist

POWER PROCEDURE TABLE

## Advanced Safety

- Safety release actuators
- Safe extra-low voltage actuators
- 19" (48.3 cm) low height
- 40" (101.6 cm) maximum height
- Leg rest extends an additional 5" (12.7 cm)
- Adjustable three-point pivot headrest
- Easy-to-clean seamless upholstery
- 18" (45.7 cm) or 21" (53.3 cm) paper dispenser

## Table Dimensions

- Height: 19" – 40" (48.3 cm – 101.6 cm)
- Length: 69" (75.3 cm)
- Length with headrest & leg rest extended: 80" (203.2 cm)
- Length in upright position: 40" (101.6 cm)
- Width: 28" (71.1 cm)
- Leg rest: 16.6" x 11.25" (42.2 cm x 27.9 cm)
- Headrest: 13" x 11" (33 cm x 27.9 cm)
- Tilt range / Trendelenburg: to 30°
- Base: 24" W x 36" L (61 cm x 91.4 cm)
- Power cord: 70" (177.8 cm)

## Other Specifications

- 450 lb (204 kg) weight capacity
- Standard three-year warranty
- Standard adjustable stirrups
- Standard foot control
- Hand pendant standard on 7500 models
- 115 VAC, 60 Hz
- 12 standard upholstery color options
- CAL-133 fire resistant upholstery available in 4 colors
- Premium, plush upholstery available in 5 colors
- Custom colors available



### models

#7000	Assist Power Procedure Table
#7500	AssistPRO Programmable Power Procedure Table

Make procedures easier, safer and more comfortable for staff and patients.

Brewer's Power Procedure Tables are infinitely adjustable to meet any procedure needs.



Supine



Left Lateral Plane




Lithotomy



Proctology

## Floor Plans & Work Flow



Room  
14 ft. x 12 ft.

# BrewerAccess High-Low EXAM TABLE

## Advanced Safety

- Safety grab bars for patient stability
- Return-to-chair feature
- On/off switch prevents accidental movement
- Sensors prevent injuries during table movement
- 18" low height, 28" when flat (45.7cm,71.1cm)
- Easy on/off for patients
- Powered backrest adjusts for easy positioning

## Abundant Storage

- Pass-through work surface on side of table
- Patented pass-through side drawers accessible from right or left side of the table
- Front drawer: 19.8" x 19.2" x 3.8" (50.5 cm x 48.9 cm x 9.8 cm)
- Side drawers: 22.6" x 20.2" x 7" (57.5 cm x 51.4 cm x 17.8 cm)
- Covered storage for up to 6 paper rolls
- Paper storage capacity: 23"x 23"x 3.25" (58 cm x 58 cm x 8 cm)
- Slide through holder accepts 21" & 18" rolls (53.3 cm & 45.7 cm)

## Table Dimensions

- Height: 18" - 37" (45.7 cm - 94 cm)
- Length: 57" extends to 78" (144.8cm-198.1cm)
- Upholstery width: 28" (71.1 cm)
- Leg extension length: 21" (53.3 cm)
- Goes flat at 28" (71.1 cm)

## Other Specifications

- 450 lb (204 kg) & 700 lb (317 kg) weight capacities
- Standard patient assist grab bars
- Standard adjustable stirrups
- Standard foot control
- Standard three-year warranty
- 115VAC, 60 Hz, 5.8 amps max
- Safe, low voltage actuators
- Safe, single power cord design, detachable
- Easy to clean, seamless upholstery
- 15 standard upholstery color options
- CAL-133 fire resistant upholstery available in 4 colors
- Premium, plush upholstery available in 5 colors
- Custom colors available

## 6801 model only:

- Hospital grade outlet standard
- 115VAC, 60 Hz, 9.3 amps max
- 5° pelvic tilt, 60-Watt drawer heater in front top drawer

## models


#6000	Access High-Low with pneumatic back
#6001	Access High-Low with pneumatic back, pelvic tilt, and front drawer warmer
#6500	Access High-Low with power back
#6501	Access High-Low with power back, pelvic tilt, front drawer warmer and outlet
#6800	Access High-Low 700 with power back
#6801	Access High-Low 700 with power back, pelvic tilt, front drawer warmer and outlet



Advance staff safety and care for more of your patients.



## Floor Plans & Work Flow



Room  
10 ft. x 10 ft.



# BrewerAccess

EXAM TABLE

## Advanced Safety

- Stable four-leg patient step
- Largest step surface area available
- Simple pneumatic backrest for easy patient positioning

## Abundant Storage

- Largest storage capacity available
- Patented pass-through side drawers accessible from right or left side
- Front drawer: 19.8" x 19.2" x 3.8" (50.5 cm x 48.9 cm x 9.8 cm)
- Side drawers: 22.6" x 20.2" x 7" (57.5 cm x 51.4 cm x 17.8 cm)
- Slide through holder accepts 21" & 18" rolls (53.3 cm & 45.7 cm)
- Additional paper storage under backrest

## Table Dimensions

- Height: 32" (81.3 cm)
- Length: 57.3" extends to 71" (145.5 cm - 180.3 cm)
- Width: 28" (71.1 cm)
- Leg extension length: 15" (38.1 cm)
- Step Dimensions: 20.4"W x 17.7"L x 8"H (51.8cm x 45 cm x 20.3cm)



## Other Specifications

- 500 lb (226 kg) weight capacity
- Standard five-year warranty
- Standard adjustable stirrups
- 115VAC, 60Hz, 5.8 amps max - outlet
- Easy-to-clean seamless upholstery
- 15 standard upholstery color options
- CAL-133 fire resistant upholstery available in 4 colors
- Custom colors available
- Optional 5° pelvic tilt, 60-Watt drawer heater in front top drawer


## models

- |       |  |
|-------|--|
| #5000 | Access Exam Table with pneumatic back                                |
| #5001 | Access Exam Table with pneumatic back, pelvic tilt and drawer heater |

Drive safety, efficiency and patient satisfaction in the examination room.



## Floor Plans & Work Flow



# LiftMate™

## Brewer Low/High Patient Lift

### Advanced Safety

- Large non-slip platform surface
- Two emergency stop buttons
- Grab bars for patient stability and positioning
- Four auto-lock rubber feet for stability
- Four retractable swivel casters prevent rolling during use
- Smooth movement for enhanced patient security
- Quiet actuators for patient comfort
- Large platform with ample room for patient positioning
- Easy access for walkers
- Designed to prevent patient falls
- Reduces need to assist patients onto exam tables

### Dimensions

- Height: 38" - 47" (96.5 cm-119.4 cm)
- Depth (outer dimensions): 21.5" (54.6 cm)
- Width (outer dimensions): 34.75" (88.3 cm)
- Depth (platform): 18" (45.7 cm)
- Width (platform): 21.5" (54.6 cm)
- Width (between grab bars): 27.75" (70.5 cm)

### Other Specifications

- Standard three-year warranty
- 500 lb (226 kg) weight capacity
- 115 VAC, 60 Hz
- Safe, extra-low voltage actuators
- Convenient hand pendant for easy operation
- Powered by standard outlet plug or rechargeable battery
- Fits through standard doorways

- **Helps prevent patient falls**
- **Improves clinical staff safety**
- **Removes patient mobility barriers**
- **Reduces capital investment**
- **Mobile, maneuverable and easy to use**

*Lift with your finger, not your back.*



model

#3000 Brewer LiftMate™ Low/High  
Mobile Patient Lift

Protect your patients as well as your staff. Lift with your finger, not your back.



HELPS PREVENT  
PATIENT FALLS

MOBILE,  
MANEUVERABLE AND  
EASY TO USE

IMPROVES CLINICAL  
STAFF SAFETY

REDUCES CAPITAL  
INVESTMENT

REMOVES PATIENT  
MOBILITY BARRIERS



www.brewerliftmate.com

# Additional Brewer Product Lines

## Lighting

## Medical Seating

Millennium Series

Century Series

Value Plus Series

Tradition Series

Side Chairs

Blood Draw Chairs

## IV/Infusion Pump Stand

## Hamper Stands

## Utility Carts


## Instrument Stands

## Kick Bucket

## Overbed Tables

## Step Stools

## Waste Cans



## LED ECO-SERIES SERIES

### Light

- Bright, best-in-class white light drives clinical efficiency
- Cool to the touch and easy to maneuver
- 6" long, 5" dia., tri-lobed face
- Maximum 170° left/right head rotation
- 12-Watt LED bulb included
- Color temperature - 4000° Kelvin
- 8800 Lumens; 5" Beam Dia. at 18"
- Bulb life: Approx. 50,000 hrs

### Arm

- Flexible, "no-drift," 23" (58.4 cm) gooseneck
- Maximum Height: 66" (167.6 cm)
- Smooth, easy-to-clean

### Base

- Internal cord exits weighted base at floor
- Optional five-caster base
- Base Diameter: 13" (33 cm)
- Cord length: 10' (304.8 cm)

## HALOGEN 35 SERIES

### Light

- Crisp white light
- Non-metallic head eliminates surface heat
- Dichroic film-coated reflectors filter out 80% of beam heat
- Heat sink technology removes bulb socket heat and keeps the housing cool
- 35-Watt bulb included
- Color temperature - 3000° Kelvin
- 1500 Footcandles @ 24" (61 cm)
- Bulb life: 4000 hrs

### Arm

- Flexible, "no-drift," 24" (61 cm) gooseneck
- 540° rotation of light head and arm
- Maximum Height: 65.5" (166.4 cm)

### Base

- Internal power cord exits base at floor
- Optional five-caster base (shown)
- Base Diameter: 13" (33 cm)
- Cord length: 9' 3" (281.9 cm)

## INCANDESCENT SERIES

### Light

- Advanced reflector system provides shadow-free, crisp white light
- Non-metallic head eliminates surface heat
- 25-Watt incandescent bulb included
- Bulb life: 500 hrs

### Arm

- Flexible, "no-drift," 18" (45.7 cm) gooseneck offers superior beam alignment
- 540° rotation of light head and arm
- Maximum Height: 62" (157.5 cm)

### Base

- Internal power cord exits base at floor
- Optional five-caster base
- Base Diameter: 13" (33 cm)
- Cord length: 9' 3" (281.9 cm)

## CHROME SERIES

### Light

- 10 to 40-Watt incandescent bulb (not included)

### Base

- Stable chrome plated base
- Base Diameter: 11" (27.9 cm)
- Grounded 3-wire plug
- Cord length: 12' (365.8 cm)

### Arm

- Durable chrome finish
- 11" (27.9 cm) gooseneck
- Maximum Height: 72" (182.9 cm)

Some assembly is required



# MILLENNIUM SERIES

Premium specialty seating with ergonomic design.

## Stool Specifications

- 18.5" (47 cm) wide contoured seat cushion
- 2" (5.1 cm) dual-wheel nylon casters with polyurethane tread; locking caster optional
- 23" (58.4 cm) five leg powder coated black aluminum base
- Weight capacity: 250 lb (113 kg)
- Standard five-year warranty



## Task Stool

- Pneumatic height adjustment: 18" - 24.75" (48.7 cm - 62.9 cm)
- Adjustable backrest with synchronized traversing movement of seat
- Adjustable armrests optional

## Lab Stool

- Pneumatic height adjustment: Glides: 20.75" - 31" (52.7 cm - 78.7 cm)
- Casters: 21.75" - 32" (55.2 cm - 81.3 cm)
- 18" (48.7 cm) foot ring standard
- Adjustable backrest

## Surgeon Stool

- Hands-free foot pneumatic height adjustment: 19" - 24" (48.3 cm - 61 cm)
- 18.5" (47 cm) poly-foam seat cushion on 21340B only; 16" (40.6 cm) wide on 21340
- Adjustable backrest on 21340B

# CENTURY SERIES

Preferred by medical facilities worldwide.

## Stool Dimensions

- 6" (40.6 cm) wide, 4" (10.2 cm) thick poly-foam seat cushions
- Five-leg, 23" (58.4 cm) brushed aluminum accent base
- Pneumatic height adjustment; two heights available:
  - 17.75" - 24.5" (45.1 cm - 62.2 cm)
  - 20" - 27.5" (50.8 cm - 69.9 cm)

## Other Specifications

- Optional adjustable contoured backrest with black plastic shroud
- Adjustable 18" (48.7 cm) foot ring standard on higher height stool models
- Weight capacity: 250 lb (113 kg)
- Standard five-year warranty



# VALUE PLUS SERIES

Affordable without compromising style and reliability.

## Stool Specifications

- 15" (38.1 cm) wide, 4" (10.2 cm) thick poly-foam seat cushion
- Pneumatic or spinlift height adjustment: 17" - 21.25" (43.2 cm - 54 cm)
- Five-leg, 21" (53.3 cm) black ABS plastic composite base
- 2" (5.1 cm) black plastic dual-wheel nylon casters
- Optional adjustable contoured backrest with black plastic shroud
- Standard five-year warranty



# TRADITION SERIES

Economical exam stools that stand the test of time.

## Tradition Series Stool Specifications

- 14" (38.1 cm) wide, 3.5" (10.2 cm) thick poly-foam seat cushion
- Spinlift height adjustment: 19.5" - 27.5" (49.5 cm - 69.9 cm)
- Sturdy four-leg, 19.5" (49.5 cm) chrome base
- 2" (5.1 cm) rubber wheel, ball bearing casters
- 47 standard upholstery color options
- Standard five-year warranty



# SIDE CHAIRS SERIES

Comfort. Reliability. Durability.

## Chair Specifications

- Kick-out leg design helps prevent wall dings and dents
- Comfortable, 2" (5.1 cm) thick poly-foam seat cushion
- Seat 18" x 18" (45.7 cm x 45.7 cm)
- Contoured backrest 18" x 13.38" (45.7 cm x 34 cm)
- Durable steel tube frame
- Weight capacity: 250 lb (113 kg)
- Standard five-year warranty
- Meets BIFMA (Business and Institutional Furniture Manufacturer's Association) Standards

## Model with arms:

- Inside arm width: 20" (50.8 cm)
- Outside arm width: 23" (58.4 cm)



# BLOOD DRAWING CHAIR SERIES

Comfort. Adjustability. Durability.

## Chair Specifications

- Weight capacity: 350 lb (158 kg) weight capacity
- Standard five-year warranty
- 2" (5.1 cm) thick poly-foam seat cushion
- Two standard add-on work surfaces: 8" L x 12" W (20.3 cm x 30.5 cm)
- Work surfaces are easily removable to accommodate space constraints
- Optional drawer attaches on either side of chair
- Work surfaces covered in a bacteriological, chemical, and fungal resistant material
- Meets BIFMA standards

## Chair Dimensions

- Height: 37" (94 cm)
- Length: 23.5" (59.7 cm)
- Width: 20" (50.8 cm), 29.25" (74.3 cm) with armrests, 43.25" (109.9 cm) with both work surfaces
- Floor to top of seat: 22" (55.9 cm)
- Chair depth: 23.5" (59.7 cm)
- Armrest Angle: 15°
- Armrest Height: 29.5" - 38" (74.9 cm - 96.5 cm)
- Seat cushion: 18" x 18" (45.7 cm x 45.7 cm)
- Backrest cushion: 18" x 12.5" (45.7 cm x 31.8 cm)

## OPTIONAL

## BLOOD DRAWING CHAIR DRAWER

### Drawer Specifications

- Easily attaches to either side of chair
- Interior plastic lining for easy cleaning

### Drawer Dimensions

- Length: 18" (45.7 cm)
- Height: 4" (10.2 cm)
- Width: 8" (20.3 cm)



# SUPPORT PRODUCTS

Short Wheel Infusion Pump Stand



Large Infusion Pump Stand



Heavy Base Four-Leg IV Pole



Aluminum Base Five-Leg IV Pole



Two-Hook Four-Leg IV Pole



18" Round Linen Hamper



18" Tilt-Top Round Linen Hamper



Tilt-Top Square Linen Hamper



Adjustable Folding Linen Hamper



Stainless Steel All Purpose Cart



Plastic Utility Cart



Multi Purpose Cart



Mobile Equipment Cart



All Support Products • 5 Year Warranty

Mobile Instrument Stand



Mayo Instrument Stand



Mayo Instrument Stand



Stainless Steel Mayo Stand



Kick Bucket w/Frame



"U" Base Table



"H" Base Table



Low "H" Base Table



Step Stool



Step Stool with Hand Rail



Bariatric Step Stool with Hand Rail



Two-Step, Step Stool with Hand Rail



Large Square Steel Waste Can



Small Square Steel Waste Can



Large Round Steel Waste Can



Large Rectangular Plastic Waste Can



Small Rectangular Plastic Waste Can



All Support Products • 5 Year Warranty


## Standard Upholstery Color Selections

For Brewer Tables, Side Chairs, Blood Drawing Chairs and Exam Stools



## California Technical Bulletin 133 Upholstery Color Selections (CAL-133)


For Brewer Tables, Century Series Stools and Value Plus Series Pneumatic Stools



CAL TB 133 is a written legal standard outlining flammability test procedures for upholstered furniture. These vinyl color selections contain a flame blocking barrier that helps to prevent ignition, reflects heat and is flame retardant.

## Premium Plush Upholstery Color Selections

For Brewer Assist Procedure Tables, Access High-Low Exam Tables, Side Chairs and Century Series Stools




All Premium Plush Upholstery Colors are Ultrafabrics ULTRALEATHER™.

**Note:** Tables featuring a plush upholstery top are two-toned, stitched, with Stucco and one of the premium plush upholstery colors. Exam stools and side chairs are one solid color, stitched. (-UL)

For Actual Material Swatches call 1•888•Brewer•1

## Additional Standard Upholstery Color Selections

For Brewer Exam Stools Only



For Actual Material Swatches call 1•888•Brewer•1

# Brewer

*The Power to Advance*

The Brewer Company, LLC  
N88W13901 Main Street, Suite 100  
Menomonee Falls, WI 53051

**P** 1.888.Brewer.1

**F** 262.251.2332

[www.brewercompany.com](http://www.brewercompany.com)

