

# Brewer

## Assist PRO & Assist Power Procedure Tables®

### Specifications

#### Brewer Differentiators

- *Adjustable headrest design accommodates a broad range of clinical procedures*
- *Enables quick transitioning from one type of procedure to the next for increased throughput*
- *Tilt range of up to 30° allows true Trendelenburg positioning*
- *ADA-compliant, 19" (48.3 cm) seat height helps enable safe lateral wheelchair transfers*
- *Standard 3-year warranty*
- *Made in USA*

#### Product Specifications

- Weight capacity: 450 lbs. (204 kg)
- Seat Height: 19" - 40" (48.3 cm - 101.6 cm)
- Length: 69" (75.3 cm) when flat
- Length with headrest & leg rest extended: 80" (203.2 cm)
- Length in upright position: 40" (101.6 cm)
- Width: 28" (71.1 cm)
- Leg rest: 17" x 12" (42.2 cm x 30.5 cm)
- Headrest: 13" x 11" (33 cm x 27.9 cm)
- Tilt range / Trendelenburg: to 30°
- Base: 24" W x 36" L (61 cm x 91.4 cm)
- Power cord: 70" (177.8 cm)
- 115VAC, 60Hz
- Model 7000 also available in 230 Volt



Assist Power  
Procedure Table®

AssistPRO Programmable  
Power Procedure Table®

#### Standard Features

- Industry-best 3-year warranty
- Adjustable ergonomic stirrups
- Foot control
- Hand pendant standard on AssistPRO (7500 models)
- Safe extra-low voltage actuators
- ADA-compliant 19" (48.3 cm) low height
- 40" (101.6 cm) maximum height
- Leg rest extends an additional 5" (12.7 cm)
- Adjustable three-point pivot headrest
- Easy-to-clean seamless upholstery
- Multi-density, bolstered support upholstery
- 18" (45.7 cm) or 21" (53.3 cm) paper dispenser

## Assist PRO & Assist Power Procedure Tables®

### Specifications



#### Certifications

- UL 60601-1 compliant (120 Volt models only)

#### Additional Specifications for AssistPRO Model 7500

- Intuitive hand pendant and foot control included
- Four preferred positions can be programmed into the table's memory
- Simple two-step programming process
- Auto-run toggle for safety



*Hand Pendant*



*Foot Control*

#### Models

- 7500 - Brewer AssistPRO Power Procedure Table
- 7000 - Brewer Assist Power Procedure Table

#### Optional Accessories

- 99497 - Grab bars (pair)\*
- 99499 - Articulating armboard\*
- 99500 - Hand pendant kit for Assist 7000
- 99501 - Side seat rails (pair of 2)
- 99502 - One universal socket clamp
- 99504 - Deluxe articulating knee crutches (left & right)\*\*
- 99505 - Stainless steel treatment pan
- 100253 - Urology drain pan system
- 101726 - Hand pendant replacement kit for AssistPRO 7500

\* Needs one order of #99501

\*\* Needs one order of #99501 and two orders of #99502

#### Optional Accessories



*Safety Grab Bars*



*Articulating Armboard*



*Deluxe Articulating Knee Crutches*

The Power to Advance. Safely. Financially. Competitively.

www.brewercompany.com | P: 1.888.273.9371

The Brewer Company is an ISO 13485:2003 certified company.

LIT-M-009-Rev A