

What can Brewer's latest
beyond-the-box thinking do for you?



By delivering all the functionality of a power table in a remarkably
efficient solution, Brewer redefines exam table flexibility & value.

The Power to Advance



FLEX[™] Access Exam Table

Brewer

Step up to power exam table performance for less than the price of a standard, fixed-height box table... all in an incredibly flexible package.

The Brewer FLEX™ Access Exam Table revolutionizes your exam room opportunities.

FLEX transforms your practice by delivering the kind of flexibility you can't get with a fixed-height exam table... and the standard value-added features you won't find on competitors' offerings.

Provide a better patient experience.

Concerned about the growing challenges of safe, efficient patient handling in the exam room? FLEX makes it easy and affordable to offer patients a higher quality care experience.

- **Safely handle patients** with the touch of the foot-operated controls.
- **Achieve ADA compliance** even in exam rooms with very small footprints.
- **Be recognized for bariatric excellence** thanks to an industry-leading 700-lb. patient weight capacity.
- **Enhance patient satisfaction** by offering them increased comfort and confidence.

Drive your bottom-line success.

FLEX offers distinct competitive advantages designed to optimize your practice revenue.

- **Expand your patient population** by accommodating a broader range of patients.
- **Perform a wider range of examinations** than would be possible with non-powered tables.
- **Improve clinical throughput** so you can see more patients per day.
- **Unmatched ergonomic patient access** helps keep you and your staff healthy and more productive.



Available in your choice of models: 5700 Series with manual backrest and 5800 Series with powered backrest.




SafeGlide™ leg extension enhances patient & staff safety.

Easy-to-manuever design accommodates patients weighing up to 700 lbs, while allowing caregivers to slide the extension in and out with minimal effort.



Compact footprint still offers plenty of storage capacity.

Holds up to six 18" paper rolls or other frequently needed supplies so they are always within easy reach. Helps drive clinical efficiency.



Widest variety of standard color choices.

Conveniently match our standard, seamless upholstery to your color scheme.

Add flexibility with options to suit your practice:



Patient Safety Grab Bars

Provide stable and secure grab point to enhance patient safety during transfers. Conveniently swing out and are easily removable.



Stirrups, Pelvic Tilt & Drawer Heater

Cost-effectively include ergonomic options as needed for OB/GYN, urology and other specialties.



Clinical Work Surface

Keep laptops and instruments within easy reach for greater clinical efficiency. Pass-through design allows access on either side.

By redefining exam table flexibility and value, Brewer delivers the Power to Advance. Safely. Financially. Competitively.

Brewer

Visit brewercompany.com

With *FLEX*™, less is so much more.

With Brewer's *FLEX* Access Exam Table, you wind up paying less than the price of a traditional, fixed-height exam table and you get so much more in the process.

Receive substantial tax credits.

Unlike traditional, fixed-height exam tables, power exam tables qualify for significant tax credits. As a result, the Brewer *FLEX* Access Exam Table provides you with industry-leading, practice-enhancing power table performance for less than the cost of a typical exam table.

Optimize reimbursement.

With the implementation of the Affordable Care Act, reimbursement dollars are now tied directly to patient safety and satisfaction ratings. By enhancing patient comfort, confidence and dignity, *FLEX* helps you improve the overall patient experience, which can help you maximize your reimbursements.

Enhance your competitive position.

FLEX enables you to attract a broader patient population, perform a wider range of examinations and see more patients per day than would be possible with fixed-height exam tables.

And as practice needs change, *FLEX* can conveniently and cost effectively adapt with your practice.

The Ultimate Exam Table Value

<i>FLEX</i> Access Exam Table	
Typical Pricing (Includes resellers markup)	\$4,500-\$5,000
Section 44 Disabled Access Tax Credit	– \$2,250-\$2,500*
Section 179 Deduction for Certain Capital Expenditures	– \$630-\$700*
Final Customer Price	\$1,620-\$1,800

Standard Fixed-Height
Exam Table \$1,700-\$2,000**

* Visit <http://brewercompany.com/wp-content/uploads/Tax-Credit-Flyer.pdf> for more detailed information regarding qualifying tax credits.

**Includes resellers markup.



Brewer offers the only power exam table line capable of accommodating bariatric patients weighing up to 700 lbs. 5700 Series with manual backrest shown.

FLEX also comfortably accommodates patients over 6' 6" tall. 5800 Series with powered backrest shown.

Brewer

The Brewer Company, LLC
N88 W13901 Main Street, Suite 100
Menomonee Falls, WI 53051

1.888.273.9371
www.brewercompany.com

