

- Sloped and cushioned baby-cradle edges keep younger children safely in place during exams

- Large side drawer offers 3.7 cubic feet of space to comfortably store pediatric scale and clinician essentials



- Padded leg extension accommodates infants to older children and teen

Product Specifications

- Weight capacity: 500 lbs. (226 kg)
- Seat height: 32.5" (82.55 cm)
- Pediatric table length: 51"; expandable to 65" (129.5 cm-165.1 cm)
- Pediatric table width: 26" (66 cm)
- Leg extension length: 14" (35.5 cm)
- Front storage drawer: 21.3" x 19.75" x 4.25" (59.9 cm x 50.1 cm x 9.7 cm)
- Side storage drawer: 23.6" x 35.75" x 7.6" (57.5 cm x 90.8 cm x 19.3 cm)
- Paper roll holder: accommodates up to 21" paper rolls
- UL certified: IEC60601-1 Ed. 3.1

Standard Features

- Industry-best 5-year warranty
- Easy-connect, removable upholstery top supports infection control
- Large, sturdy patient step with integral anti-slip surface offers easy cleanability and infection control
- Easy-to-clean, seamless upholstery available in 12 standard color choices
- Numerous additional custom colors available to choose from

Models

- Pediatric model without stirrups: 1000PX-XXDX-XXX-__ [2 digit color code]
- Pediatric model with stirrups: 1000PX-XSDX-XXX-__ [2 digit color code]

Optional Accessories

- 3500288 – Front Drawer Dividers
- 3500347 – Brewer Welch Allyn® Exam Light / IV Bracket (factory installed)
- 3500398 – Brewer Welch Allyn® Exam Light / IV Bracket (field installed)

* Product specification measurements are approximate and may vary +/- 0.5"



We Help Healers Heal

Learn more at brewercompany.com