

# Access<sup>®</sup> 1000 Fixed Height Exam Table



## Specifications

### Highly Customizable

- Easy-connect, removable upholstery top is available in standard & pediatric configurations
- Add options as needed to economically meet practice requirements

### Improved Patient Experience

- Large, sturdy step with integral, anti-slip surface enhances patient safety
- Stable, cushioned platform with 500-lb patient weight capacity comfortably supports patients of all builds
- Optional multi-position patient-assist handles sized for more confident transitions and easier patient self-positioning
- Seamless upholstery & integral step surface offer easy cleanability and infection control

### More Ergonomic Access

- Rounded corners and softened edges throughout safeguard patients and caregivers, while improving patient access
- Pass-through side drawer purposefully positioned at waist height to reduce bending and twisting
- Front drawer strategically positioned so caregivers do not infringe on a patient's personal space

### Table Specifications

- **Weight Capacity:** 500 lbs (227 kg)
- **Seat height:** 32.5" (82.6 cm)
- **Length:** 58"; expandable to 72" (147.3 cm - 182.9 cm)
- **Pediatric table length:** 51"; expandable to 65" (129.5 cm - 165.1 cm)
- **Leg extension length:** 14" (35.6 cm)
- **Upholstery width:** 27.25" (69.2 cm)
- **Pediatric table width:** 26" (66 cm)
- 4.2 cubic feet (119 liters) of storage
- **Front drawer:** 20.9" x 18.8" x 3.8" (53.1 cm x 47.8 cm x 9.7 cm)
- **Side drawer:** 23" x 34.8" x 7" (58.4 cm x 88.4 cm x 17.8 cm)
- Paper roll storage capacity: accommodates up to 21" paper rolls
- Easy connect, removable upholstery top

- 120VAC, 60Hz outlet
- Industry-best 5 year warranty
- Easy-to-clean, seamless upholstery
- Two side drawer dividers for organization
- Large, sturdy patient step with integral anti-slip surface

### Optional Features

- Hospital-grade electric receptacle
- 60-Watt front drawer heater
- 5° pelvic tilt

### Certifications

- All Brewer products conform to US & Canadian electrical & medical industry codes & regulations and meet or exceed industry requirements.
- AAMI ES60601-1 (2005) + AMD 1 (2012)
  - CAN/CSA C22.2 No. 60601-1 (2014)
  - IEC60601-1 Ed. 3.1

## Models:

### Recommended Models:

- **1000A**  
Standard upholstered top, pneumatic back, pass-through rear drawer, ergonomic stirrups and 120V electrical receptacle
- **1000A-NR**  
Standard upholstered top, pneumatic back, pass-through rear drawer and ergonomic stirrups
- **1001A**  
Standard upholstered top, pneumatic back, pass-through rear drawer, ergonomic stirrups, pelvic tilt, 120V electrical receptacle and drawer warmer
- **1001A-NR**  
Standard upholstered top, pneumatic back, pass-through rear drawer, ergonomic stirrups, pelvic tilt and drawer warmer
- **1000P**  
Pediatric upholstered top and pass-through rear drawer

## Accessories:

- **3500920:**  
Assist Handles (factory installed)
- **98464:**  
Stainless Steel Treatment Pan
- **3500286:**  
Deluxe articulating knee crutches with assembly kit
- **100253-TG:**  
Urology Drain Pan System
- **3500288:**  
Front Drawer Dividers
- **3500347:**  
Brewer Welch Allyn<sup>®</sup> Exam Light / IV Bracket (factory installed)
- **3500398:**  
Brewer Welch Allyn<sup>®</sup> Exam Light / IV Bracket (field installed)
- **3500305:**  
Paper Fastener Strap
- **3500077:**  
Footrest Pad

# Create an elevated patient-provider experience.

Recommended Models	1000A	1000A-NR	1001A	1001A-NR	1000P
Weight Capacity	500	500	500	500	500
Standard or Pediatric	Standard	Standard	Standard	Standard	Pediatric
Patient Assist Handles	Optional	Optional	Optional	Optional	
Stirrups	☑	☑	☑	☑	
Drawer Warmer			☑	☑	
Pelvic Tilt			☑	☑	
Storage Capacity	4.2 cu ft	4.2 cu ft	4.2 cu ft	4.2 cu ft	4.2 cu ft
Side Drawer Dividers	☑	☑	☑	☑	☑
Electrical Receptacle	☑		☑		

## Available Colors.



Clamshell  
21



Gunmetal  
22



Feather  
23



Black Satin  
24



Saddle  
25



Tapestry  
26



Cabernet  
27



Cocoa  
28



Deep Sea  
29



Ivy  
30



Deep Sapphire  
31



Blue Fog  
32

Need a closer look? Contact us to receive physical samples of our vinyl or base materials.

*Custom colors available.*



1-888-273-9371 | [BrewerCompany.com](http://BrewerCompany.com)

## The Brewer Difference.

Industry's most  
safety features.

Industry's only  
5-year warranty.

USA-made in  
the Midwest.