

# FLEX™ High-Low PLUS Power Exam Table



## Specifications

### Advanced Safety

- Exclusive SafeGlide™ leg extension increases safety and satisfaction
- ComfortGrip™ Patient Assist Handles facilitate easy transfers and adjustments, enhancing patient independence and safety during examinations (*optional*)
- Return to chair feature on power back models
- Anti-tamper on/off switch limits accidental motion
- Two safety zones prevent injuries during table movement
- 18" (45.7 cm) low height that raises to 37" (94 cm), both ADA compliant
- Easy on/off for patients
- Powered backrest adjusts for easy patient positioning\*

\* Pneumatic back models are also available

### Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density support upholstery
- Adjustable ergonomic stirrups (*optional*)
- Toe/leg extension & plastic debris tray

### Clinician Convenience

- Mayo Pass-Through work surface on side of table (*optional*)
- **Front drawer:**  
19.75" x 19.25" x 4" (50.2 cm x 63.5 cm x 10.2 cm)
- **Paper roll holder accepts:**  
18" & 21" rolls

### Table Specifications

- 700 lb (317 kg) weight capacity
- ComfortGrip™ patient assist handles (*optional*)
- Ergonomic adjustable stirrups come standard (*5701PLUS and 5801PLUS models*)
- Foot control comes standard
- Standard 5-year warranty
- **Height:** 18"-37" (45.7 cm – 94 cm) — Both heights are ADA compliant
- **Length:** 60", extends to 78" (152.4 cm – 198.1 cm)
- **Upholstery Width:** 28" (71.1 cm)
- **SafeGlide™ Leg extension length:** 18" (45.7 cm)
- 115VAC, 60Hz outlet
- Safe, low voltage actuators
- Up to 6° pelvic tilt
- 60-Watt front drawer heater
- Easy to clean, seamless upholstery

### Certifications

All Brewer products conform to US & Canadian electrical & medical industry codes & regulations and meet or exceed industry requirements.

- UL 60601-1
- CAN-CSA C22.2 No. 601.1
- ISO 13485:2016

## Accessories:

- **3500687PLUS:**  
EasyGlide™ Mobility System
- **99505:**  
Stainless Steel Treatment Pan
- **99502:**  
One Universal Socket Clamp
- **99504:**  
Deluxe Articulating Knee Crutches\*
- **100253-TG:**  
Urology Drain Pan System
- **100407PLUS:**  
Front Drawer Dividers
- **3500077:**  
Footrest Pad
- **100540:**  
Two Support Adapters for Knee Crutches (factory installed)
- **100550:**  
Two Support Adapters for Knee Crutches (field installed)
- **102031PLUS:**  
Brewer Welch Allyn® IV Bracket (factory installed)
- **102032PLUS:**  
Brewer Welch Allyn® IV Bracket (field installed)  
\*Needs two 99502 and one of either 100540 or 100550

## Optional Kits








- **2102894PLUS:**  
Stirrups for 5700PLUS/5800PLUS (set of two factory installed)
- **2102898PLUS**  
Stirrups for 5700PLUS/5800PLUS (set of two field installed)
- **2102891PLUS**  
ComfortGrip™ Patient Assist Handles (set of two, factory installed)
- **2102895PLUS**  
ComfortGrip™ Patient Assist Handles (set of two, field installed)
- **2102892PLUS**  
Mayo Pass-Through Work Surface (factory installed)
- **2102896PLUS**  
Mayo Pass-Through Work Surface (field installed)



# Create an elevated patient-provider experience.

Models	5700PLUS	5701PLUS	5800PLUS	5801PLUS
Weight Capacity	700	700	700	700
ADA Compliant 18" Seat Height	✓	✓	✓	✓
Power Base	✓	✓	✓	✓
Power Back			✓	✓
Manual Back	✓	✓		
Patient Assist Handles	Optional	Optional	Optional	Optional
Stirrups	Optional	✓	Optional	✓
Pelvic Tilt & Drawer Warmer		✓		✓
Outlet		✓		✓
Work Surface	Optional	Optional	Optional	Optional
Mobility	Optional	Optional	Optional	Optional

## Available Colors.

 Chateau Grey AF906	 Flannel AF945	 Horizon AF928	 Mist AF938	 Black Satin PR40
 Cool Grey MV821	 Eclipse AF919	 Blue Fog PR46	 Cocoa PR51	 Cabernet PR36
 Feather PR37	 Deep Sapphire PR55	 Gunmetal PR92	 Clamshell PR57	 Ivy PR53

Need a closer look? Contact us to receive physical samples of our vinyl or base materials.

*Custom colors available.*



1-888-273-9371 | [BrewerCompany.com](https://www.BrewerCompany.com)

## The Brewer Difference.

Industry's most  
safety features.

Industry's only  
5-year warranty.

USA-made in  
the Midwest.

LIT-M-068 RevA