

Brewer

Clinical & Exam Room Solutions

We Help Healers Heal.



Introducing Versa:
the versatile
fixed-height exam
table designed
to serve your
entire population
(See page 14)

Making a difference by design

Relentlessly adding value to your clinical care.

We understand that today's healthcare reality presents a host of challenges to your clinical, financial and operational success. That's why we work closely with medical professionals to develop the clinical and exam room advancements that matter.

Redefining patient & staff safety.

Our solutions set the standards for safe, confident handling of even your most problematic patients, including elderly, limited mobility, and bariatric populations. At the same time, they minimize caregiver handling, safeguarding the health of your staff and your bottom line.

Improving your competitive position.

Innovative designs and capabilities enable you to see the broadest possible range of patients and provide a wider range of services. They also help preserve patient dignity and enhance patient satisfaction – from the exam room to the procedure suite.



Model 6801 in Clamshell



Optimizing your return on investment.

At Brewer, expert ergonomic design, quality materials and made-in-USA craftsmanship come together to create lasting values. Versa is a perfect example. It's the first major design innovation of a fixed-height exam table in the industry in more than two decades. And it redefines patient and caregiver access, safety and clinical efficiency. Plus, you can count on Brewer for helpful touches, like easy-to-clean surfaces to expedite infection control and industry-leading warranties for greater peace of mind.

So be sure to take a good look at all the ways Brewer adds value to your clinical care.

TABLE OF CONTENTS

Power Exam Tables

Access® High-Low 700 Exam Table	8
Access® High-Low Exam Table	8
FLEX™ Access High-Low Exam Table ...	10

Power Procedure Tables

Assist PRO Power Procedure Table® ...	12
Assist Power Procedure Table®	13

Fixed-Height Exam & Treatment Tables

Versa Exam Table	14
Access® Exam Table	16
Element Treatment Table	17

Medical Seating

135 Saddle Seat Series	18
Polyurethane Seating	19
Millennium Series	20
Century Series	22
Value Plus Series	23
Side Chairs	24
Blood Drawing Chair	24

Task Lighting

LED Eco Series	24
----------------------	----

Support Products

Waste Cans	25
Utility & Equipment Carts	26
Hampers	27
Step Stools	28
Instrument Stands	29
Mayo Stands	29
Over Bed Tables	29
Kick Buckets	29
Infusion Pump Stands	30
IV Poles	30

Upholstery Color Options.....	31
-------------------------------	----

Doing more to safeguard your patients & your practice

Leading the way in safety.

Patient handling can be a risky proposition for patients and staff. That's why Brewer is proud to take the lead in introducing the advancements that help assure patient and staff safety throughout so many exam and procedure environments.

More confident patient transfers

Only Brewer offers a power exam table that allows the seat back to be raised in the chair position when at its lowest seat height. The upright seat back, combined with patient safety grab bars – standard, provides three points of contact during wheelchair transfers. This unique capability helps assure more confident patient transfers, enhances patient independence, and minimizes physical demands on caregivers.





More convenient positioning for exams

Brewer's SafeGlide Leg Extension system enables caregivers to easily position the extension with their fingertips. heavy-duty design supports patients up to 700 pounds.

Also, the patient safety grab bars rotate nearly 180 degrees to facilitate patient positioning during exams.

Unmatched accident protection

Only Brewer's Access High-Low Exam Table features three safety switches, providing the most thorough protection available against accidental injury.



Creating a culture of safety for your practice

Safer products are just the beginning.

Healthcare professionals want the safest environment possible for their patients and staff. And Brewer understands just how important the team behind the products can be to making that goal a reality.

Safe patient handling expertise

Our Safe-Patient-Handling Specialists are trained to assess your care requirements in order to match the right equipment solutions to your particular needs.

On-site training

Our specialists also ensure that care teams get the expert training needed to make effective use of all the safety features Brewer incorporates into our products.

Responsive support

We also are there for you after the sale, with industry-leading warranties and a team of experts dedicated to helping you maximize the return on your safety investment.





Brewer

Brewer

Visit brewercompany.com

Power Exam Tables

Advancing patient & staff safety.

Brewer Access® High-Low Exam Tables set the standard for safe patient handling. Our 18" (45.7 cm) profile in the seated position offers level transfers from wheelchairs while the rotating patient safety grab bars and upright seat back enable more confident transfers, with less caregiver interaction. Brewer's exclusive SafeGlide™ leg extension allows caregivers to slide the heavy-duty extension in and out with minimal effort. Seat height ranges from 18"-39" (45.7 cm-94 cm) accommodate various examination needs, while our innovative pass-through drawers and work surface give you the flexibility to make the most of compact exam room footprints. Add in more standard features and industry-best, 3-year warranties and it's easy to see how Brewer offers you superior value.

Access® High-Low Exam Table 700

Meet a wider range of patient needs with an industry-best 700-lb (317 kg) weight capacity. Includes unmatched storage capacity and swing-out patient safety grab bars standard.

Models

6800* - Access High-Low 700 with power back

6801 - Access High-Low 700 with power back, pelvic tilt, front drawer heater and hospital-grade outlet



Model 6801 in Cabernet

Access® High-Low Exam Table

A choice of pneumatic or power back models enables you to select the design that best suits your practice needs. Includes 450-lb weight capacity, unmatched storage capacity and swing-out patient safety grab bars standard.

Models

6000* - Access High-Low with pneumatic back

6001 - Access High-Low with pneumatic back, pelvic tilt, and front drawer heater

6500* - Access High-Low with power back

6501 - Access High-Low with power back, pelvic tilt, front drawer heater and hospital grade outlet

* This model is also available in 230 Volt.



Model 6501 in Cocoa



Model 6501 in Black Satin

Standard Features



Patient safety grab bars



Pass-through Work Surface



SafeGlide Leg Extension



Optional footrest pad adds comfort and support for enhanced patient experience

**Accessories
(All Access® High-Low Exam Tables)**

- 99505 - Stainless steel treatment pan
- 99502 - One universal socket clamp
- 99504 - Deluxe articulating knee crutches**
- 100253 - Urology drain pan system
- 100407 - Front drawer dividers
- 100540 - Two support adapters for knee crutches - *(factory installed)*
- 100550 - Two support adapters for knee crutches - *(field installed)*
- 102031 - Brewer Welch Allyn® exam light/IV bracket - *(factory installed)*
- 102032 - Brewer Welch Allyn® exam light/IV bracket - *(field installed)*
- 101800 - Electrical outlet (6001 only) - plant-installed - must order w/table
- 3500077 - Footrest pad

**Needs two orders of #99502 and one order of either #100540 or #100550



Visit brewercompany.com

FLEX™ Access High-Low Exam Tables

Step up to power exam table performance for less.

The FLEX™ Access High-Low Exam Table combines the productivity advantages of a power table with an industry-leading 700-lb. (317 kg) weight capacity in one affordable package. As a result, FLEX enables you to attract a broader patient population, perform a wider range of examinations and see more patients per day than would be possible with a fixed-height exam table. Its compact footprint and unmatched ergonomic patient access make it ideal for the confines of small exam rooms. Seat height ranges from 18"-39" (45.7 cm - 94 cm) accommodate your various examination needs. Available in two upholstery top choices, Premium and Euro Top. Plus, optional accessories – like patient safety grab bars, a pass-through clinical work surface and stirrups, pelvic tilt and a drawer heater – allow you to easily adapt FLEX to your clinical needs.

FLEX™ Access High-Low Exam Table – Pneumatic Back

Pneumatic lever allows you to manually raise and lower seat back to transition from consultation to examination.

Models

5700* - (Euro Top) | 5700-P* (Premium)
FLEX Access Exam Table with
pneumatic back

5701 - (Euro Top) | 5701-P (Premium)
FLEX Access Exam Table with pneumatic
back, pelvic tilt, front drawer heater,
ergonomic stirrups, and outlet



Model 5700 in Gunmetal

FLEX™ Access High-Low Exam Table – Power Back

Foot operated controls allow hands free raising and lowering of the seat back to optimize productivity as well as patient and staff safety.

Models

5800* - (Euro Top) | 5800-P* (Premium) FLEX Access Exam Table
with power back

5801 - (Euro Top) | 5801-P (Premium) FLEX Access Exam Table
with power back, pelvic tilt, front drawer heater, ergonomic
stirrups, and outlet

* This model is also available in 230 Volt.



Model 5801-P in Deep Sapphire

Available in two attractive upholstery tops



FLEX Premium Top



FLEX Euro Top



Pass-through Work Surface



Optional footrest pad (3500077) adds comfort and support for enhanced patient experience

Optional Kits



Stirrups, Pelvic Tilt & Drawer Heater



Patient safety grab bars

Optional Kits

- 2102894 - Ergonomic stirrup for 5700 and 5800
(factory installed)
- 2102898 - Ergonomic stirrup for 5700 and 5800
(field installed)
- 2102897 - Ergonomic stirrup for 5701 and 5801
(field installed)
- 2102891 - Safety grab bars *(factory installed)*
- 2102895 - Safety grab bars *(field installed)*
- 2102892 - Pass-through laptop work surface *(factory installed)*
- 2102896 - Pass-through laptop work surface *(field installed)*

Accessories

- 99505 - Stainless steel treatment pan
- 3500077 - Footrest pad
- 99502 - One universal socket clamp
- 99504 - Deluxe articulating knee crutches**
- 100253 - Urology drain pan system
- 100407 - Front drawer dividers
- 100540 - Two support adapters for knee crutches -
(factory installed)
- 100550 - Two support adapters for knee crutches -
(field installed)
- 102031 - Brewer Welch Allyn® exam light/
IV bracket - *(factory installed)*
- 102032 - Brewer Welch Allyn® exam light/
IV bracket - *(field installed)*

**Needs two orders of #99502 and one order of either #100540 or #100550

Power Procedure Tables

Boost your profitability through Brewer flexibility.

Now, you can accommodate a broad range of clinical procedures with a single, versatile Brewer power procedure table. That means added patient peace of mind – knowing a procedure can be performed in familiar surroundings.

All Brewer power procedure tables are infinitely adjustable to quickly transition from one type of procedure to the next for maximum throughput. A seat height of 19" (48.3 cm) helps enable safe patient transfers. The ability to extend the table from 70"-89.5" (177.8 cm - 227.3 cm) when headrest is fully extended allows you to handle a diverse range of patients. And multi-density, bolstered support upholstery and adjustable, three-point pivot headrest help optimize patient comfort.

Plus, all our power procedure tables feature a standard 3-year warranty.

AssistPRO Programmable Power Procedure Table®

Four preferred positions can be pre-programmed for quick, precise push-button transitions. Features 450-lb. weight capacity and simple two-step programming process. Intuitive hand pendant and foot control included.

Model

7500 - AssistPRO Programmable Power Procedure Table



Model 7500 in Cocoa with patient safety grab bars accessory



Multiple positions



Supine



Left Lateral Plane



Lithotomy



Proctology




Model 7500 in Blue Fog upholstery with footrest and pendant control

Optional Accessories



Patient safety grab bars



Articulating Armboard



Deluxe Articulating Knee Crutches

Assist Power Procedure Table®

All our power procedure tables feature a 450-lb. (204-kg) weight capacity; adjustable stirrups; a tilt range of up to 30° for true Trendelenburg positioning; and a standard three-year warranty.



Model 7000 in Black Satin

Model

7000* - Assist Power Procedure Table

Assist & AssistPro Accessories

- 99497 - Patient safety grab bars (pair)**
- 99499 - Articulating armboard**
- 99500 - Hand pendant kit for Assist 7000
- 99501 - Side seat rails (pair of 2)
- 99502 - One universal socket clamp
- 99504 - Deluxe articulating knee crutches (left & right)***
- 99505 - Stainless steel treatment pan
- 100253 - Urology drain pan system
- 101726 - Hand pendant replacement kit for AssistPRO 7500
- 102024 - Brewer Welch Allyn® exam light/IV bracket (factory installed)
- 102025 - Brewer Welch Allyn® exam light/IV bracket (field installed)

* This model is also available in 230 Volt.

** Needs one order of #99501

*** Needs one order of #99501 and two orders of #99502

NEW *Versa Fixed-Height Exam Table*

One exam table for your entire patient population.

We understand that fixed-height exam tables continue to be a mainstay of care facilities. That's why we felt they deserved to incorporate the very latest in ergonomic design concepts. The result is Versa – the innovative solution that redefines care quality and efficiency for your entire patient population.

- **Improved Patient Experience** – thoughtful design touches include integral non-slip surface of the large, sturdy step and more strategic design and placement of multi-position patient assist handles
- **Highly Customizable** – cost-effectively meet your care needs now and over time
- **More Convenient Accessibility** – rounded corners and softened edges throughout help safeguard patients and caregivers while drawers are positioned to reduce caregiver bending and twisting.

Versa® Exam Table

Versa gives you the ability to readily and cost-effectively adapt to your changing needs. When it comes to care, the standard model provides a new standard of quality and efficiency for most practices. Easily convert to a pediatric top to safely and effectively care for children from newborns to teenagers.



Standard model in Maverick Gray

Recommended Models

Standard model: 1000AX-XXDX-XXX-__ [2 digit color code]

Standard model with pelvic tilt, front drawer heater and outlet:

1000AX-XSDT-RHB-__ [2 digit color code]

Standard model with assist handles:

1000AX-HXDX-XXX-__ [2 digit color code]

Pediatric model: 1000PX-XXDX-XXX-__ [2 digit color code]



Standard top in Maverick Gray



Standard top in Maverick Gray with patient assist handles, stirrups and pelvic tilt

Easily adapt to changing practice needs



Standard Top



Pediatric Top

Optional Accessories



Knee Crutches



Footrest Pad adds comfort and support for enhanced patient experience

Pediatric top in Maverick Gray



Patient Assist Handles for more confident transitions



Women's Health Package includes warming tray, pelvic tilt and ergonomic stirrups

Accessories

98464 - Stainless Steel Treatment Pan

3500286 - Deluxe articulating knee crutches
with assembly kit

100253 - Urology Drain Pan System

3500288 - Front Drawer Dividers

3500347 - Brewer Welch Allyn® exam light/IV Bracket
(factory installed)

3500398 - Brewer Welch Allyn® exam light/IV bracket
(factory installed)

3500237 - Paper fastener strap

3500077 - Footrest pad

3500369 - Standard to Pediatric Conversion Kit

3500408 - Pediatric to Standard Conversion Kit

Brewer

Visit brewercompany.com

Fixed-Height Exam & Treatment Tables

Traditional designs infused with Brewer's value-added touches.

Our fixed height exam tables feature the most generous storage space in their class to keep your clinic organized and running efficiently. Exclusive pass-through drawers offer easy access from either side. Our treatment table includes all steel frame construction for added dependability. All exam and treatment tables are outfitted with easy-to-clean seamless upholstery.

Access Fixed-Height Exam Table

Offers largest step surface in the industry for added safety.

Multi-density support upholstery enhances patient comfort.

Patented pass-through drawer system – standard – provides more than 6 ft³ of storage capacity as well as

convenient access from either side to easily adapt to exam room environment. Simple pneumatic backrest enables easy patient positioning on this 500-lb. (226-kg) capacity unit. All backed by a standard 5-yr. warranty.

Height: 32" (81.3 cm). Width: 28" (71.1 cm).

Models

5000 - Access Exam Table **w/electrical outlet**

5000-SP-CN - Access Exam Table **w/o electrical outlet**

5001 - Access Exam Table with pelvic tilt, front drawer heater and **electrical outlet**

5001-SP-CN - Access Exam Table with pelvic tilt and front drawer heater **w/o electrical outlet**

Accessories

98464 - Stainless steel treatment pan

100253 - Urology drain pan

100407 - Front drawer dividers

102022 - Brewer Welch Allyn® Exam Light IV bracket
(factory installed)

102023 - Brewer Welch Allyn® Exam Light IV bracket
(field installed)



Model 5000 in Blue Fog



Model 5001 in Gunmetal upholstery



Element Treatment Table

Features a 500-lb. (226-kg) weight capacity and a standard 3-yr. warranty. Two-piece top design allows single person delivery for added freight savings. Adjustable back models have 6 positions.



Model 2240 in Deep Sea

Models

	Flat Top	Flat Top Adjustable	Pillowtop Adjustable
Basic Frame	- 2010	2030	2040
w/Shelf	- 2110	2130	2140
w/Shelf & Drawers	- 2210	2230	2240

Accessories

10196 - Table paper strap (Qty 2)

100407 - Drawer separators (Drawer models only)



Visit brewercompany.com

NEW *Dynamic-Motion Saddle-Seat Technology*

Exclusive capability promotes back health.

Brewer is the only U.S. manufacturer to offer 360° Dynamic Motion technology to healthcare professionals. This ergonomic capability is designed to engage core muscles in order to help strengthen back muscles over time. As a result, it helps keep care professionals more comfortable and productive throughout the day.

Amazon Series

Features a narrower seat and a choice of two height ranges specifically designed to accommodate smaller-framed care professionals.



Model 135ASB Amazon in Raven Wing

Models

135 Amazon Models*	Adjustable Height Range	Recommended User Height
■ 135ASG	24"-32"	5' 3"-5' 11"
■ 135ASB	24"-32"	5' 3"-5' 11"
■ 135ASP	24"-32"	5' 3"-5' 11"
■ 135ASG-L	21"-27"	Under 5' 3"
■ 135ASB-L	21"-27"	Under 5' 3"
■ 135ASPL	21"-27"	Under 5' 3"

Jumper Series

Features a wider seat and a choice of two height ranges specifically designed to accommodate taller, larger-framed care professionals.



Model 135JSG Jumper in Dove Gray

Models

135 Jumper Models*	Adjustable Height Range	Recommended User Height
■ 135JSG	24"-32"	5' 3"-5' 11"
■ 135JSB	24"-32"	5' 3"-5' 11"
■ 135JSP	24"-32"	5' 3"-5' 11"
■ 135JSG-H	27"-37"	Over 5' 11"
■ 135JSB-H	27"-37"	Over 5' 11"
■ 135JSP-H	27"-37"	Over 5' 11"

"3 Ultraleather options to choose from"

* G - Dove Gray B - Raven Wing P - Papyrus



Model 135ASP Amazon in Papyrus upholstery
Saddle seats feature a 5-year warranty and a 250-lb. weight capacity.

Comfortable and durable.

Choice of styles and configurations lets you match the design to your care needs for the best possible fit. All feature rugged polyurethane construction for long life and easy cleaning.



Model PS-2

Round Stool

Features a 14" polyurethane seat and 5-leg poly base.

Models

- PR-1 - with stationary glides (16.5"-21.5")
- PR-1-C - with casters
- PR-2 - with adjustable foot ring and stationary glides (19.5"-27")
- PR-2-C - with adjustable foot ring and casters
- PR-3 - with high height cylinder (22.5"-32.5"), adjustable foot ring and stationary glides
- PR-3-C - with high height cylinder (22.5"-32.5"), adjustable foot ring and casters



Model PR-2

Task Stool

Features an 18" contoured polyurethane seat, 16" adjustable backrest and 5-leg poly base.

Models

- PS-1 - with stationary glides (16.5"-21.5")
- PS-1-C - with casters
- PS-2 - with adjustable foot ring and stationary glides (19.5"-27")
- PS-2-C - with adjustable foot ring and casters
- PS-3 - with high height cylinder (22.5"-32.5"), adjustable foot ring and stationary glides
- PS-3-C - with high height cylinder (22.5"-32.5"), adjustable foot ring and casters



Model PS-2

Task Stool with Seat Tilt

Features seat tilt for added comfort and supported range of motion, an 18" contoured polyurethane seat, 16" adjustable backrest and 5-leg poly base.

Models

- PT-1 - with stationary glides (17"-22.5")
- PT-1-C - with casters
- PT-2 - with adjustable foot ring and stationary glides (20.5"-28")
- PT-2-C - with adjustable foot ring and casters
- PT-3 - with high height cylinder (22.5"-32.5"), adjustable foot ring and stationary glides
- PT-3-C - with high height cylinder (22.5"-32.5"), adjustable foot ring and casters



Model PT-3

Medical Seating

Ergonomically designed to meet your practice's needs.

Seating tends to be one of the most under appreciated aspects of a clinic. Yet it can dramatically impact the productivity and health of your entire staff. That's why our ergonomic experts work closely with clinicians to develop specific comfort solutions for their specialized requirements. The Millennium Series offers the thoughtful design and premium features that help make every day on the job seem less like work. Plus, a standard five-year warranty means added peace of mind.

Millennium Series

Premium specialty seating features seamless upholstery, an 19" (48.3 cm) wide, contoured seat cushion and adjustable contoured backrest unless otherwise noted. For stability, all incorporate a 23" (58.4 cm) five-leg, powder-coated black aluminum base and 2" (5.1 cm) dual-wheel nylon casters with polyurethane tread unless otherwise noted. Weight capacity: 250-lbs. (113 kg).

Task Models

Fully adjustable forward seat tilt and synchronized seat and backrest tilt support neutral posture to reduce fatigue and back pain. Pneumatic height adjustment: 18"-24" (45.7cm-61 cm).

21435BA - w/ arm rests

21435B - w/o arm rests



*Model 21435BA
in Azure Blue*

A - Arm rests B - Backrest

Heights are measured from top of seat cushion to floor. Because our seating is custom built, heights and widths could vary within 1/2 of an inch of product listing.



Task Model 21435BA in Azure Blue



*Lab Model 21521B
in Blue Fog (DUR-026)*

Lab Models

Features 18" (45.7 cm) foot ring standard and pneumatic height adjustment:

- w/ glides: 21"-31" (53.3 cm-78.7 cm)
- w/ casters: 22"-32" (55.9 cm-81.3 cm)

21520B - w/ glides

21521B - w/ casters



*Model 21521B
in Blue Fog (DUR-026)*

Surgical Models

Offers hands-free, foot-operated pneumatic adjustment:

- 18.5"-23.5" (47 cm-59.7 cm) with back
- 20"-25" (50.8 cm-63.5 cm) without back

21340B - 18.5" (47 cm)

21340BV - 18.5" (47 cm)

vacuum - formed seat

21340 - 16" (40.6 cm)

seamless upholstery
not standard

21340V - 18.5" (47 cm)

Vacuum-formed seat



*Model 21340
in Black (PR-40)*



*Model 21340B
in Wedgewood
(DAV-002)*

Accessories

2102265 - set of 5 pressurized locking casters (sit to lock casters)

3102277 - replacement casters (set of 5)

B - Backrest V - Vacuum-formed/Seamless

*Surgeon Model 21340 in
Black (PR-40)*

Medical Seating

An ideal blend of quality, comfort and value.

Brewer's ability to work quality craftsmanship and comfort-driven features into affordable seating solutions has made our Century Series and Value Plus Series popular choices at medical facilities around the world.

Century Series

Offers the comfort of a 16" (40.6 cm) wide, 4" (10.2 cm) thick poly-foam seat cushion with the dependability of a five-leg, 23" (58.4 cm) brushed aluminum accent base. Features a 5-year warranty and a 250-lb. weight capacity.

21"-28"

(53.3 cm - 71.1 cm)

Adjustable Height

Models with footring

11001BFR

11001BVFR

11001FR

11001VFR



Model 11001BFR in Cabernet (PR-36) Model 11001FR in Cabernet (PR-36)

18"-24"

(45.7 cm - 61 cm)

Adjustable Height

Models

11001B

11001BV

11001BLCD

11001BG*

11001B-UL

11001

11001V

11001LCD

11001G*

11001-UL



Model 11001B in Black (PR-40) Model 11001 in Black (PR-40)



Model 11001BFR in Cabernet (PR-36)

B - Backrest FR - Footring V - Vacuum-formed/Seamless
LCD - Locking casters G - Glides UL - Ultraleather™ upholstery

* Height range 17"-23" (43.2 cm-58.4 cm)



*Model 11001 in Black Satin
(PR-40)*

Value Plus Series

Designed for affordability, each features a 15" (38.1 cm) wide, 4" (10.2 cm) thick poly-foam seat cushion and five-leg, 21.5" (54.6 cm) black ABS plastic composite base.

17.5" - 22.5"
(44.5 - 57.2 cm)

**Pneumatic Height
Adjustment**

Models

- 22500B
- 22500BV
- 22500
- 22500V



*Model 22500B in Ash
(DUR-002)*

*Model 22500 in Ash
(DUR-002)*

17" - 22"
(43.2 - 55.9 cm)

**Spinlift Height
Adjustment**

Models

- 22400B
- 22400BV
- 22400
- 22400V



*Model 22400B in Ash
(DUR-002)*

*Model 22400 in Ash
(DUR-002)*

B - Backrest FR - Footring V - Vacuum-formed/Seamless
LCD - Locking casters G - Glides UL - Ultra leather™ upholstery

Heights are measured from top of seat cushion to floor.

Because our seating is custom built, heights and widths could vary within 1/2 of an inch of product listing.

Task Lighting & Medical Seating

Reliable solutions for a wide range of needs.

Round out your clinical care needs with lasting Brewer values. It's a choice your staff and patients will appreciate.



LED ECO-Series Task Lighting

This eco-friendly LED light consumes significantly less energy than halogens and its 50,000-hour bulb life is twice that of other LED exam lights available. Efficient power consumption also means it remains cooler to the touch.

Yet it provides an industry-best, bright white light (more than 20,000 lux / 880 Lumens with a 5" diameter beam at 18") to drive clinical efficiency. Color temperature is 4000° Kelvin.

Conveniently position that industry-best, bright white light precisely where you need it with the flexible, "no-drift" gooseneck arm. Features a 1-year warranty.

Model 19100

Accessories

2102755 - Replacement bulb
(50,000-hour bulb life)



Side Chair Series

Perfect for exam room, front office or waiting room environments, these durable chairs feature 250-lb. weight capacity steel tube frames and kick-out rear leg design to help prevent wall dings and dents. Their contoured backrests and 2" thick poly-foam seat cushions offer plenty of comfort. Features a 5-year warranty.

Models

1200 - w/ sturdy plastic arm rests

1250 - w/o arm rests



Model 1200 in Ivy



Model 1250 in Saddle

Blood Drawing Chair Series

Seat alignment positions patients securely on the chair, while innovative angled armrests and easy-to-use knob height adjustments provide optimal arm position for each patient's blood draw. Features include two standard add-on work surfaces (to accommodate prep materials, paperwork and vials) and a 350-lb. weight capacity. (Model # 101059) can be attached on either side of chair to keep prep materials nearby yet out of sight. Features a 5-year warranty.

Model 1500

Accessories

101059 - Drawer option



Model 1500 in Cabernet

Support Products

A single convenient source for so many quality solutions.



Square Steel Waste Can
32 Quart (8 Gallon) Models
12" D x 12" W x 20" H
35266 - White enamel finish
35267 - Red enamel finish
35268 - Beige enamel finish

Rectangular Plastic Waste Can
52 Quart (13 Gallon) Models
12" D x 15" W x 26" H
25266 - White plastic
25267 - Red plastic
25268 - Beige plastic

20 Quart (5 Gallon) Models
12" D x 12" W x 16" H
35269 - White enamel finish
35270 - Red enamel finish
35271 - Beige enamel finish

32 Quart (8 Gallon) Models
12" D x 15" W x 19" H
25269 - White plastic
25270 - Red plastic
25271 - Beige plastic



Visit brewercompany.com

Support Products

Ergonomic design solutions.

Attention-to-detail ergonomic design and quality construction are hallmarks of every product Brewer offers, whether it's a power table or a hamper stand.

NEW 5-Drawer Mobile Medical Cart

Premium construction, security and ease-of-use design make this the perfect addition to your practice.

Features a universal key lock to secure all 5 drawers at once. Offers 5.57 ft³ of storage. Premium, full-swivel 5" diameter casters with parking locks and an ergonomic push handle assure exceptional maneuverability and secure placement. Includes a seamless, easy-clean plastic work surface, a strong steel frame and base. Plus, soft-close return glides offer smooth, silent drawer closure and help prevent finger pinching. Drawer Dimensions: 15.5" L x 23" W; 2 drawers 3" H, 2 drawers 6" H and 1 Drawer 9" H.

Model 3500101



Stainless Steel All Purpose Cart

Features non-corrosive stainless steel construction & three 15.5" W x 24" L (39.4 cm x 61 cm) shelves.

Total weight capacity:

300 lbs (136 kg).

Dimensions with handle: 27" L x 16" W x 32" H (68.6 cm x 40.6 cm x 81.3 cm).

Model 63500



5-Drawer Mobile Medical Cart

Plastic Utility Cart

Includes three plastic shelves, 16" W x 21.75" L (40.6 cm x 55.3 cm) in light gray. Total weight capacity: 200 lbs (90.7 kg). Dimensions with handle: 28" L x 17.5" W x 35" H (71.1 cm x 44.5 cm x 94 cm).

Model 63530



Multi-Purpose Cart

Strong steel frame with durable white powder coat finish. Total weight capacity: 50 lbs (22.7 kg). Dimensions with handle: 37.5" H x 19.5" W x 29" D (95.3 cm x 49.5 cm x 73.6 cm). Baskets sold separately.

Models

- 63400 - Cart frame only
- 63021 - 12" White powder coat wire basket
- 63020-1 - 6" White powder coat wire basket (1 per carton)
- 63020-2 - 6" white powder coat wire basket (2 per carton)



Adjustable Folding Linen Hamper

Holds 18", 22" or 25" diameter bags (45.7 cm, 55.9 cm or 63.5 cm). Features durable white powder-coat finish. Folds for compact storage. Dimensions: 37" H x 20" x 20" D (93.9 cm x 50.8 cm x 50.8 cm).

Model 33395

Accessories

- 31331 - 18" Rope style 10 oz. cotton bag
- 31341 - 22"-25" Rope style 10 oz. cotton bag



18" Tilt-Top Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Step-on pedal for hands-free use. Dimensions: 35.75" H x 18" outside diameter (90.6 cm x 45.7 cm).

Model 33334

Accessories

- 31331 - 18" Rope style 10 oz. cotton bag
- 34000 - Bag holder



Tilt-Top Square Hamper

Open design for easy unloading. Step-on pedal for hands-free use. White powder-coat finish. Built-in bag holder supports heavy-duty bags. Dimensions: 38" H x 18.5" W x 19.5" D (91.4 cm x 47 cm x 49.5 cm).

Model 11410

Accessories

- 31331 - 18" Rope style 10 oz. cotton bag



18" Round Linen Hamper (45.7 cm)

Durable steel construction with white powder-coat finish. Dimensions: 35.75" H x 18.25" outside diameter (90.8 cm x 46.3 cm).

Model 33330

Accessories

- 31331 - 18" Rope style 10 oz. cotton bag
- 34000 - Bag holder



Support Products

More quality solutions for your care needs.



Step Stool

Features rugged, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety.

Dimensions: 9" H x 11.25" W x 14.25" L (22.8 cm x 28.6 cm x 36.2 cm). With handrail, total height is 35.25" (89.5 cm).

Weight capacity: 350 lbs (158 kg).

Models

- 11200-1 - 1 Stool per box
- 11220 - w/Hand rail

Accessories

- 99121 - Replacement crutch tips (set of 4)



Two-Step, Step Stool

Features heavy-duty, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety.

Dimensions (1st step): 9" H x 12" W x 16" L (22.9 cm x 30.5 cm x 40.6 cm).

Dimensions (2nd step): 7" H x 12" W x 16" L (17.8 cm x 30.5 cm x 40.6 cm). With handrail, total height is 44" (111.8 cm).

Weight capacity: 600 lbs (272 kg).

Models

- 31200 - Step stool
- 31220 - Step stool w/hand rail

Accessories

- 50016-1 - Hand rail
- 99121 - Replacement Crutch Tips (set of 4)



Bariatric Step Stool

Features heavy-duty, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety.

Dimensions: 9.25" H x 11.25" W x 14.25" L (23.5 cm x 28.6 cm x 36.2 cm). With handrail, total height is 34" (86.4 cm).

Weight capacity: 600 lbs (272 kg).

Models

- 21200 - Step stool
- 21220 - Step stool w/hand rail

Accessories

- 99121 - Replacement crutch tips (set of 4)
- 50015-1 - Hand rail



Mobile Instrument Stand

Adjustable height from 29"-47.5" (73.6 cm-120.7 cm) for comfortable use sitting or standing. Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 43465

Accessories

2102568 - Replacement tray




Stainless Steel Mayo Stand

Features non-corrosive stainless steel construction. Adjustable height from 31"-48.5" (78.7 cm-123.2 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 48461

Accessories

2102568 - Replacement tray



Mayo Instrument Stand

U-shaped base minimizes interference during procedures. Adjustable height from 31"-48.5" (78.7 cm-123.2 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 43460

Accessories

2102568 - Replacement tray




"U" Base Table

"U" base fits easily in restricted spaces. Offers convenient, spring-loaded height adjustment from 26.75"-43.25" (68 cm-109.9 cm). Top measures 15" x 30" (38.1 cm x 76.2 cm).

Models

11620 - w/Walnut grain laminate top & chrome base
11630 - w/Gray laminate top & gray epoxy finish on base



Mayo Instrument Stand

"California" style base for easy positioning. Adjustable height from 34"-52" (86.4 cm-132.1 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 11100

Accessories

2102568 - Replacement tray



Kick Bucket w/Frame

Features non-corrosive stainless steel construction. Heavy-duty design includes 360° rubber bumper on frame 13 quart size (3.25 gallon).

Model 48711

48712 - Kick bucket only
48714 - Kick bucket frame only



Low "H" Base Table

Designed for low care beds. Features spring-loaded height adjustment. With casters: 19.5"-26.75" (49.5 cm-67.9 cm). With glides: 17.625"-24.5" (44.8 cm-62.2 cm). Top measures 15" x 30" (38.1 cm x 76.2 cm).

Model

11640 - w/Walnut grain laminate top & chrome base



"H" Base Table

Offers convenient, spring-loaded height adjustment from 26.75"-43.25" (67.9 cm-109.9 cm). Top measures 15" x 30" (38.1 cm x 76.2 cm).

Models

11600 - w/Gray laminate top & gray epoxy finish on base
11610 - w/Walnut grain laminate top & chrome base




Support Products

All support products are backed by a standard 5-year warranty.

Large Infusion Pump Stand

Our 26" diameter (66 cm) design delivers the most stable base available, handling up to four infusion pumps and/ or IV controllers (w/ accessories 43433 & 43430). Exclusive "No Loss" knob allows easy height adjustment from 54"-96.5" (137.2 cm-245.1 cm).



Models

43408 - 2 Ram's
horn hook

43409 - 4 Ram's
horn hook

Short Wheel Infusion Pump Stand

Stable, space-saving design features heavy center weight with only a 16" diameter (40.6 cm and exclusive "No Loss" knob for easy height adjustment from 54"-96.5" (137.2 cm-245.1 cm).



Models

43416 - 2 Ram's
horn hook


43417 - 4 Ram's
horn hook

Accessories for 43416, 43417, 43408, 43409

The following accessories all attach with the #43433 clamp, which is sold separately.

Universal Clamp


Offers fast, flexible attachment & positioning of accessories with set-screw locks.



Model 43433

Add-A-Tray

Removable stainless steel tray measures 9.75" x 13.5" x 0.75" (48.9 cm x 32.1 cm x 1.9 cm).



Model 43431

Heavy Base Four-Leg IV Pole

Features tough, four-leg epoxy-coated steel base with 22" (55.9 cm) diameter, wall-protecting rubber bumpers on the legs, and exclusive "No Loss" knob for easy height adjustment from 51.75"-93" (131.5 cm-236.2 cm).



Models

43403 - 2 Ram's
horn hook

43406 - 4 Ram's
horn hook

Aluminum Base Five-Leg IV Pole

Includes five-leg cast aluminum base with 25" (63.5 cm) diameter, wall-protecting rubber bumpers on the legs, and exclusive "No Loss" knob for easy height adjustment from 52.25"-92.75" (132.7 cm-235.6 cm).



Models

11350 - 2 Ram's
horn hook

11360 - 4 Ram's
horn hook

Two-Hook Four-Leg IV Pole

Offers four-leg, epoxy-coated steel base with 22" (55.9 cm) diameter and easy, twist-lock height adjustment from 47"-92" (119.4 cm-233.6 cm).




Model

11300 - 2 Ram's
horn hook

Upholstery Color Options

Standard upholstery color selections.

For Brewer Tables, Side Chairs, Blood Drawing Chairs and Exam Stools.



Clamshell (21) - PR 57



Gunmetal (22) - PR 92
Vacuum Upholstery -
DAV 016



Feather (23) - PR 37



Black Satin (24) - PR 40



Saddle (25) - PR 52



Tapestry Red (26) -
PR 47



Cabernet (27) - PR 36



Cocoa (28) - PR 51



Deep Sea (29) - PR 49



Ivy (30) - PR 53



Deep Sapphire (31) -
PR 55



Blue Fog (32) - PR 46

Dozens more upholstery options available, contact Brewer customer service or go to brewercompany.com

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1

Clinical & Exam Room Solutions



Access® High-Low Exam Table FLEX™ Access High-Low Exam Table Assist Power Procedure Table® Versa Fixed-Height Exam Table

Our process-oriented solutions are designed to help you improve operational efficiencies while more effectively protecting patients, staff and your practice. From exclusive safety advances for confident transfers to features like our large, signature pass-through drawers and work surfaces. Brewer helps you meet the challenges of today's healthcare reality.

What's more, our tables and stools are proudly made in the USA. Our extensive line of solutions is manufactured to Brewer's rigorous quality standards and backed by best-in-class standard warranties. It's all part of our commitment to help healers heal.

For more information about the many ways we help healers heal, contact your Brewer representative today.

Brewer

The Brewer Company, LLC
N88 W13901 Main Street, Suite 100
Menomonee Falls, WI 53051

1.888.273.9371
www.brewercompany.com



We Help Healers Heal.