

Assist[®] Power Procedure Table

Specifications

Advanced Safety

- ADA-compliant 19" seat height helps enable safe lateral wheelchair transfers.
- Return to chair feature on power back models
- Anti-tamper on/off switch limits accidental motion
- 19" (48.3 cm) low height that raises to 40" (101.6 cm), both ADA compliant
- Easy on/off for patients

Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density, easy-to-clean seamless upholstery
- Standard adjustable ergonomic stirrups
- Tilt range of up to 30 degrees allows true Trendelenburg positioning.

Clinician Convenience

- Enables quick transitioning from one type of procedure to the next for increased throughput

Table Specifications

- 450 lb (204 kg) weight capacity
- Ergonomic adjustable stirrups come standard
- Foot control comes standard
- Industry best 3-year warranty
- **Height:** 19"-40" (48.3 cm – 101.6 cm)
— Both heights are ADA compliant
- **Length:** 69", extends to 80" with headrest and leg rest (175.3 cm – 203.2 cm)
- **Length in upright position:** 40" (101.6 cm)
- **Upholstery Width:** 28" (71.1 cm)
- **Leg rest:** 16.5" x 11.25" (41.9 cm x 28.6 cm)
- **Headrest:** 13.75" x 11.75" (34.9 cm x 29.8 cm)
- **Base:** 24"W x 36"L (61 cm x 91.4 cm)
- 115VAC, 60Hz
- Safe, low voltage actuators
- **Slide-through paper roll holder** accepts 21" & 18" rolls
- Tilt range/Trendelenburg: to 30 degrees
- Easy to clean, seamless upholstery
- Hand pendant standard on AssistPRO (7500 model)

Certifications

All Brewer products conform to US & Canadian electrical & medical industry codes & regulations and meet or exceed industry requirements.

- UL 60601-1
- CAN-CSA C22.2 No. 601.1
- ISO 13485:2016



Models:

Recommended Models:

- **7000**
Brewer Assist Power Procedure Table
- **7500**
Brewer AssistPRO Power Procedure Table

Additional Specifications for AssistPRO Model 7500:

- Intuitive hand pendant and foot control included
- Four preferred positions can be programmed into the table's memory
- Auto-run toggle for safety

Accessories:

- **99497**
Grab bars (qty 2)*
- **99499:**
Articulating armboard*
- **99500:**
Hand pendant kit for Assist 7000
- **99501:**
Side seat rails (qty 2)
- **99502:**
Universal socket clamp (qty 1)
- **99504:**
Deluxe articulating knee crutches (left and right) **
- **100253:**
Urology drain pan system
- **102024:**
Welch Allyn Bracket (Factory Install)
- **102025:**
Welch Allyn Bracket (Field Install)

* Needs one order of #99501

** Needs one order of #99501 and two orders of #99502

Create an elevated patient-provider experience.

Models	7000	7500
Weight Capacity	450	450
ADA Compliant Seat Height	☑	☑
Adjustable Headrest	☑	☑
Tilt Range up to 30°	☑	☑
Seamless Upholstery	☑	☑
Ergonomic Stirrups	☑	☑
Standard 3-Year Warranty	☑	☑
Programmable Positions		☑
Hand Pendant Included		☑



Trendelenburg Position



Lithotomy Position



Proctology Position

Available Colors.



Clamshell
21



Gunmetal
22



Feather
23



Black Satin
24



Saddle
25



Tapestry
26



Cabernet
27



Cocoa
28



Deep Sea
29



Ivy
30



Deep Sapphire
31



Blue Fog
32

Need a closer look? Contact us to receive physical samples of our vinyl or base materials.
Custom colors available.



BREWER

1-888-273-9371 | BrewerCompany.com

The Brewer Difference.

Industry’s most safety features. | Industry’s only 5-year warranty. | USA-made in the Midwest.