

FLEX™ High-Low PLUS Power Exam Table with 17" Low Seat Height



Specifications

Advanced Safety

- Exclusive SafeGlide™ leg extension increases safety and satisfaction
- ComfortGrip™ Patient Assist Handles facilitate easy transfers and adjustments, enhancing patient independence and safety during examinations (*optional*)
- Return to chair feature on power back models
- Two safety zones offer protection against pinching and impact-related injuries.
- 17" (43.2 cm) low height that raises to 37" (94 cm), both ADA compliant
- Easy on/off for patients
- Powered backrest adjusts for easy patient positioning*

* Pneumatic back models are also available

Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density support upholstery
- Adjustable ergonomic stirrups (*optional*)
- Toe/leg extension & plastic debris tray

Clinician Convenience

- **Front drawer:**
19.75" x 19.25" x 4" (50.2 cm x 48.9 cm x 10.2 cm)
- **Paper roll holder accepts:**
18" (45.7 cm) & 21" (53.3 cm) rolls

Table Specifications

- 700 lb (317 kg) weight capacity
- ComfortGrip™ patient assist handles (*optional*)
- Ergonomic adjustable stirrups come standard
(*standard on 5701PLUS-17 and 5801PLUS-17 models; optional on 5700PLUS-17 and 5800PLUS-17*)
- Foot control comes standard
- Standard 5-year warranty
- **Height:** 17"–37" (43.2 cm – 94 cm) — Both heights meet the 2024 Access Board recommendation.
- **Length:** 60", extends to 78" (152.4 cm – 198.1 cm)
- **Upholstery Width:** 27.5" (69.9 cm)
- **SafeGlide™ Leg extension length:** 18" (45.7 cm)
- 115VAC, 60Hz outlet
- Safe, low voltage actuators
- Up to 6° pelvic tilt
- 60-Watt front drawer heater
- Easy to clean, seamless upholstery

Certifications

All Brewer products conform to US & Canadian electrical & medical industry codes & regulations and meet or exceed industry requirements.

- UL 60601-1
- CAN-CSA C22.2 No. 601.1
- ISO 13485:2016



1-888-273-9371 | [BrewerCompany.com](https://www.BrewerCompany.com)

17" Seat Height Option:

- **5700PLUS-17:**
Pneumatic back & 17" low-transfer height
- **5701PLUS-17:**
Pneumatic back, stirrups, pelvic tilt, front drawer heater, outlet & 17" low-transfer height
- **5800PLUS-17:**
Power back & 17" low-transfer height
- **5801PLUS-17:**
Power back, stirrups, pelvic tilt, front drawer heater, outlet & 17" low-transfer height

Accessories:

- **3500687PLUS:**
EasyGlide™ Mobility System
 - **99505:**
Stainless Steel Treatment Pan
 - **99502:**
One Universal Socket Clamp
 - **99504:**
Deluxe Articulating Knee Crutches*
 - **100253-TG:**
Urology Drain Pan System
 - **100407PLUS:**
Front Drawer Dividers
 - **100540:**
Two Support Adapters for Knee Crutches (factory installed)
 - **100550:**
Two Support Adapters for Knee Crutches (field installed)
 - **102031PLUS:**
Brewer Welch Allyn® IV Bracket (factory installed)
 - **102032PLUS:**
Brewer Welch Allyn® IV Bracket (field installed)
- *Needs two 99502 and one of either 100540 or 100550

Optional Kits

- **2102894PLUS:**
Stirrups for 5700PLUS-17/5800PLUS-17 (set of two factory installed)
- **2102898PLUS**
Stirrups for 5700PLUS-17/5800PLUS-17 (set of two field installed)
- **2102891PLUS**
ComfortGrip™ Patient Assist Handles (set of two, factory installed)
- **2102895PLUS**
ComfortGrip™ Patient Assist Handles (set of two, field installed)

Create an elevated patient-provider experience.

Models	5700PLUS-17	5701PLUS-17	5800PLUS-17	5801PLUS-17
Weight Capacity	700	700	700	700
Access Board 17" Seat Height Recommendation	✓	✓	✓	✓
Power Base	✓	✓	✓	✓
Power Back			✓	✓
Manual Back	✓	✓		
Patient Assist Handles	Optional	Optional	Optional	Optional
Stirrups	Optional	✓	Optional	✓
Pelvic Tilt & Drawer Warmer		✓		✓
Outlet		✓		✓
EasyGlide Mobility	Optional	Optional	Optional	Optional
Mobility	Optional	Optional	Optional	Optional

Available Colors.

 Chateau Grey AF906	 Flannel AF945	 Horizon AF928	 Mist AF938	 Black Satin PR40
 Cool Grey MV821	 Eclipse AF919	 Blue Fog PR46	 Cocoa PR51	 Cabernet PR36
 Feather PR37	 Deep Sapphire PR55	 Gunmetal PR92	 Clamshell PR57	 Ivy PR53

Need a closer look? Contact us to receive physical samples of our vinyl or base materials.

Custom colors available.



1-888-273-9371 | [BrewerCompany.com](https://www.BrewerCompany.com)

The Brewer Difference.

Industry's most
safety features.

Industry's only
5-year warranty.

USA-made in
the Midwest.

LIT-M-074 RevB