



Clinical & Exam Room Solutions

Choose different, Choose better.

From patented safety features to whisper-quiet mobility and intuitive ergonomics, every Brewer product is designed to support your team and protect your patients. Backed by the industry's only 5-year warranty and proudly made in the USA, we help you deliver better care with tools you can trust.

Industry's most
safety features.

Industry's only
5-year warranty.

USA-made in
the Midwest.

FULLY ADA COMPLIANT

Every Brewer power exam table meets current ADA regulations, including the required 17"-19" seat height range, to support safe, dignified patient transfers. From ComfortGrip™ Assist Handles to wheelchair-accessible design, our tables make accessibility standard.



POSITIONING MADE EASY

From fingertip leg extension adjustments to supportive assist handles and high weight capacities, every Brewer exam table is designed for safer, more comfortable positioning. Whether fixed-height or powered, our exam tables help patients feel secure and help caregivers work with confidence.

The industry's only one-touch, silent mobility system has arrived.



A simple turn of a switch engages
and disengages the EasyGlide
Mobility System.

INTRODUCING THE BREWER EASYGLIDE™ MOBILITY SYSTEM.

The only one to offer a one-touch, whisper-quiet mobility system, Brewer makes it possible to reconfigure and clean exam rooms quickly and with zero disruption to patients and providers. Our advanced EasyGlide Mobility System is an available option on our Access® PLUS and FLEX™ PLUS Exam Tables.



See our EasyGlide Mobility System in action.

Safer products are just the beginning.

Healthcare professionals want the safest environment possible for their patients and staff. Brewer understands just how important the team behind the products can be to making that goal a reality.



EXPERT GUIDANCE

Our knowledgeable product consultants are here to help you find the right exam table for your space, staff, and patient needs. Whether you're focused on accessibility, ergonomics, or budget, we'll help guide you to the best-fit solution.



ON-SITE TRAINING

We provide hands-on support to ensure your team understands how to safely and effectively use every feature built into your Brewer equipment.



RESPONSIVE SUPPORT

Our relationship doesn't end after delivery. With the industry's only 5-year warranty and a dedicated support team, we're here to ensure your equipment continues to perform day after day, year after year.

Table of Contents

01

POWER EXAM TABLES

Access® High-Low PLUS Power Exam Table	8
FLEX™ High-Low PLUS Power Exam Table	10

02

POWER PROCEDURE TABLE

Assist™ Power Procedure Table	12
-------------------------------------	----

03

FIXED-HEIGHT EXAM & TREATMENT TABLES

Access® 1000 Medical Exam Table	14
Access® Medical Exam Table	16
Element Treatment Table	18

04

MEDICAL SEATING & CLINICAL SUPPORT PRODUCTS

Millenium Series	21
Century Series	24
Value Plus Series	25
Side Chairs	26
Blood Drawing Chair	27
135AS and 135JS Saddle Stool	28
Eco Series LED Exam Light	29

05

ACCESSORIES

Accessories	30
-------------------	----

01



5-YEAR
WARRANTY



POWER EXAM TABLES

Access® High-Low PLUS
Power Exam Table



The one that has it all.

Access[®] High-Low PLUS Power Exam Table



450 LB MODELS

6000PLUS - Pneumatic back with stirrups, assist handles & work surface

6001PLUS - Pneumatic back with stirrups, assist handles, work surface, pelvic tilt & drawer warmer

6500PLUS - Power back with stirrups, assist handles & work surface

6501PLUS - Power back with outlet, stirrups, assist handles, work surface, pelvic tilt & drawer warmer

700 LB MODELS

6800PLUS - Power back with stirrups, assist handles & work surface

6801PLUS - Power back with outlet, stirrups, assist handles, work surface, pelvic tilt & drawer warmer

01



5-YEAR
WARRANTY



POWER EXAM TABLES

FLEX™ High-Low PLUS

Power Exam Table



Power exam table tailored to your practices needs.

FLEX™ High-Low PLUS Power Exam Table



700 LB MODELS

5700PLUS - Pneumatic back

5701PLUS - Pneumatic back, stirrups, pelvic tilt, front drawer heater and outlet

5800PLUS - Power back

5801PLUS - Power back, stirrups, pelvic tilt, front drawer heater and outlet

02



POWER PROCEDURE TABLES

Assist™ Power Procedure Table



Boost your profitability with Brewer flexibility.

Assist™ Power Procedure Table

MODEL

7000 - Assist Power Procedure Table

7500 - AssistPRO 4-Position Programmable Power
Procedure Table



TRENDELENBURG
POSITION



ERGONOMIC STIRRUPS



LITHOTOMY POSITION



PROCTOLOGY POSITION

Scan the QR code to view full product specs and available accessories.

03



FIXED-HEIGHT EXAM & TREATMENT TABLES

Access® 1000 Medical
Exam Table



One exam table for your entire patient population.

Access[®] 1000 Medical Exam Table



ADULT MODELS

1000A - Pneumatic back, stirrups and outlet

1000A-NR - Pneumatic back and stirrups

1001A - Pneumatic back, stirrups, pelvic tilt, outlet and drawer warmer

1001A-NR - Pneumatic back, stirrups, pelvic tilt, and drawer warmer

PEDIATRIC MODEL

1000P - Cushioned pediatric top

*The table shown is fully equipped with optional accessories.
Scan the QR code to view full product specs and available accessories.*

03



FIXED-HEIGHT EXAM & TREATMENT TABLES

Access[®] Medical Exam Table



Traditional designs infused with Brewer's value-added touches.

Access[®] Medical Exam Table



MODELS

5000 - Pneumatic back, stirrups and outlet

5000-SP-CN - Pneumatic back and stirrups

5001 - Pneumatic back, stirrups, outlet, pelvic tilt and drawer heater

5001-SP-CN - Pneumatic back, stirrups, pelvic tilt and drawer heater

03



5-YEAR
WARRANTY



FIXED-HEIGHT EXAM & TREATMENT TABLES

Element Treatment Table



Fundamental features for a growing number of healthcare practices.

Element Treatment Table






2110



2040

MODELS

	flat top	flat top adjustable	pillow top adjustable
basic frame	 #2010	 #2030	 #2040
frame w/ shelf	 #2110	 #2130	 #2140
frame w/ shelf & drawers	 #2210	 #2230	 #2240

04



MEDICAL SEATING & CLINICAL SUPPORT PRODUCTS



Premium specialty seating with an ergonomic design.

Millennium Series Lab Stool

PART NUMBER	HEIGHT RANGE	CASTERS / GLIDES	FOOTRING	BACKREST	SEAMLESS	STND. COLORS
21520B	21"-31"	Glides	Yes	Yes	Yes	15
21521B	22"-32"	Casters	Yes	Yes	Yes	15



21520B



21521B

Hands-free height adjustable for greater flexibility.

Millennium Series Surgeon Stool

POPULAR MODELS	HEIGHT RANGE	FOOT OPERATED	BACKREST	SEAMLESS	STND. COLORS
21340	20"-25"	Yes	No	No	15
21340B	18.5"-23.5"	Yes	Yes	Yes	15



21340



21340B

The most sophisticated ergonomic design in its class.

Millennium Series Task Stool

POPULAR MODELS	HEIGHT RANGE	ARM RESTS	FOOTRING	BACKREST	SEAMLESS	STND. COLORS
21435B	18"-24"	No	No	Yes	Yes	15
21435BA	18"-24"	Yes	No	Yes	Yes	15



21435B



21435BA

Preferred by medical facilities
worldwide.

Century Series

POPULAR MODELS	HEIGHT RANGE	FOOTRING	BACKREST	SEAMLESS	STND. COLORS
11001	18"-24"	No	No	No	15
11001FR	21"-28"	Yes	No	No	15
11001B	18"-24"	No	Yes	No	15
11001BFR	21"-28"	Yes	Yes	No	15



11001



11001B



11001BFR

Affordable without compromising style and reliability.

Value Plus Series

POPULAR MODELS	HEIGHT RANGE	CASTERS / GLIDES	BACKREST	SEAMLESS	STND. COLORS
22500	17.5"-22.5"	Casters	No	No	15
22500B	17.5"-22.5"	Casters	Yes	No	15
22400	17"-22"	Casters	No	No	15
22400B	17"-22"	Casters	Yes	No	15



22500



22500B

Perfect for the exam room, front office or waiting room environment.

Side Chairs

MODELS

1200 – Brewer Side Chair With Arms

1250 – Brewer Side Chair Without Arms



Accommodates prep materials, paperwork, and vials.

Blood Drawing Chair

MODEL

1500 - Brewer Blood Drawing Chair

101059 - Drawer Accessory



Dynamically different for comfort and performance.

135AS AND 135JS Saddle Stool

PART NUMBER	ADJUSTABLE HEIGHT RANGE	RECOMMENDED USER HEIGHT
135AS	24" - 32"	5' 3" - 5' 11"
135ASB-LS	24" - 32"	5' 3" - 5' 11"
135AS-L	21" - 27"	Under 5' 3"

PART NUMBER	ADJUSTABLE HEIGHT RANGE	RECOMMENDED USER HEIGHT
135JS	24" - 32"	5' 3" - 5' 11"
135JSB-LS	24" - 32"	5' 3" - 5' 11"
135JS-H	27" - 37"	Over 5' 11"



135ASP



135ASB-LS



135JSG

**AVAILABLE
UPHOLSTERY COLORS:**

- G- DOVE GRAY
- B- RAVEN WING
- P- PAPYRUS


Best-in-class illumination and a long bulb life.

ECO Series LED Exam Light

MODELS

19100 - ECO-Series LED Light - Beige

19100-230 - Eco-Series LED Light-230 Volt
*(only sold in USA)



05

ACCESSORIES

Access® High-Low PLUS Power Exam Table

HLTPLUS17-STD: 17IN Low Transfer Upholstery Top

6115PLUS17-STD: 17IN Low Transfer Upholstery Top Service Kit

3500687PLUS: EasyGlide™ Mobility System

99502: One Universal Socket Clamp

99504: Deluxe Articulating Knee Crutches*

100253-TG: Urology Drain Pan System

100407PLUS: Front Drawer Dividers

3500077: Footrest Pad

100540: Two Support Adapters for Knee Crutches (factory installed)

100550: Two Support Adapters for Knee Crutches (field installed)

102031PLUS: Brewer Welch Allyn® IV Bracket (factory installed)

102032PLUS: Brewer Welch Allyn® IV Bracket (field installed)

**Needs two 99502 and one of either 100540 or 100550*

FLEX™ High-Low PLUS Power Exam Table

FLEXPLUS17-STD: 17IN Low Transfer Upholstery Top

5715PLUS17-STD: 17IN Low Transfer Upholstery Top Service Kit

3500687PLUS: EasyGlide™ Mobility System

99502: One Universal Socket Clamp

99504: Deluxe Articulating Knee Crutches*

100253-TG: Urology Drain Pan System

100407PLUS: Front Drawer Dividers

3500077: Footrest Pad

100540: Two Support Adapters for Knee Crutches (factory installed)

100550: Two Support Adapters for Knee Crutches (field installed)

102031PLUS: Brewer Welch Allyn® IV Bracket (factory installed)

102032PLUS: Brewer Welch Allyn® IV Bracket (field installed)

**Needs two 99502 and one of either 100540 or 100550*

Assist™ Power Procedure Table

99497: Grab Bars (Qty. 2) **

99499: Articulating Arm Board **

99500: Hand Pendant

99501: Side Seat Rails (Qty. 2)

99502: Universal Socket Clamp (Qty. 1)

99504: Deluxe Articulating Knee Crutches (Pair – Right & Left) ***

100253: Urology Drain Pan Kit

102024: Brewer Welch Allyn® Exam Light IV Bracket (plant installed)

102025: Brewer Welch Allyn® Exam Light IV Bracket (field installed)

*** Requires two side seat rails (one order of 99501)*

**** Requires two side seat rails (one order of 99501) and two universal socket clamps (two orders of 99502)*

Access® 1000 Medical Exam Table

3500929: Assist Handles (factory installed)

98464: Stainless Steel Treatment Pan

100253-TG: Urology Drain Pan System

3500288: Front Drawer Dividers

3500447: Brewer Welch Allyn® Light IV Bracket (plant installed)

3500398: Brewer Welch Allyn® Light IV Bracket (field installed)

3500286: Deluxe articulating knee crutches with assembly kit

3500305: Paper Fastener Strap

3500077: Footrest Pad

Access® Medical Exam Table

98464: Stainless Steel Treatment Pan

100253: Urology Drain Pan System

100407: Front Drawer Dividers (can order up to 7)

102022: Brewer Welch Allyn® Light IV Bracket (plant installed)

102023: Brewer Welch Allyn® Light IV Bracket (field installed)

3500077: Footrest Pad

Element Treatment Table

10196: Table Paper Straps

100407: Drawer Dividers (Models 2210, 2230 & 2240 only)

Medical Seating

2102265: Set of 5 pressurized locking casters

101059: Blood Drawing Chair Drawer

For Healers. For Life.

Brewer's exam tables and stools are thoughtfully designed to help clinicians work smarter, safer, and with greater confidence. From exclusive safety features that support more secure transfers to workflow-enhancing touches like our spacious pass-through drawers and integrated work surfaces, every detail reflects our commitment to today's healthcare realities.

All of our solutions are manufactured in the USA, backed by best-in-class warranties, and built to Brewer's rigorous quality standards. Because when equipment works better, caregivers can focus on what matters most: delivering exceptional care.

We Help Healers Heal.

To learn more, contact your Brewer representative today.



The Brewer Company, LLC
N88 W13901 Main Street, Suite 100
Menomonee Falls, WI 53051
1.888.273.9371
www.brewercompany.com