

Brewer

DENTAL SEATING

Ergonomically designed to support the health of your practice



CATALOG



Redefining comfort with the latest and greatest ergonomic advances.

From carefully researched design refinements to entirely new concepts in seating, Brewer is constantly working to offer you the most advanced and comfortable seating solutions in the industry. All models feature a 250-lb weight capacity and a standard 5-year warranty.

- 4 SADDLE SERIES**
Amazone + Jumper Models
- 6 9100 SERIES**
Dentist/Hygienist + Assistant Models
- 7 9500 SERIES**
Dentist/Hygienist + Assistant Models
- 8 3300 SERIES**
Dentist/Hygienist + Assistant Models
- 9 3100 SERIES**
Dentist/Hygienist + Assistant Models
- 10 2000 SERIES**
Dentist/Hygienist + Assistant Models
- 11 11001 SERIES**
Dentist/Hygienist + Assistant Models
- 12 ROOM SEATING**
Side Chairs
- 13 SUPPORT PRODUCTS**
Task Lighting + Mobile Cart

TABLE OF CONTENTS

SADDLES

AMAZONE SERIES

Features narrower seat and choice of two height ranges specifically designed to accommodate smaller-framed dental professionals.

135 Amazone Models*	Adjustable Height Range	Recommended User Height
■ 135ASG	24" - 32"	5'3" - 5'11"
■ 135ASB	24" - 32"	5'3" - 5'11"
■ 135ASP	24" - 32"	5'3" - 5'11"
■ 135ASG-L	21" - 27"	Under 5'3"
■ 135ASB-L	21" - 27"	Under 5'3"
■ 135ASP-L	21" - 27"	Under 5'3"

JUMPER SERIES

Features wider seat and choice of two height ranges specifically designed to accommodate taller, larger-framed dental professionals.

135 Jumper Models*	Adjustable Height Range	Recommended User Height
■ 135JSG	24" - 32"	5'3" - 5'11"
■ 135JSB	24" - 32"	5'3" - 5'11"
■ 135JSP	24" - 32"	5'3" - 5'11"
■ 135JSG-H	27" - 37"	Over 5'11"
■ 135JSB-H	27" - 37"	Over 5'11"
■ 135JSP-H	27" - 37"	Over 5'11"

3 Ultraleather options to choose from

* G - Dove Gray B - Raven Wing P - Papyrus



Exclusive capability promotes back health

Brewer is the only U.S. manufacturer to offer 360° Dynamic Motion technology. This ergonomic capability is designed to engage core muscles and help strengthen back muscles over time. As a result, dentists and hygienists can stay more comfortable and productive throughout the day.

**360°
Dynamic
Motion**

*Model 135ASP Amazone
in Papyrus*



*Model 135JSB Jumper
in Raven Wing*



9100 SERIES

DENTIST / HYGIENIST MODEL

Features synchronized seat / back tilt and independent backrest adjustment at the touch of your fingertips to comfortably allow greater range of movement. Plus, all models feature our exclusive HybriGEL™ memory foam for unmatched pressure management and cooling comfort.

9100 Series	Adjustable Height Range	Seamless Upholstery
9100B	19" - 25"	-
9100BV	19" - 25"	Yes



Model 9100B in Burgundy

ASSISTANT MODELS

Features synchronized seat / back tilt and independent backrest adjustment, ratchet-adjustable contoured body support with HybriGEL™ memory foam and an adjustable foot ring.

9100 Series	Adjustable Height Range	Seamless Upholstery	Footring	Body Support
9120BL	24.5" - 32.5"	-	Yes	Left
9120BLV	24.5" - 32.5"	Yes	Yes	Left
9120BR	24.5" - 32.5"	-	Yes	Right
9120BRV	24.5" - 32.5"	Yes	Yes	Right



Model 9120BL in Walnut

Accessories

101282 - Optional stool arms (sold separately)



101282 Stool Arms

Cushion Sizing:

Backrest 17.5" W x 11.5" H, Seat 19.5" W x 15" D



9500 SERIES

DENTIST / HYGIENIST MODEL

Features synchronized seat / back tilt and independent backrest adjustment at the touch of your fingertips to comfortably allow greater range of movement. Plus, all models feature our exclusive HybriGEL™ memory foam for unmatched pressure management and cooling comfort.

9500 Series	Adjustable Height Range	Seamless Upholstery
9500B	19.5" - 25.5"	-
9500BV	19.5" - 25.5"	Yes



Model 9500B in Blue Fog

ASSISTANT MODELS

Features synchronized seat / back tilt and independent backrest adjustment, ratchet-adjustable contoured body support with HybriGEL™ memory foam and adjustable foot ring.

9500 Series	Adjustable Height Range	Seamless Upholstery	Footring	Body Support
9520BL	25" - 33"	-	Yes	Left
9520BLV	25" - 33"	Yes	Yes	Left
9520BR	25" - 33"	-	Yes	Right
9520BRV	25" - 33"	Yes	Yes	Right



Model 9520BL in Charcoal

Accessories

101282 - Optional stool arms (sold separately)



101282 Stool Arms

Cushion Sizing:

Backrest 17" W x 9.5" H, Seat 19.5" W x 15" D



Scan for video

3300 SERIES

DENTIST / HYGIENIST MODEL

Features synchronized seat / back tilt and comes equipped with 24" diameter, powder-coated, 5-leg cast-aluminum base and hooded, dual-wheel casters.

3300 Series	Adjustable Height Range	Arms
3335B	18" - 25"	-
3335BV	18" - 25"	-
3335BA	18" - 25"	Yes
3335BVA	18" - 25"	Yes



Model 3335B
in Shetland



Model 3335BA
in Fawn

ASSISTANT MODELS

Features 24" diameter, powder-coated, 5-leg cast-aluminum base and hooded, dual-wheel casters and patented ratchet body support.

3300 Series	Adjustable Height Range	Footring	Body Support	Backrest
3345L	23" - 31"	Yes	Left	-
3345LV	23" - 31"	Yes	Left	-
3345R	23" - 31"	Yes	Right	-
3345RV	23" - 31"	Yes	Right	-
3345BL	23" - 31"	Yes	Left	Yes
3345BLV	23" - 31"	Yes	Left	Yes
3345R	23" - 31"	Yes	Right	Yes
3345BRV	23" - 31"	Yes	Right	Yes



Model 3345BL
in Admiral

Cushion Sizing:

Backrest 13" W x 15" H, Seat 17" W x 19" D



Scan for video

3100 SERIES

DENTIST / HYGIENIST MODEL

Features contoured tilting backrest, 24" diameter, powder-coated, 5-leg cast-aluminum base and hooded, dual-wheel casters.

3100 Series	Adjustable Height Range	Backtilt
3125B	18" - 24"	Yes
3125BV	18" - 24"	Yes



Model 3125B
in Blue Fog

ASSISTANT MODELS

Features 24" diameter, powder-coated, 5-leg cast-aluminum base and hooded, dual-wheel casters and patented ratchet body support.

3100 Series	Adjustable Height Range	Footring	Body Support
3145L	23" - 31"	Yes	Left
3145LV	23" - 31"	Yes	Left
3145R	23" - 31"	Yes	Right
3145RV	23" - 31"	Yes	Right



Model 3145L
in Wedgewood



Model 3145R
in Gunmetal

Cushion Sizing:

Backrest 13" W x 15" H, Seat 17" W x 19" D



2000 SERIES

DENTIST / HYGIENIST MODEL

Contemporary and Traditional Series provides practices with cost-effective solutions without compromises in comfort or quality.

2000 Series	Adjustable Height Range	Backrest
2020	18.5" - 24.5"	-
2020B	18.5" - 24.5"	Yes



Model 2020
in Deep Sapphire

ASSISTANT MODELS

Contemporary and Traditional Series provides practices with cost-effective solutions without compromises in comfort or quality.

2000 Series	Adjustable Height Range	Ratcheted Body Support	Backrest	Footring	Body Support
2042L	22" - 28"	-	-	Yes	Left
2042R	22" - 28"	-	-	Yes	Right
2052L	22" - 28"	Yes	-	Yes	Left
2052R	22" - 28"	Yes	-	Yes	Right
2052BL	22" - 28"	Yes	Yes	Yes	Left
2052BR	22" - 28"	Yes	Yes	Yes	Right



Model 2020B
in Dove



Model 2052BL
in Feather



Model 2042L
in Slate Blue

Cushion Sizing:

Backrest 14" W x 10" H, Seat 16" W x 4" thick



Scan for video

11001 SERIES

DENTIST / HYGIENIST MODEL

Contemporary and Traditional Series provides practices with cost-effective solutions without compromising in comfort or quality.

11001 Series	Adjustable Height Range	Backrest
11001-D	18" - 24"	-
11020-D	21" - 28"	-
11001B-D	18" - 24"	Yes



Model 11001-D
in Fawn

ASSISTANT MODELS

Contemporary and Traditional Series provides practices with cost-effective solutions without compromises in comfort or quality.

11001 Series	Adjustable Height Range	Foot Ring	Swing Arm
11001FRA-D	21" - 28"	Yes	Yes
11001VFRA-D	21" - 28"	Yes	Yes



Model 11001B-D
in Velvet Brown



Model 11001FRA-D
in Stone

Cushion Sizing:

Backrest 14" W x 10" H, Seat 16" W x 4" thick



ROOM SEATING

SIDE CHAIR SERIES

Perfect for the exam room, front office or waiting room environment. These affordable and easy-to-clean chairs are excellent in any healthcare setting.

Models	Contoured Backrest	Arms
1200	Yes	Yes
1250	Yes	-



*Model 1250 in
Blue Fog*



*Model 1200 in
Black Satin*

Cushion Sizing:

Backrest 14" W x 10" H, Seat 16" W x 2" thick

SUPPORT PRODUCTS

Task Lighting

This eco-friendly LED light outperforms halogens, with significantly lower energy consumption and bulb life of 50,000 hours, twice that of other LED exam lights. Efficient power usage ensures it stays cool to the touch while providing an industry-best bright white light, exceeding 20,000 lux / 880 lumens with a 5" diameter beam at 18". Easily position the light with the flexible, "no-drift" gooseneck arm.

5-Drawer Mobile Medical Cart

Features a universal key lock to secure all five drawers simultaneously. Equipped with premium, full-swivel 5" diameter casters with parking locks and an ergonomic push handle for exceptional maneuverability and secure placement. Includes seamless, easy-to-clean plastic work surface, sturdy steel frame and base.

Dimensions: 40"H x 27"W x 18"D



Model 3500101

Accessories:

2102755 Replacement Bulb (50,000-hour bulb life)



Model 19100



Scan for video

Item Number	Height Range	Description
9100B-L	15.5" - 20"	Operator, Low Cylinder
9100BV-L	15.5" - 20"	Operator, Low Cylinder
9100B-H	20.5" - 26.5"	Operator, High Cylinder
9100BV-H	20.5" - 26.5"	Operator, High Cylinder
9100BFR	22" - 31"	Operator, Adjustable Foot Ring
9100BVFR	22" - 31"	Operator, Adjustable Foot Ring
9500B-L	15.5" - 20"	Operator, Low Cylinder
9500BV-L	15.5" - 20"	Operator, Low Cylinder
9500B-H	20.5" - 26.5"	Operator, High Cylinder
9500BV-H	20.5" - 26.5"	Operator, High Cylinder
9500BFR	22" - 31"	Operator, Adjustable Foot Ring
9500BVFR	22" - 31"	Operator, Adjustable Foot Ring
3335BFR	20" - 26"	Operator, Adjustable Foot Ring
3335BVFR	20" - 26"	Operator, Adjustable Foot Ring
3355BL	22" - 31"	Assistant, Backrest, Seat/Back Tilt
3355BLV	22" - 31"	Assistant, Backrest, Seat/Back Tilt
3355BR	22" - 31"	Assistant, Backrest, Seat/Back Tilt
3355BRV	22" - 31"	Assistant, Backrest, Seat/Back Tilt
3125BFR	20" - 26"	Operator, Back Tilt, Foot Ring
3125BFRV	20" - 26"	Operator, Back Tilt, Foot Ring
2020V	18" - 24"	Operator
2020BV	18" - 24"	Operator, Backrest
2020G	18" - 24"	Operator w/ Glides
2021	20" - 26"	Operator, High Height
2021V	20" - 26"	Operator, High Height
2021VFR	20" - 26"	Operator, High Height, Foot Ring

Item Number	Height Range	Description
2021B	20" - 26"	Operator, High Height, Backrest
2021BV	20" - 26"	Operator, High Height, Backrest
2021BFR	20" - 26"	Operator, High Height, Backrest, Foot Ring
2021BVFR	20" - 26"	Operator, High Height, Backrest,
2042LV	20" - 26"	Assistant, Adjustable Body Support on Left
2042RV	20" - 26"	Assistant, Adjustable Body Support on Right
2052LV	20" - 26"	Assistant, Ratcheted Body Support on Left
2052RV	20" - 26"	Assistant, Ratcheted Body Support on Right
2052BLV	20" - 26"	Assistant, Backrest, Ratcheted Body Support, Left
2052BRV	20" - 26"	Assistant, Backrest, Ratcheted Body Support, Right
11001V-D	18" - 24"	Operator
11001BV-D	18" - 24"	Operator, Backrest
11001VFR-D	20" - 27.5"	Operator, Foot Ring
11001FR-D	20" - 27.5"	Operator, Foot Ring
11001BFR-D	20" - 27.5"	Operator, Backrest, Foot Ring
11001BVFR-D	20" - 27.5"	Operator, Backrest, Foot Ring
11020-D	20" - 27.5"	Operator
11001FRA-D	20" - 27.5"	Assistant, Swing Arm
11001VFRA-D	20" - 27.5"	Assistant, Swing Arm
101282	-	Set of Stationary Arms
3102858	-	Casters - 135 Series Stools - Set of 5
3102827	-	Casters - 9XXX Series Stools - Set of 5
3102278	-	Casters - 33XX Series Stools - Set of 5
3102275	-	Casters - 31XX and 22XX Series Stools - Set of 5
3102274	-	Casters - 11XX Series Stools - Set of 5



Through better ergonomic seating solutions for you and your staff, Brewer helps practices achieve smarter economics.

In addition, our extensive line of solutions is manufactured to our rigorous quality standards and backed by best-in-class standard warranties. It's all part of our commitment to provide you with The Power to Advance.

For more information about the many ways we deliver The Power to Advance, contact your Brewer representative today!



The Brewer Company, LLC
N88 W13901 Main Street, Suite 100
Menomonee Falls, WI 53051

1.888.273.9371
www.brewercompany.com

