



Brewer

Dental Seating

Ergonomically  
designed to  
support the health  
of your practice.

*Delivering the Power of Neutral Posture*

## Musculoskeletal Disorders are a pervasive problem for dental professionals.

*MSDs can cost a practice more than \$67,000 per incident.<sup>1</sup> They have been documented to affect more than 80% of all dental professionals and are known to shorten careers and erode a practice's bottom line. Prolonged sitting in poorly designed stools will result in muscular fatigue and can contribute to musculoskeletal pain and injuries.*

### Proper posture is critical to the long-term health and productivity of your dental team.


The risk of injuries involving the back, neck and shoulders can be greatly reduced by working in Neutral Posture – a balanced, upright position with your weight evenly distributed. In addition, appendages are neither moved away from the body's midline, nor laterally turned or twisted.



*Model 9000BST shown in Black Satin upholstery with arms (part #101282) sold separately*

### Designed to make a difference.

Working in conjunction with Marquette University, Brewer has analyzed working conditions in the dental operator. Through extensive testing and observation, we have pinpointed the proper seating mechanics needed to more effectively meet the strenuous demands of dental professionals. The results of these efforts are a line of highly advanced, ergonomic seating solutions specially designed to help safeguard the very backbone of your practice.



*According to the ADA, over 300,000 injuries can be prevented each year with ergonomic seating.*

### So be sure to take a good look at all the ways Brewer adds value to your practice.

#### TABLE OF CONTENTS

<b>Premium Ergonomic Seating</b>	
9500 Series .....	6
9000 Series .....	8
9100 Series .....	8
9200 Series .....	9
<b>Dynamic Saddle Seat</b>	
135AS Series .....	11
135JS Series .....	11
<b>Performance Ergonomic Seating</b>	
3100 Series .....	12
3300 Series .....	12
<b>Contemporary Seating</b>	
2000 Series .....	14
<b>Traditional Seating</b>	
11001 Series .....	15

Design of a stool to improve the working posture of dentists and dental hygienists – Marquette University grant #59121, December 1995.

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1

<sup>1</sup> Electro Kinetic Technologies; How Much Does a Musculoskeletal Injury Cost?; Jan. 2016.

*Advancing Ergonomic Innovation to Enhance Your Comfort*

## Specializing in breakthroughs that matter.

*From carefully researched design refinements to entirely new concepts in seating, Brewer is constantly working to offer you the most advanced, most comfortable seating solutions in the industry.*

*Most are features we introduced to the market, and many remain exclusives to this day. One thing they all have in common is our commitment to provide you with the ultimate dental seating experience to keep you more productive all day, every day.*

### **Dual Dynamic Motion Technology.**

Brewer is the first in the industry to incorporate Dynamic Motion Technology into both the seat and backrest. Now, you get two ways to engage your core as you work to promote back health.

### **Patented, Dynamic Saddle Technology.**


This Brewer exclusive combines a seat design that naturally supports neutral posture with full 360° tilt to help strengthen your core as you work.

### **Seat tilt and lumbar support.**

Helps maintain the ideal low-back curve, decrease disc pressure\*, enable closer positioning to the patient with less muscle strain and may help reduce low-back pain.


### **Unmatched pressure management & comfort.**

HybriGel™ Foam is a Brewer exclusive in the industry. This gel-infused memory foam supports the body for superior pressure relief, while its open cell structure dramatically increases airflow for greater cooling comfort.




*Model 9120BL shown in Deep Sapphire upholstery*

*Pressure Mapping*



*Without HybriGel - Avg. Pressure = 49*



*With HybriGel - Avg. Pressure = 36.9*

\*Journal of Medicine and Life: The dentist's operating posture – ergonomic aspects. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4151237/>

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1

Premium Ergonomic Seating – 9500 Series

Redefining comfort with latest & greatest ergonomic advances.

The 9500 Series integrates the best of Brewer’s proven ergonomic expertise with the latest design advancements in the industry. It includes the introduction of our third generation of Dynamic Motion Technology and an innovative wrap-around backrest. Everything about this stool contributes to facilitating freedom of movement – for enhanced comfort, reduced fatigue and optimum productivity. All models feature a 250-lb weight capacity and a standard 5-year warranty.

**9500 Series**

The evolution of ergonomic comfort continues with the state-of-the-art in the fight against back pain.

**Dentist/Hygienist Models**

- 9500B - Stitched upholstery – Standard height
- 9500BV - Seamless vacuum-formed upholstery – Standard height
- 9500B-L - Stitched upholstery – Low height
- 9500BV-L - Seamless vacuum-formed upholstery – Low height
- 9500B-H - Stitched upholstery – High height
- 9500BV-H - Seamless vacuum-formed upholstery – High height
- 9500BFR - Stitched upholstery – Adjustable foot ring
- 9500BVFR - Seamless vacuum-formed upholstery – Adjustable foot ring

**Assistant Models**

Feature contoured, ratchet-adjustable body support with HybriGel™ Foam and adjustable foot ring standard.

- 9520BL - Stitched upholstery – Left support
- 9520BLV - Seamless vacuum-formed upholstery – Left support
- 9520BR - Stitched upholstery – Right support
- 9520BRV - Seamless vacuum-formed upholstery – Right support

**Accessories**

- 101282 - Optional stool arms (sold separately) – shown on pages 2-3

**Cushions**

Backrest 17" W x 9.5" H, Seat 19.5" W x 15" D

**Chair Heights**

Standard, 19.5"-25.5" Low, 17"-21" High, 22"-29"

Dentist/Hygienist foot ring, 23"-33",

Assistant foot ring, 25"-33"



Model 9500B in Wedgewood



Model 9520BL in Black Satin



Model 9500B shown in Black Satin upholstery

Exclusive use of HybriGel™ Foam in seat (and body support of assistant models), delivers unmatched pressure management and cooling capabilities to optimize comfort.

Flexible, waterfall seat front for greater ease of movement and improved comfort.

Easily adjust desired level of synchronized seat & backrest tilt to help maintain neutral posture.

Pneumatic height adjustment makes it easy to get yourself into the ideal neutral posture stance.

New, low-profile wrap-around lumbar support cradles you in comfort, while allowing a generous range of motion.

Durable, high-compression rubber spacers strategically placed in seat & backrest to enable micro-movements.

Backrest adjusts independently at the touch of your fingertips to comfortably allow greater range of movement.

B - Backrest V - Vacuum Formed Seamless Seat & Backrest FR - Adjustable Foot Ring

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1



Visit [brewercompany.com](http://brewercompany.com)

Premium Ergonomic Seating – 9000, 9100 & 9200 Series

In sync with your comfort and productivity needs.

Unlike most dental stools, our Premium Ergonomic Series seating features continuous, synchronized movement of the saddle-style seat and lumbar support. This combination allows for proper weight distribution and encourages neutral posture throughout your typical range of motion to help prevent fatigue and promote muscular health. All our Premium Ergonomic Seating now features independent backrest adjustment at the touch of your fingertips to comfortably allow a greater range of movement. Plus, all models feature our exclusive HybriGel™ Foam for unmatched pressure management and cooling comfort, a 250-lb weight capacity, and a standard 5-year warranty.

**9000 Series**

The original ergonomic stool offers thick, plush upholstery with a deeper, cut out seat and full-size backrest.

**Dentist/Hygienist Model**

Available in your choice of column heights to fit the seat to your stature for maximum comfort.

**9000B** - Seamless upholstery – Standard height

**9000BST** - Stitched upholstery – Standard height

**9000B-L** - Seamless upholstery – Low height

**9000BST-L** - Stitched upholstery – Low height

**9000B-H** - Seamless upholstery – High height

**9000BST-H** - Stitched upholstery – High height

**9000BFR** - Seamless upholstery –

Adjustable foot ring

**9000BFRST** - Stitched upholstery –

Adjustable foot ring

**Assistant Models**

Feature contoured, ratchet-adjustable body support with HybriGel™ Foam and adjustable foot ring standard.

**9020BL** - Seamless upholstery – Left support

**9020BLST** - Stitched upholstery – Left support

**9020BR** - Seamless upholstery – Right support

**9020BRST** - Stitched upholstery – Right support

**Accessories**

**101282** - Optional stool arms

(sold separately) – shown on pages 2-3

**Cushions**

9000 Series – Backrest 19.5" W x 15" H, Seat 19" W x 16.5" D 9100 Series – Backrest 11.5" W x 17.5" H, Seat 19.5" W x 15" D

**Chair Heights**

Standard, 19"-25" Low, 16.5"-20.5" High, 21.5"-28.5" Dentist/Hygienist foot ring, 22.5"-32.5", Assistant foot ring, 24.5"-32.5"

**9100 Series**

Features all the benefits of our 9000 Series while offering a slightly less deep seat and sleeker backrest designed to allow a freer range of movement.

**Dentist/Hygienist Model**

Available in your choice of column heights to fit the seat to your stature for maximum comfort.

**9100B** - Stitched upholstery – Standard height

**9100BV** - Seamless vacuum-formed upholstery – Standard height

**9100B-L** - Stitched upholstery – Low height

**9100BV-L** - Seamless upholstery – Low height

**9100B-H** - Stitched upholstery – High height

**9100BV-H** - Seamless upholstery – High height

**9100BFR** - Stitched upholstery –

Adjustable foot ring

**9100BFRST** - Seamless upholstery –

Adjustable foot ring

**Assistant Models**

Feature contoured, ratchet-adjustable body support with HybriGel™ Foam and adjustable foot ring standard.

**9120BL** - Stitched upholstery – Left support

**9120BLV** - Seamless upholstery – Left support

**9120BR** - Stitched upholstery – Right support

**9120BRV** - Seamless upholstery – Right support

**Accessories**

**101282** - Optional stool arms

(sold separately) – shown on pages 2-3

**9200 Series**

Offers the best of both worlds, pairing the deeper seat from the 9000 Series with the smaller, sleek backrest from the 9100 Series.

**Dentist/Hygienist Model**

Available in your choice of column heights to fit the seat to your stature for maximum comfort.

**9200B** - Seamless upholstery – Standard height

**9200BST** - Stitched upholstery – Standard height

**9200B-L** - Seamless upholstery – Low height

**9200BST-L** - Stitched upholstery – Low height

**9200B-H** - Seamless upholstery – High height

**9200BST-H** - Stitched upholstery – High height

**9200BFR** - Seamless upholstery –

Adjustable foot ring

**9200BFRST** - Stitched upholstery –

Adjustable foot ring

**Assistant Models**

Feature contoured, ratchet-adjustable body support with HybriGel™ Foam and adjustable foot ring standard.

**9220BLST** - Stitched upholstery – Left support

**9220BL** - Seamless vacuum-formed upholstery – Left support

**9220BRST** - Stitched upholstery – Right support

**9220BR** - Seamless vacuum-formed upholstery – Right support

**9220BRV** - Seamless upholstery – Right support

**Accessories**

**101282** - Optional stool arms

(sold separately) – shown on pages 2-3

9200 Series – Backrest 11.5" W x 17.5" H, Seat 19" W x 16.5" D

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1



Model 9000B in Purple Gray



Model 9020BL in Desert Rose



Model 9100B in Cocoa



Model 9120BL in Azure Blue



Model 9200B in Dusty Jade



Model 9220BL in Graystone

B - Backrest V - Vacuum Formed Seamless Seat & Backrest FR - Adjustable Foot Ring



Visit [brewercompany.com](http://brewercompany.com)

## Dynamic-Motion Saddle-Seat Technology

### Exclusive capability promotes back health.

Brewer is the only U.S. manufacturer to offer 360° Dynamic Motion technology. This ergonomic capability is designed to engage core muscles in order to help strengthen back muscles over time. As a result, it helps keep dentists and hygienists more comfortable and productive throughout the day. We're such big believers in the power of 360° Dynamic Motion that we've expanded our offering so more dental professionals can take advantage of this transformative technology.

**360° Dynamic-Motion technology** allows the freedom of movement that helps engage core muscles. Adjustable design enables you to limit the range of motion initially. Then, gradually increase the range of motion over time to continue strengthening core muscles, supporting enhanced comfort, daily productivity and career longevity.

All Brewer saddle seats feature a 250-lb weight capacity and a standard 5-year warranty.

#### Create an ideal personal fit.

Our latest saddle seats are available in two seat sizes and a choice of seating height ranges to assure a comfortable, personal fit for virtually any dental professional. All feature a comfy, ultra-cushioned seat in your choice of Ultraleather™ upholstery in Raven Wing, Dove Gray or Papyrus\*. The ergonomically shaped saddle tilts the pelvis naturally to improve spinal posture and enhance long-term comfort.



Model 135ASP Amazone shown in Papyrus upholstery



Model 135JSB Jumper shown in Raven Wing upholstery

#### Amazone Series

Features a narrower seat and a choice of two height ranges specifically designed to accommodate smaller-framed dental professionals.



Model 135ASB Amazone in Raven Wing

#### Models

135 Amazone Models*	Adjustable Height Range	Recommended User Height
■ 135ASG	24"-32"	5' 3"-5' 11"
■ 135ASB	24"-32"	5' 3"-5' 11"
■ 135ASP	24"-32"	5' 3"-5' 11"
■ 135ASG-L	21"-27"	Under 5' 3"
■ 135ASB-L	21"-27"	Under 5' 3"
■ 135ASP-L	21"-27"	Under 5' 3"

#### Jumper Series

Features a wider seat and a choice of two height ranges specifically designed to accommodate taller, larger-framed dental professionals.



Model 135JSG Jumper in Dove Gray

#### Models

135 Jumper Models*	Adjustable Height Range	Recommended User Height
■ 135JSG	24"-32"	5' 3"-5' 11"
■ 135JSB	24"-32"	5' 3"-5' 11"
■ 135JSP	24"-32"	5' 3"-5' 11"
■ 135JSG-H	27"-37"	Over 5' 11"
■ 135JSB-H	27"-37"	Over 5' 11"
■ 135JSP-H	27"-37"	Over 5' 11"

\*3 Ultraleather options to choose from"

\* G - Dove Gray B - Raven Wing P - Papyrus

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1

**Brewer**

Visit [brewercompany.com](http://brewercompany.com)

Premium Ergonomic Seating – 3100 & 3300 Series

A practical blend of ergonomics, comfort and durability.

All our Performance Ergonomic Seating is designed to be a welcome addition to your practice for years to come. Our contoured 17" wide x 19" deep seat cushion features a flexible, waterfall front edge for enhanced comfort. Dentist's & hygienist's models include an adjustable, contoured 13" high x 15" wide backrest – standard. To hold up to demanding work days, all come equipped with a 24" diameter, powder-coated, 5-leg cast-aluminum base and hooded, dual-wheel casters – complete with polyurethane tread to protect your floor surfaces. Each offers a 250-lb weight capacity and is backed by a standard 5-year warranty for added peace of mind.

**3100 Series**

Available in three design configurations and three different adjustable pneumatic height ranges to meet the needs of your team.

**Dentist/Hygienist Models**

**3125B** - Pneumatic height adjustment from 18"-24"; dual adjustable, contoured, tilting backrest

**3125BV** - Vacuum-formed upholstery

**3125BFR** - Pneumatic height adjustment from 20.5"-27.5"; dual adjustable backrest & adjustable foot ring

**3125BFRV** - Vacuum-formed upholstery

**Assistant Models**

**3145L** - Adjustable left body support with patented ratchet; pneumatic height adjustment from 23"-31"; and adjustable foot ring standard

**3145LV** - Vacuum-formed upholstery

**3145R** - Adjustable right body support with patented ratchet; pneumatic height adjustment from 23"-31"; and adjustable foot ring standard

**3145RV** - Vacuum-formed upholstery



Model 3125B in Dusty Jade



Model 3125BFR in Royal



Model 3145L in Sage

**3300 Series**

Available in five different design configurations and three different adjustable pneumatic height ranges to meet the needs of your team.

**Dentist/Hygienist**

**3335B** - Pneumatic height adjustment from 18.5"-24.5"; synchronized seat & backrest tilt

**3335BV** - Vacuum-formed upholstery

**3335BA** - Pneumatic height adjustment from 18.5"-24.5"; synchronized seat & backrest tilt; stationary armrests standard

**3335BVA** - Vacuum-formed upholstery

**3335BFR** - Pneumatic height adjustment from 21"-28"; synchronized seat & backrest tilt; adjustable foot ring standard

**3335BFRV** - Vacuum-formed upholstery

**Assistant Models**

**3345L** - Adjustable left body support with patented ratchet; pneumatic height adjustment from 23"-31"; and adjustable foot ring standard

**3345LV** - Vacuum-formed upholstery

**3345R** - Adjustable right body support with patented ratchet; pneumatic height adjustment from 23"-31"; and adjustable foot ring standard

**3345RV** - Vacuum-formed upholstery



Model 3335B in Mushroom



Model 3335BA in Black Satin



Model 3345L in Sand



Model 3335B shown in Graystone upholstery

**Assistant Models cont.**

**3345BL** - Contoured backrest; adjustable left body support with patented ratchet; pneumatic height adjustment from 23"-31"; and adjustable foot ring standard

**3345BLV** - Vacuum-formed upholstery

**3345BR** - Contoured backrest; adjustable right body support with patented ratchet; pneumatic height adjustment from 23"-31"; and adjustable foot ring standard

**3345BRV** - Vacuum-formed upholstery

**3355BL** - Synchronized seat and backrest tilt; adjustable left body support with patented ratchet; pneumatic height adjustment from 24"-32"; and adjustable foot ring standard

**3355BLV** - Vacuum-formed upholstery

**3355BR** - Synchronized seat and backrest tilt; adjustable right body support with patented ratchet; pneumatic height adjustment from 24"-32"; and adjustable foot ring standard

**3355BRV** - Vacuum-formed upholstery



Model 3345BL in Dove



Model 3355BL w/ seat tilt in Royal

B - Backrest V - Vacuum Formed Seamless Seat & Backrest FR - Adjustable Foot Ring

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1

**Brewer**

Visit [brewercompany.com](http://brewercompany.com)

Contemporary Seating – 2000 Series

Brewer quality & comfort in our most economical offerings.

Our Contemporary & Traditional Series offer budget-conscious practices dependable solutions to their comfort and quality needs. Both feature round, 16" wide, 4" thick poly-foam seats for comfort, a durable 5-leg base and hooded, dual-wheel. Each offers a 250-lb. weight capacity and is backed by a standard 5-year warranty.

**2000 Series**

Feature a 24" diameter, powder-coated base and polyurethane treads on the casters.

**Dentist/Hygienist Models**

- 2020 - Pneumatic height adjustment from 18.5"-24.5"
- 2020V - Pneumatic height adjustment from 18.5"-24.5", with seamless upholstery
- 2021 - Pneumatic height adjustment from 21"-28"
- 2021V - Pneumatic height adjustment from 21"-28", with seamless upholstery
- 2021FR - Pneumatic height adjustment from 21"-28", with adjustable foot ring
- 2021VFR - Pneumatic height adjustment from 21"-28", with seamless upholstery & adjustable foot ring
- 2020B - Pneumatic height adjustment from 18.5"-24.5", with contoured backrest
- 2020BV - Pneumatic height adjustment from 18.5"-24.5", with seamless upholstery & contoured backrest
- 2021BFR - Pneumatic height adjustment from 21"-28", with adjustable foot ring & contoured backrest
- 2021BVFR - Pneumatic height adjustment from 21"-28", with seamless upholstery, adjustable foot ring & contoured backrest



Model 2020 in Azure Blue



Model 2020B in Feather

**Assistant Models**

- All feature pneumatic height adjustment from 22"-28" & adjustable foot ring standard
- 2042L - 360° adjustable left body support arm
  - 2042LV - 360° adjustable left body support arm, with seamless upholstery
  - 2042R - 360° adjustable right body support arm
  - 2042RV - 360° adjustable right body support arm, with seamless upholstery
  - 2052L - Patented ratcheting left body support arm
  - 2052LV - Patented ratcheting left body support arm, with seamless upholstery
  - 2052R - Patented ratcheting right body support arm
  - 2052RV - Patented ratcheting right body support arm, with seamless upholstery
  - 2052BL - Patented ratcheting left body support arm, with contoured backrest
  - 2052BLV - Patented ratcheting left body support arm, with seamless upholstery & contoured backrest
  - 2052BR - Patented ratcheting right body support arm, with contoured backrest
  - 2052BRV - Patented ratcheting right body support arm, with seamless upholstery & contoured backrest



Model 2042L in Dove



Model 2052L in Dusty Jade



Model 2052BR in Dusty Jade

Traditional Seating – 11001 Series



Model 2020B shown in Taupe upholstery

**11001 Series**

Feature a 23" diameter, brushed cast aluminum base.

**Dentist/Hygienist Models**

- 11001-D - Pneumatic height adjustment from 18"-24"
- 11001V-D - Pneumatic height adjustment from 18"-24", with seamless upholstery
- 11020-D - Pneumatic height adjustment from 21"-28"
- 11001FR-D - Pneumatic height adjustment from 21"-28" & adjustable foot ring standard
- 11001VFR-D - Pneumatic height adjustment from 21"-28", with seamless upholstery & adjustable foot ring standard
- 11001B-D - Pneumatic height adjustment from 18"-24" with contoured backrest
- 11001BV-D - Pneumatic height adjustment from 18"-24", with seamless upholstery & contoured backrest
- 11020B-D - Pneumatic height adjustment from 21"-28" with contoured backrest
- 11001BFR-D - Pneumatic height adjustment from 21"-28", adjustable foot ring standard & contoured backrest
- 11001BVFR-D - Pneumatic height adjustment from 21"-28", with seamless upholstery, adjustable foot ring standard & contoured backrest



Model 11001-D in Black Satin



Model 11001FR-D in Cloud Gray



Model 11001B-D in Bone



Model 11001BV-D in Bone



Model 11001FRA-D in Dove

**Assistant Models**

- 11001FRA-D - Pneumatic height adjustment from 21"-28" with adjustable foot ring standard
- 11001VFRA-D - Pneumatic height adjustment from 21"-28", with seamless upholstery & adjustable foot ring standard

Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1



Visit [brewercompany.com](http://brewercompany.com)

B - Backrest V - Vacuum Formed Seamless Seat & Backrest FR - Adjustable Foot Ring





Through better ergonomic seating solutions for you and your staff, Brewer helps practices achieve smarter economics.

In addition, our extensive line of solutions is manufactured to our rigorous quality standards and backed by best-in-class standard warranties. It's all part of our commitment to provide you with The Power to Advance.

For more information about the many ways we deliver The Power to Advance, contact your Brewer representative today.

**Brewer**

The Brewer Company, LLC  
N88 W13901 Main Street, Suite 100  
Menomonee Falls, WI 53051

1.888.273.9371  
[www.brewercompany.com](http://www.brewercompany.com)

