

9500 Series Premium Ergonomic Seating

Specifications

Discover greater freedom of movement feels with the all-new 9500 Series from the ergonomic experts at Brewer. This evolution in comfort has been developed to enable you to customize your seating around both preference and ergonomic needs.

Brewer Differentiators

- 3rd generation Dynamic Motion Technology
- Low-profile, wrap-around backrest
- Separate, independent backrest tilt
- Flexible waterfall seat front for added comfort and circulation
- Exclusive use of HybriGel[™] Foam

Product Benefits

- Our 3rd generation Dynamic Motion Technology in the seat & backrest engages your core, helping strengthen back muscles as you work to enhance comfort and stamina
- Low-profile, wrap-around support cradles you in comfort, while allowing you to enjoy a greater range of motion
- Exclusive use of HybriGel[™] Foam for unmatched pressure management and cooling comfort
- Designed to provide correct weight distribution, pelvic tilt, and comfort throughout a broad range of routine motion



- Continuous, synchronized movement of the saddle-style seat and lumbar support allows for proper weight distribution, encourages neutral posture and promotes muscular health
- Adjust desired level of seat and back tilt together to accommodate front-to-back reaching without sacrificing lumbar support
- Backrest adjusts independently at the touch of your fingertips, allowing more comfort, improved sight lines and patient access
- Flexible waterfall seat front features less depth than standard seats to meet individual comfort preferences



9500 Series Premium Ergonomic Seating

Specifications

All 9500 Series Premium Ergonomic Seating Feature:

- Less deep, saddle-style seat dimensions: 18" x 16"
- Includes exclusive use of HybriGel[™] Foam in seat (and body supports in assistant's models)
- Wrap-around backrest with ratcheted height adjustment
- Paddle control for seat height adjustment
- Paddle control enables back support to move independently of the seat

- Paddle control adjusts level of synchronized tilt for seat and back
- Backrest dimensions: 9" x 16"
- 25" diameter powder-coated, 5-leg cast aluminum base with hooded, dual wheel, polyurethanetread casters
- 250 lb weight capacity
- Standard 5 year warranty
- Made in the USA

Dentist/Hygienist Models:

- Four optional column height ranges:
 - 15.5" 20"
 - -18.5'' 24.5''
 - 20.5" 26.5"
 - -22" 31" (w/footring)



Part Number	Height Range	Seamless Upholstery	Foot Ring	Body Support
9500B	18.5"-24.5"	-	_	_
9500BV	18.5"-24.5"	Yes	-	_
9500B-L	15.5"-20"	_	-	_
9500BV-L	15.5"-20"	Yes	_	_
9500B-H	20.5"-26.5"	-	_	_
9500BV-H	20.5"-26.5"	Yes	_	_
9500BFR	22"-31"	_	Yes	_
9500BVFR	22"-31"	Yes	Yes	_

Stationary arms (item number 101282) sold separately

Assistant Models:

- Contoured, ratcheted, adjustable body support with HybriGelTM Foam
- Left and right body support models available
- Adjustable footring standard
- Height: 24" 32"



Part Number	Height Range	Seamless Upholstery	Foot Ring	Body Support
9520BL	24"-32"	-	Yes	Left
9520BLV	24"-32"	Yes	Yes	Left
9520BR	24"-32"	_	Yes	Right
9520BRV	24"-32"	Yes	Yes	Right