

# Brewer

## 135 Series Saddle Seats


Exclusive 360° Dynamic Motion technology optimizes comfort

### Designed to provide unmatched comfort

- Innovative 360° Dynamic Motion capability engages core muscles, helping strengthen back muscles over time – reducing back pain and sitting fatigue
- Saddle design eases load on back muscles by tilting pelvis forward to enhance posture, naturally align vertebrae and disks, and relieve shoulder and neck tension
- Enables dentists and hygienists to be more comfortable and productive throughout their workday

### Personalized comfort

- Choice of saddle seat sizes and cylinder heights enables you to customize the fit for all users under 5' 3", those that are 5' 3"-5' 11", and those over 5' 11"
- Conveniently adjust cylinder height to exact user preference



## 135 Series Saddle Seat

Brewer exclusives help optimize a practice's bottom-line

### Brewer Differentiators

- *360° Dynamic Motion technology promotes back health*
- *Conveniently adjust the level of dynamic motion desired to suit personal preference*
- *Standard 5-year warranty*

### 135 Amazone Series

#### All 135 Amazone Seats feature:

- Narrower 13" upholstered seat design
- Press on ring for easy height adjustment
- Turn ring to conveniently adjust level of dynamic motion
- UltraLeather upholstery available in Dove Gray, Raven Wing and Papyrus
- 23" diameter cast aluminum base w/large castors
- 250 lb. weight capacity



135 Amazone Models	Adjustable Height Range	Recommended User Height
135AS	24"-32"	5' 3"-5' 11"
135AS-L	21"-27"	Under 5' 3"

### 135 Jumper Series

#### All 135 Jumper Seats feature:

- Standard upholstered seat design
- Press on ring for easy height adjustment
- Turn ring to conveniently adjust level of dynamic motion
- UltraLeather upholstery available in Dove Gray, Raven Wing and Papyrus
- 23" diameter cast aluminum base w/large castors
- 250 lb. weight capacity



135 Jumper Models	Adjustable Height Range	Recommended User Height
135JS	24"-32"	5' 3"-5' 11"
135JS-H	27"-37"	Over 5' 11"

### 135DSS Series

#### All 135DSS Saddle Seats feature:

- Innovative split-saddle-seat design with adjustable dynamic movement offering full 360° tilt
- Available in a wide variety of color options
- 23" diameter cast aluminum base w/large casters
- 250 lb. weight capacity



135DSS Models	Adjustable Height Range	Recommended User Height
135DSS	22.5"-29.75"	5' 3"-5' 11"
135DSS-L	18.75"-23.75"	Under 5' 3"
135DSS-H	24"-33.75"	Over 5' 11"