

# Brewer

## Access<sup>®</sup> High-Low Exam Tables

### Specifications

#### Brewer Differentiators

- *Industry-leading 700-lb. (317 kg) weight capacity for Model 6800/6801*
- *Combination of 18" (45.7 cm) low height in seated position and safety grab bars that rotate nearly 180° enables more confident transfers*
- *Exclusive SafeGlide™ leg extension allows caregivers to slide the heavy-duty extension in and out with minimal effort*
- *Innovative pass-through drawers and work surface offer access from either side*
- *Made in USA*

#### Product Specifications

- Weight capacity: 700 lbs. (317 kg) for Model 6800/6801, 450 lbs. (204 kg) for Model 6000/6001/6500/6501
- Seat Height: 18"-37" (45.7cm-94 cm)
- Length: 57"; extendable to 78" (144.8 cm-198.1 cm)
- SafeGlide™ leg extension length: 21" (53.3 cm)
- Lies flat at 28" (71.1 cm)
- Upholstery width: 28" (71.1 cm)
- Front storage drawer: 19.7" x 19.2" x 4" (50.3 cm x 48.8 cm x 9.7 cm)
- Side storage drawer: 22.6" x 20.2" x 7" (57.5cm x 51.4cm x 17.8cm)
- Paper storage capacity: 23"x 23"x 3.25" (58 cm x 58 cm x 8 cm)
- Slide through paper holder accepts 21" & 18" rolls (53.3 cm & 45.7 cm)



Model 6801 in Slate Blue

#### Standard Features

- Industry-best 3-year warranty
- Patient assist grab bars
- Adjustable stirrups
- Laptop work surface
- Easy on/off for patients
- Easy-to-clean, seamless upholstery
- Return-to-chair feature on power back models
- Adjustable backrest for easy patient positioning (pneumatic or power back available)
- Sensors to minimize possibility of injuries during table movement
- Anti-tamper on/off switch prevents accidental movement
- Foot control
- Low voltage actuators
- Single power cord design

