

Brewer

Assist & AssistPRO

POWER PROCEDURE TABLE

ADVANCED PATIENT CARE FOR
THE HEALTH OF YOUR PRACTICE.

- Preserve patient dignity and comfort during procedures, increasing your patient satisfaction rankings.
- Maximize throughput by having accessories and simple one-touch programming at clinician fingertips.
- ADA compliant to accommodate multiple patient populations.
- Reduce your total cost of ownership with our industry-leading **three year warranty** for peace of mind.



Assist & AssistPRO

POWER PROCEDURE TABLE

Make procedures easier, safer and more comfortable for staff and patients.



PROGRAMMABILITY

- AssistPRO - simple two-step programming process
- Memory for storing up to four preferred positions
- Visual and audible signals with auto-run for safety

MORE FLEXIBILITY

- Supine, left lateral plane, lithotomy and proctology positions standard
- Patented 3-point pivot headrest system
- Extendable leg rest

ADVANCED SAFETY

- 19" low height for easy wheelchair transfers
- 40" height and positioning for efficient procedures
- Pinch-free design and safety switches

PATIENT COMFORT

- Seamless, multi-density, bolstered support
- ADA compliant to accommodate many patient populations

SUPERIOR QUALITY

- Industry-leading **three year warranty**, standard
- UL 60601 certified
- **Made in USA**

models

#7000	Assist Power Procedure Table
#7500	AssistPRO Programmable Power Procedure Table

standard features

Three-year warranty
450 pound weight capacity
Height: 19" - 40"
Length: 69", extends to 80"
Width: 28"
Tilt range/Trendelenburg: to 30°
UL 60601-1 Certified

accessories

#99497	grab bars**
#99499	articulating armboard**
#99501	side seat rails (pair)
#99502	one universal socket clamp
#99503	standard arm board**
#99504	delux articulating knee crutches (pair) *
#99505	stainless steel treatment pan
#100253	urology drain pan system
#102025	Brewer Welch Allyn® IV bracket (plant install)
#102025	Brewer Welch Allyn® IV bracket (field install)

* Needs one order of #99501

** Needs one order of #99501 and two of #99502

