

# Brewer

# *FLEX*

## Access Exam Table



### ADVANCED STAFF SAFETY AND CARE FOR MORE OF YOUR PATIENTS

- Perform a wider range of examinations than would be possible with non-powered tables and see more patients per day.
- Improve clinical throughput and achieve bariatric excellence with FLEX's industry leading 700lbs. patient weight capacity.
- Meet ADA standards with our 18" low height for safe lateral wheelchair transfers.
- Unmatched ergonomic patient access helps keep you and your staff healthy and more productive.



# FLEX

## Access Exam Table

The Brewer **FLEX** Access Exam Table revolutionizes your exam room opportunities.

### CHAIR TO TABLE DESIGN

- Rapidly adjusts from consultation to examination
- 18" low height for easy and safe patient transfers
- 37" maximum seat height

### SUPERIOR FUNCTIONALITY

- Compact, highly flexible power-exam table footprint
- Unmatched ergonomic patient access increases comfort and efficiency
- Made in USA



### ADVANCED SAFETY

- 700lb (317kg) weight capacity
- Anti-tamper on/off switch
- Standard foot control
- UL60601 certified

### ERGONOMIC UPHOLSTERY

- Seamless & easy to clean
- 12 standard colors available
- 4 CAL-133 fire resistant upholstery available

### Models

- #5700 - FLEX Access Exam Table with pneumatic back
- #5701 - FLEX Access Exam Table with pneumatic back, pelvic tilt, front drawer heater, and outlet
- #5800 - FLEX Access Exam Table with power back
- #5801 - FLEX Access Exam Table with power back, pelvic tilt, front drawer heater, and outlet

### Kits and Accessories

#### Kits:

- #2102894 Stirrup for 5700 and 5800 (factory installed)
- #2102898 Stirrup for 5700 and 5800 (field installed)
- #2102893 Stirrup for 5701 and 5801 (factory installed)
- #2102897 Stirrup for 5701 and 5801 (field installed)
- #2102891 Safety grab bars (factory installed)
- #2102895 Safety grab bars (field installed)
- #2102892 Pass-through laptop work surface (plant installed)
- #2102896 Pass-through laptop work surface (field installed)

#### Accessories:

- #98464 Stainless steel treatment pan
- #99502 One universal socket clamp
- #99504\* Deluxe articulating knee crutches
- #100253 Urology drain pan system
- #100407 Front drawer dividers
- #100540 Two support adapters for knee crutches (plant installed)
- #1000550 Two support adapters for knee crutches (field installed)
- #102031 Brewer Welch Allyn® IV bracket (plant installed)
- #102032 Brewer Welch Allyn® IV bracket (field installed)



The Brewer Company LLC  
N88W13901 Main Street, Ste 100  
Menomonee Falls, WI 53051  
1.888.Brewer.1  
brewercompany.com

Brewer is an ISO 13485:2003 certified company.