

Brewer Power Table Highlights & Clinical Differentiators

We welcome the opportunity to showcase the strengths of the Brewer offering. For your convenience, we have provided this brief summary of the distinct advantages our exam and procedure table designs deliver.



Access® High-Low Exam Table
6800/6801 Series & 6500/6501 Series



FLEX™ Access® Exam Table
5700/5701 Series & 5800/5801 Series



Assist Power Procedure Table
7000/7500 Series

Patient Positioning

Safety grab bars — standard — facilitate patient positioning.

Infinitely adjustable foot controls allow caregiver to position patient precisely where needed.

Powered backrest easily adjusts for patient positioning.

Exclusive SafeGlide™ leg extension system extends to 21" to comfortably accommodate taller patients.

Infinitely adjustable foot controls allow caregiver to position patient precisely where needed.

Available with pneumatic backrest or powered backrest to easily adjust for patient positioning.

Exclusive SafeGlide™ leg extension system extends to 21" to comfortably accommodate taller patients.

Safety grab bars available as an option to further facilitate patient positioning.

Infinitely adjustable to meet any procedure needs, including supine, left lateral plane, lithotomy and proctology.

Exclusive, easy-to-adjust 3-point pivot headrest.

Infinitely adjustable foot controls allow caregiver to position patient precisely where needed.

Powered footrest extends an additional 5".

Patient Comfort & Satisfaction

Ergonomically designed to preserve patient dignity, comfortably accommodates the widest range of patients.

Seamless, multi-density support upholstery is standard, enhancing patient comfort.

Ergonomically designed to preserve patient dignity, comfortably accommodates the widest range of patients.

Designed to preserve patient dignity, comfortably supporting patients weighing up to 450 lbs. and it is ADA compliant.

Seamless, multi-density support upholstery is standard, enhancing patient comfort.

Healthcare Provider Ergonomics

Minimize physical caregiver interaction required while bolstering patient independence:

- Convenient, low-effort foot-operated controls.
- Unique, 3-sided support during patient transitions enables confident, well-supported transfers.
- Exclusive SafeGlide™ leg extension system delivers enhanced strength and ease of use.

Offers unmatched access to the patient for increased physician comfort.

Minimize physical caregiver interaction required while bolstering patient independence:

- Convenient, low-effort foot-operated controls.
- Exclusive SafeGlide™ leg extension system delivers enhanced strength and ease of use.
- Optional safety grab bars help create Brewer's unique, 3-sided support for more confident, well-supported transfers.

Standard intuitive foot controls enable caregivers to effortlessly position patient for any procedure.

Bariatric Excellence

Only Brewer offers industry-leading 700-lb. weight capacities to accommodate the broadest range of patients.

Equipped with rotating safety grab bars — standard — to create a highly stable platform for independent patient transfers.

SafeGlide™ leg extension features a heavy-duty design.

Offers industry-leading 700-lb. weight capacities to accommodate the broadest range of patients.

Features an industry-leading 450-lb. weight capacity.

Available with patient safety grab bars to help preserve patient independence and dignity.

Access® High-Low Exam Table
6800/6801 Series & 6500/6501 Series

FLEX™ Access® Exam Table
5700/5701 Series & 5800/5801 Series

Assist Power Procedure Table
7000/7500 Series

Storage Capacity

Convenient covered paper roll storage holds up to six rolls.

Patented pass-through drawer design provides more than 60ft³ of storage capacity as well as convenient access from either side to easily adapt to exam room environment.

Spacious front drawer offers additional storage space.

Front storage drawer can accommodate up to six 18" x 3" paper rolls.

Includes convenient table paper storage.

Accessories/Customization

Available with ergonomic stirrups, 5° pelvic tilt, stainless steel debris tray and 60-Watt drawer heater for Ob/Gyn and Urology specialties.

Further customize the unit with the purchase of the following accessories: knee crutches and stainless steel treatment pan as well as a Welch Allyn® light bracket.

Available with pelvic tilt and drawer heater, front drawer heater and outlet.

Options include safety grab bars, a pass-through laptop work surface and ergonomic stirrups.

Available accessories include stainless steel treatment pan, urology drain pan system, a Welch Allyn® light bracket, deluxe articulating knee crutches and a universal socket clamp.

Adjustable ergonomic stirrups are standard.

Also available with patient safety grab bars, standard or articulating arm board side seat rails, universal accessory clamp deluxe articulating knee crutches, stainless steel treatment pan, urology drain pan system and a Welch Allyn® light bracket.

Patient Safety

Anti-tamper on/off switch prevents accidental movement.

Sensors prevent injury during table movement.

Sensors prevent injury during table movement.

Features pinch-free design and safety switches.

Includes visual and audible signals with auto-run for safety.

Wheel Chair Transfers

Provides ADA-compliant 18" low seat height with unique 3-sided support capability to allow more confident transfers while minimizing the need for direct caregiver interaction.

Provides ADA-compliant 18" low seat height for confident transfers.

Meets the 19", ADA-compliant low seat height, allowing easy lateral wheelchair transfers.

Care Efficiency

Motor speed allows caregivers to safely and efficiently transition patient from consultation to examination mode.

Return to chair feature easily takes table back to chair position with one touch of the controls.

Proprietary, 3-sided-support-configuration capability also helps expedite safe patient transitioning.

Pass-through work surface and generous storage capacity drive clinical efficiency by assuring caregivers have needed supplies and instruments within easy reach.

Delivers the efficiency of power table performance for about the cost of a fixed-height exam table.

Unmatched patient access supports efficient care with less fatigue.

Motor speed allows caregivers to safely and efficiently transition patient from consultation to examination mode.

Return to chair feature easily takes table back to chair position with one touch of the controls.

Maximize throughput with pre-programmable procedure configuration settings for one-touch set-up.

Full 40" height and positioning supports efficient procedures.

Table speed is dependent on the choice of powered actuators selected.

Infection Control

Features easy-to-clean seamless, medicinal stain resistant upholstery.

All exterior surfaces have been designed to facilitate quick, easy cleanup.

Features easy-to-clean seamless, medicinal stain resistant upholstery.

All exterior surfaces have been designed to facilitate quick, easy cleanup.

Features easy-to-clean seamless, medicinal stain resistant upholstery.

All exterior surfaces have been designed to facilitate quick, easy cleanup.

Brewer

The Power to Advance. Safely. Financially. Competitively.
www.brewercompany.com | P: 1.888.273.9371

All brand names and product names are trademarks, registered trademarks, or tradenames of their respective holders.

LIT-M-047