

Product Specifications & Standards

- Industry best 3 year warranty
- Weight capacity:
450lbs. for Models 6000PLUS/6001PLUS/6500PLUS/6501PLUS
700lbs. for Models 6800PLUS/6801PLUS
- Seat Height: 18" - 37"
ADA Access Board Compliant at the seated and table position
- Length: 60"; extendable to 78"
- SafeGlide™ leg extension length: 18"
- Upholstery width: 28"
- Front storage drawer: 19^{3/4}" x 19^{1/4}" x 4"
- Side storage drawer: 22^{1/2}" x 20^{1/4}" x 8"
- Paper storage capacity: 23" x 23^{1/4}" x 3^{1/4}"
- Slide-through paper holder accepts 18" & 21" rolls
- ComfortGrip™ Patient Assist Handles included
- Mayo Pass-Through Work Surface included
- 6 cubic ft. of storage, 3 times more than other brands



EasyGlide™ Mobility System



Urology Drain Pan System



Knee Crutches



Footrest Pad

Optional Accessories

- 3500687PLUS: EasyGlide™ Mobility System
- 99502: One Universal Socket Clamp
- 99504: Deluxe Articulating Knee Crutches*
- 100253: Urology Drain Pan System
- 100407PLUS: Front Drawer Dividers
- 3500077: Footrest Pad
- 100540: Two Support Adapters for Knee Crutches (factory installed)
- 100550: Two Support Adapters for Knee Crutches (field installed)
- 102031PLUS: Brewer Welch Allyn® IV Bracket (factory installed)
- 102032PLUS: Brewer Welch Allyn® IV Bracket (field installed)

* Needs two 99502 and one of either 100540 or 100550

ACCESS[®] HIGH-LOW PLUS

POWER EXAM TABLE



Upholstery Options

- Available in 15 standard colors
- Specialty colors available for an additional cost

Certifications

- UL 60601-1 compliant
- CAN-CSA C22.2 No. 601.1

Additional Specifications for Model 6000PLUS

- Hospital-grade outlet - standard
- 115VAC, 60Hz
- Pelvic tilt
- 60 Watt front drawer heater



Access High-Low PLUS
shown in Blue Fog

Models	Weight Capacity	ADA Compliant 18" Seat Height	Power Base	Power Back	Manual Back	Patient Assist Handles	Stirrups	Pelvic Tilt & Drawer Warmer	Outlet	Work Surface	Mobility
6000PLUS	450	•	•		•	•	•			•	optional
6001PLUS	450	•	•		•	•	•	•		•	optional
6500PLUS	450	•	•	•		•	•			•	optional
6501PLUS	450	•	•	•		•	•	•	•	•	optional
6800PLUS	700	•	•	•		•	•			•	optional
6801PLUS	700	•	•	•		•	•	•	•	•	optional