

Brewer

Versa

Exam Table



Total Storage Capacity

Versa

1 Drawer Front; 1 Side 4.7 ft³

Access 5000

2 Drawer Front; 2 Side 5.2 ft³

Competitive Table

2 Drawer Front; 3 Side 4.2 ft³

* Measurements are approximate



Leading-edge design & added versatility help you enhance care across patient populations

The Versa Exam Table readily and cost-effectively adapts to your changing patient and practice needs.



HIGHLY CUSTOMIZABLE

- Easy-connect, removable upholstery top is available in standard & pediatric configurations
- Add options as needed to economically meet practice requirements



IMPROVED PATIENT EXPERIENCE

- Large, sturdy step with integral, anti-slip surface enhances patient safety
- Stable, cushioned platform with 500-lb patient weight capacity comfortably supports patients of all builds
- Optional multi-position patient-assist handles sized for more confident transitions and easier patient self-positioning
- Seamless upholstery & integral step surface offer easy cleanability and infection control




MORE ERGONOMIC ACCESS

- Rounded corners and softened edges throughout safeguard patients and caregivers, while improving patient access
- Pass-through side drawer purposefully positioned at waist height to reduce bending and twisting
- Front drawer strategically positioned so caregivers do not infringe on a patient's personal space

Available in 12 standard colors

Recommended



24-Black Satin 29-Deep Sea 30-Ivy 32-Blue Fog 31-Sapphire

Also available



21-Clamshell 22-Gunmetal 23-Feather 25-Saddle 26-Tapestry 27-Cabernet 28-Cocoa

We Help Healers Heal

Contact Brewer at 1.888.273.9371 to get a quote or schedule an on-site demo.

Product Specifications

- Weight capacity: 500 lbs. (226 kg)
- Seat height: 32.5" (82.55 cm)
- Length: 58"; extendable to 72" (147.3 cm-182.9 cm)
- Pediatric table length: 51"; expandable to 65" (129.5 cm-165.1 cm)
- Pediatric table width: 26" (66 cm)
- Leg extension length: 14" (35.5 cm)
- Upholstery width: 27.25" (69.2 cm)
- Front storage drawer: 21.3" x 19.75" x 4.25" (59.9 cm x 50.1 cm x 9.7 cm)
- Side storage drawer: 23.6" x 35.75" x 7.6" (57.5 cm x 90.8 cm x 19.3 cm)
- Paper roll holder: accommodates up to 21" paper rolls
- UL certified: IEC60601-1 Ed. 3.1

Standard Features

- Industry-best 5-year warranty
- Easy-connect, removable upholstery top
- Large, sturdy patient step with integral anti-slip surface
- Easy-to-clean, seamless upholstery available in 12 standard color choices
- Numerous additional custom colors available to choose from
- Two side drawer dividers for organization

Optional Features

- Hospital-grade electric receptacle
- Paper fastener strap on standard top offers convenience
- Patient assist handles
- 120VAC, 60Hz
- 60-Watt front drawer heater
- 5° pelvic tilt
- Ergonomic Stirrups
- Footrest pad

Models

- Standard model:
1000AX-XXDX-XXX-__ [2 digit color code]
- Standard model with pelvic tilt, front drawer heater and outlet:
1000AX-XSDT-RHX-__ [2 digit color code]
- Standard model with assist handles:
1000AX-HXDX-XXX-__ [2 digit color code]
- Pediatric model:
1000PX-XXDX-XXX-__ [2 digit color code]

Optional Accessories

- 98464 – Stainless Steel Treatment Pan
- 3500286 – Deluxe articulating knee crutches with assembly kit
- 100253 – Urology Drain Pan System
- 3500288 – Front Drawer Dividers
- 3500347 – Brewer Welch Allyn® Exam Light / IV Bracket (factory installed)
- 3500398 – Brewer Welch Allyn® Exam Light / IV Bracket (field installed)
- 3500237 – Paper Fastener Strap
- 3500077 – Footrest Pad

* Product specification measurements are approximate and may vary +/- 0.5"



Brewer



We Help Healers Heal

Learn more at brewercompany.com

LIT-M-056