Brewer

5 Drawer Mobile Medical Cart

Premium construction, security, and ease-of-use design make this the perfect addition to your clinic.

BREWER DIFFERENTIATORS

- 5 drawer construction offers 5.57 cubic feet of of storage space for equipment and supplies.
- 5-in. diameter full-swivel locking casters make movement throughout the clinic a breeze.
- Premium soft-close return glides offer smooth, silent drawer closure.
- Seamless, easy-to-clean plastic work surface for convenience.

SPECIFICATIONS

Model #3500101 Height: 40 inches Width: 27 inches Depth: 18 inches Drawers: 5 Drawer Dimensions: 15 ½ in. x 23 in. Drawer heights: two 3 in., two 6 in., and one 9 in. Workable Surface Top: 24 in. x 14 1/2 in. Security: Central key lock Wheels: 5-inch full-swivel wheels Braking: 2 wheels equipped with parking brakes Handle: Push handle (can be assembled on either side of cart) Weight: 128 pounds Base / Frame: Steel Work surface top: plastic Warranty: 1 year

Assembly: Base, frame and top come fully assembled. No tools required.

The Brewer Company is an ISO 13485:2016 certified company